







Philips Reusable ECG Lead Sets & Trunk Cables

Instructions for Use


3 ECG	4 ECG	5 ECG	5 (3+2) ECG	6 (4+2) ECG	
					
					
AAMI IEC M1671A, M1675A*, M1672A, M1678A*	AAMI IEC M1532A, M1533A,	AAMI IEC M1976A, M1979A*, M1978A,	AAMI IEC M1968A M1973A* M1971A, M1974A*	AAMI IEC M1680A, M1681A	
					
AAMI IEC M1673A M1674A		AAMI M1602A	AAMI IEC M1644A M1645A		
					
AAMI IEC M1624A M1626A					
					M1663A (5+5) M1665A (6+4) M1667A (6) M1668A (5) M1669A (3) M1949A (5+5) 989803170171 (3)* 989803170181 (5)*




* Products intended for use in
Operating Room (OR) environments

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


ASSOCIATED DOCUMENT - Refer to the Addendum that accompanies this Instructions for Use (IFU) for validated cleaning agents and procedures (ECG Lead Sets, Cables, and Adapter Cables Care, Cleaning, and Disinfection). **INTENDED PURPOSE** - The ECG lead set and trunk cable is an accessory device to be used together with vital parameter monitoring device - as the main device intended to be used for monitoring and recording of, and to generate alarms for, multiple physiological parameters of adults, pediatrics, and neonates. **Intended Use** - Philips IntelliVue reusable ECG lead sets and trunk cables are limited by the indications for use of the connected monitoring and diagnostic equipment, and are for use only by healthcare professionals. They are designed for multiple patient use and indicated for monitoring cardiac signals for both diagnostic and monitoring purposes in adult, pediatric, and neonatal patients. Philips operating room (OR) ECG trunk cables are indicated for use to monitor patient ECG in an electrosurgery (ESU) environment when used in combination with a compatible ECG patient lead set on adult, pediatric, infant, and neonatal ECG applications. **Intended Users** - Healthcare/ Medical professionals only. **Indications for Use** - Philips IntelliVue reusable ECG lead sets/ trunk cables and Philips OR/ ECG trunk cables are indicated for continuous monitoring of cardiac signals for both diagnostic and monitoring purposes. These devices are limited by the indications for use of the connected monitoring and diagnostic equipment in healthcare facilities. These devices are intended to interact with patient intact skin only. **Clinical Benefit** - ECG trunk cables and lead sets are intended to be used with and support the function of the connected monitor. Beyond the indirect clinical benefits of supporting a diagnosis to be made through the intended use of the connected monitor, no additional clinical benefits can be claimed. **Contraindications** - There are no known contraindications. **APPLICABLE STANDARD** - All lead sets and trunk cables included in this IFU meet requirements of the standard, "ECG Trunk Cables and Patient Leads/ (ANSI/AAMI EC53)". **PRODUCT DESCRIPTIONS** - Tables on page i define physical characteristics and applications for all ECG lead set and trunk cable products included in this IFU. Column headings within those tables are defined below.

AAMI	Check mark ✓ indicates that the lead set is AAMI color coded.	IEC	Check mark ✓ indicates that the lead set is IEC color coded.
	Check mark ✓ identifies lead sets that are appropriate for chest application.		Check mark ✓ identifies lead sets that are appropriate for arm/leg application.
3 ECG	Indicates the total number of ECG leads in each lead set, and the type of electrode connector at the end of each lead:  = Grabber  = Snap  = Mini Clip		Check mark ✓ identifies lead sets and/or trunk cables that may be used in operating room environments.
4 ECG			Defines the lead set or trunk cable catalogue reference number.
5 ECG			Identifies length(s) of lead set leads, or trunk cable.
6 ECG			
AAMI + IEC (x) or (x+x)	Defines ECG trunk cables as being both AAMI and IEC compatible, with (x) or (x+x) identifying each trunk cable's lead set connector pin configuration.		

EXPECTED SERVICE LIFE - Expected service life of these ECG lead sets and trunk cables is 18 months of typical clinical use. **COMPATIBILITY** - These ECG lead sets and trunk cables can be used with any ECG monitor or defibrillator for which they are listed as accessories in that product's IFU. Ensure that lead sets are plugged into appropriate trunk cables designed for the same number of leads. Incompatible components may result in performance issues. **OPERATING ROOM ENVIRONMENTS** - Product REF numbers shown with an asterisk * on the front cover of this IFU are intended for use in Operating Room (OR) environments: • M1675SA, • M1678A, • M1979A, • M1973A, • M1974A, • 998003170711, • 998003170181.

WARNINGS - Modification of this equipment is not allowed. • Attach all applicable labels that are shipped with these products. • These products are intended for use by healthcare professionals only. • These lead sets and cable products are to be used only in approved combinations with compatible products and monitoring equipment. • Use of incompatible product combinations including those not defined by the monitoring equipment's IFU must be verified. • Use of incompatible products could lead to performance issues or result in a reduction of safety. • Visit www.philips.com/REACH for information on REACH (Registration, Evaluation, Authorization and Restriction of Chemicals) substances that may be contained in some Philips products. • The patient must be monitored continuously during the procedure. • Do not touch the patient while the system is connected to the patient. • Do not use any lead set or cable when visual inspection reveals teeth separation, label peeling, or other visible cable damage. • If visual inspection reveals that an inner wire of any patient lead is exposed, do not use and immediately discard the lead set. • Ensure that leads are securely connected to the trunk cable, and the trunk cable is properly plugged into the monitor/differential base. Otherwise incorrect patient physiological data may result. • Ensure that the patient is properly grounded during Electro-Stimulus (ESU) procedures to prevent patient/monitor injury (i.e. electrical shock). • Patient leads must be positioned carefully to avoid entanglement, choking, and strangulation. • Before use verify that there are no tears in the patient leads. • Do not use any lead set or cable when visual inspection reveals tears in the patient leads. • Do not use any lead set or cable when visual inspection reveals wires on all patient leads are not exposed; otherwise immediately discard the lead set.

CAUTIONS • Always clean patient cables when they are removed from their packaging for first-time use. • Do not use any cable that shows visual signs of contamination, including liquid or lint contamination in cable connectors. • Do not use in excessively wet environments or under massive influence of fluids (e.g., rain). • Do not immerse and do not autoclave or use ultrasonic cleaners on the cable. • Do not clean cable connector electrical contacts or connectors with bleach. These cables must be cleaned and disinfected using only the validated methods described in the accompanying IFU Addendum (see [Cable Cleaning and Disinfection](#) below). • Do not disconnect a cable by pulling on the lead wires. • These cables are not suitable for use in an MRI environment as there is the potential for patient burns. When monitoring ECG in the OR, always use an OR lead set or an OR trunk cable to provide appropriate ESU protection that prevents patient injury (i.e. burns); respiration measurement is not possible when using an OR lead set or an OR trunk cable. • Caretakers should assess patients for adverse reactions on the skin where the device has contact, such as redness (erythema), swelling (edema), irritation, sensitization (delayed Type IV hypersensitivity), allergy, immune response, or other reactions.

VISUAL INSPECTION PRIOR TO EACH USE – Before attaching any cable to a patient or monitor, perform a visual inspection to determine if the cable has reached end-of-life. Check for cracks, blisters, peeling, exposed wires, damaged connectors, peeling labels, and similar wear or damage that could compromise accurate readings or cause patient/user injury (e.g. cuts). Whenever visual inspection reveals a cable is no longer suitable for continued use, follow appropriate product disposal procedures (see *Product Disposal*). **PATIENT APPLICATION, ECG LEADS** – Refer to the compatible patient monitor I/FU for information regarding proper ECG lead/electrode placement that complies with standard AAMI or IEC practices. **CABLE CLEANING AND DISINFECTION** – Refer to the Addendum that accompanies this I/FU titled *ECG Lead Sets, Cables and Adapter Cables Care, Cleaning, and Disinfection* for validated cleaning and disinfection instructions. **STERILIZATION** – Do not attempt to sterilize any of the cables. The cables are not designed to be sterilized. **REUSE** – Do not reuse any of the cables or the disposable and/or disinfect reusable cables that are applied to any patient, if they become visibly soiled. **STERILIZATION** – All lead sets and cable products described in this I/FU are not intended to be sterilized.

PRODUCT DISPOSAL. Follow approved medical waste disposal methods as specified by your patient care facility or local regulations. **INCIDENT REPORTING.** Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey in which the user and/or patient is established. **ENVIRONMENTAL SPECIFICATIONS.** Temperature, humidity and atmospheric pressure ranges identified below must be maintained for all of the lead set and cable products described in this IFU. Otherwise product damage could result. **Temperature:** Continuous Operating Temperature Range: 0°C to 40°C (32°F to 104°F). Transient Operating Temperature Range: -20°C to 50°C (-4°F to 122°F). Storage/Transport Temperature Range: -20°C to 70°C (-4°F to 158°F). **Humidity:** Minimum humidity (operating/transport/storage): 15% (non-condensing, water vapor pressure up to 50 hPa). Maximum humidity (operating/storage/transport): 95% (non-condensing, water vapor pressure up to 50 hPa).

American Pressure: Operating: 650 hPa; Storage: 100 hPa; 100 hPa – 1080 hPa. **REORDERING INFORMATION:** The following products are CE marked and available in the European Union: **AAIMI Color-Coded ECG Lead Sets** 1.7M, 1M1674 (989803145091) 3 LEADSET, GRABBER, AAIMI, ICU, 1M1675 (989803145131) 3 LEADSET, GRABBER, AAIMI, OR, 1M1676 (989803145111) 3 LEADSET, SNAP, AAIMI, ICU, 1M1624A (989803144941) UNSHIELDED 3 L2 MINICAP AAIMI 0.7M, 1M1523A (989803144491) 4 LEAD SET GRABBER, AAIMI, ICU, 1M1677 (989803145121) 3 LEADSET, GRABBER, AAIMI, ICU, 1M1625A (989803144951) UNSHIELDED 3 L2 MINICAP AAIMI 0.7M, 1M1678 (989803145131) 3 LEADSET, GRABBER, CHEST, AAIMI, OR, 1M1734 (989803125861) 5 LEADSET, GRABBER, AAIMI, OR, 1M1626A (989803144911) 3 LEAD SET SNAP, CHEST, AAIMI, ICU, 1M1644A (989803144991) 5 LEADSET, SNAP, AAIMI, ICU, 1M1680A (989803145161) 6 LEAD SET GRABBER, AAIMI, ICU, **IEC Color-Coded ECG Lead Sets** 1M1672A (989803145101) 3 LEADSET, GRABBER, IEC, OR, 1M1678A (989803145141) 3 LEADSET, GRABBER, IEC, OR, 1M1674A (989803145121) 3 LEADSET, SNAP, IEC, ICU, 1M1626A (989803144951) UNSHIELDED 3 L2 MINICAP AAIMI 0.7M, 1M1675A (989803145131) 3 LEADSET, GRABBER, IEC, OR, 1M1794 (989803125891) 5 LEADSET, GRABBER, CHEST, IEC, ICU, 1M1645A (989803145091) 5 LEAD SET SNAP, IEC, ICU, 1M1681A (989803145171) 6 LEAD SET GRABBER, IEC, ICU, **ECG TRUNK Cables** 1M1669A (989803145071) 3 L ECG ECG TRUNK, AAIMI/IEC 2.7M, 1M1668A (989803145061) 5 LEAD ECG TRUNK, AAIMI/IEC 2.7M, 1M1676A (989803145051) 6 LEAD ECG TRUNK, AAIMI/IEC 2.7M, 1M1633A (989803144791) 10 L ECG ECG TRUNK, AAIMI/IEC 2M, 1M1494A (989803125831) 5+5 ECG TRUNK CABLE AAIMI, OR, 1M1668A (989803145041) 6+4 LEAD ECG TRUNK AAIMI/IEC 2.7M, 989803170171 OR 3-LEAD ECG TRUNK CABLE, AAIMI/IEC, OR 989803170181 OR 5-LEAD ECG TRUNK CABLE, AAIMI/IEC, OR

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Copyright © 2024
All rights reserved
Philips Medizin Systeme Böblingen GmbH