

VariAx[®] 2

Distal Radius Locking Plates

Operative technique

VariAx 2 DR Dorsal Plates

VariAx 2 DR Fragment Specific Plates

VariAx 2 DR Wrist Spanning Plates



This publication sets forth detailed recommended procedures for using Stryker devices and instruments. It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required.

 **WARNING**

- Follow the instructions provided in our cleaning and sterilization guide (OT-RG-1).
- All non-sterile devices must be cleaned and sterilized before use.

 **WARNING**

Multi-component instruments must be disassembled for cleaning. Please refer to the corresponding assembly / disassembly instructions.

Please remember that the compatibility of different product systems has not been tested unless specified otherwise in the product labeling. Consult Instructions for Use (www.ifu.stryker.com) for a complete list of potential adverse effects and adverse events, contraindications, warnings and precautions.

The surgeon must advise patients of surgical risks, and make them aware of adverse effects and alternative treatments.

 **WARNING**

- The patient should be advised that the device cannot and does not replicate a normal healthy bone, that the device can break or become damaged as a result of strenuous activity or trauma and that the device has a finite expected service life.
- Removal or revision of the device may be required sometime in the future due to medical reasons.



MRI Safety Information

MRI safety information



A patient with the VariAx 2 Distal Radius implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Device name	VariAx 2 Distal Radius
Static magnetic field strength (T)	1.5 T and 3.0 T
Maximum spatial field gradient	30 T/m (3000 gauss/cm)
RF excitation	Circularly Polarized (CP)
RF transmit coil type	Integrated Whole Body Transmit Coil
Operating mode	Normal Operating Mode
Maximum whole-body SAR (W/kg)	2 W/kg (Normal Operating Mode)
Scan duration	2 W/kg whole-body average SAR for 15 minutes of continuous RF (a sequence or back to back series/scan without breaks) followed by a wait time of 15 minutes if this limit is reached, for the total scanning session duration of up to 1 hour (or 60 minutes).
MR image artifact	The presence of this implant produced an image artifact of approximately 32 mm from the VariAx 2 Distal Radius implant when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.
Additional instructions	<p>⚠ CAUTION</p> <p>The MRI safety information provided is based on testing which did not include supplementary devices. If there are supplementary devices (i.e. plates, screws, wires, etc.) present in proximity to the VariAx 2 Distal Radius implant, this could result in additional MRI effects and the information provided above may not apply.</p>

VariAx 2

Distal Radius Locking Plates

Contents

Indications and contraindications	5
Overview	6
Implants: plate and screw platform	6
Instrumentation	7
Operative technique.....	11
Dorsal DR plates.....	11
Lateral DR plates	14
Dorsal medial DR plates.....	17
Dorsal lateral DR plates.....	20
Volar rim DR plates	23
Wrist spanning plates.....	26
System Components.....	32

Indications and Contraindications

VariAx 2 Distal Radius Plating System

The devices are non-active implants intended to provide temporary stabilization for bones or bone fragments.

Indications for use

The VariAx 2 Distal Radius implants are indicated for the treatment of fractures, non-unions, malunions and deformities of the distal radius.

Contraindications

The licensed healthcare professional's education, training and professional judgment must be relied upon to choose the most appropriate device and treatment. They should warn patients about these contraindications and limitations when appropriate.


Conditions presenting an increased risk of failure include:

- Any active or suspected latent infection or marked local inflammation in or about the affected area.
 - Compromised vascularity that would inhibit adequate blood supply to the fracture or the operative site.
 - Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and/or fixation of the devices.
 - Material sensitivity, documented or suspected.
 - Patients having inadequate tissue coverage over the operative site.
- Implant utilization that would interfere with anatomical structures or physiological performance.
 - Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
 - Other medical or surgical conditions which would preclude the potential benefit of surgery.

Overview

Implants: plate and screw platform

VariAx 2 DR Dorsal and DR Fragment Specific Plates



Wrist Spanning Plates may be used to span complex wrist fractures in a similar way to an external fixator.



Color Coding and Screw/Peg Options

CAUTION


When final tightening of the locking screw occurs, take care not to over-torque the screw. Excessive torque may damage the locking mechanism, the screw and /or the screwdriver blade.

Locking and non-locking screws can be used in any round hole.

Locking screws are laser marked with a 'dot' and 'ring' marking on the screw head to differentiate them from non-locking screws.

Pre-Angled Distal Screw Holes

The distal screw holes are angled to give a predetermined screw pattern in the distal bone block.



When drilling at a 0 degree angle relative to the plate hole, the screw trajectories relative to the plate surface will be achieved.

CAUTION

K-Wire with Olive stop (56-40281) is to be used only in the screw holes.

Overview

Instrumentation

SmartLock Polyaxial / Compression Drill Guide


Allows for ± 15 degrees of angulation. A lip on the drill sleeve will engage and allow toggling in the hole. The range in which the drill guide toggles will create a 30-degree cone and every angle in this range will be a locking position.

This may allow the surgeon to aim where the screw / peg should be placed. Also, depending on the placement of the plate, there may be a need to angle a screw / peg out of the fracture line.

The 2.0mm drill guide for T8 Screws (703684) facilitates drilling a 2.0mm pilot hole for a 2.4 or 2.7mm T8 screw or a 2.0mm peg centrally for locking or non-locking screws.

Not using a drill guide may lead to drilling out of specified locking range and compromise the locking capabilities.

Ensure the drill guide is fully engaged in the hole and is aimed in the desired direction.



CAUTION

- First fully engage the drill guide in the hole and then aim the drill in the desired direction.
- Make sure to drill perpendicular to oblong holes.
- Only use non-locking bone screws in oblong holes.

Overview

Instrumentation

Drills & Drill Guides for Lagging


In addition to the standard Drills and Drill guides, a number of solutions are also available to perform a lag screw technique independently.

Dedicated overdrills for each screw size are available for overdrilling the near cortex when placing a lag screw independently. In addition to being marked with the actual drill diameter on the AO Coupling, these overdrills are also marked with a single color ring corresponding to the desired screw diameter. This marking matches the marking on the correct side of the lagging drill guide.


CAUTION

Always match the screw color with at least one of the color ring markings.

In order to insert a lag screw independently of a plate, the Independent Lag Screw Drill Guides (703688 for 2.4mm screws and 703884 for 2.7mm screws) should be used. First, the near cortex should be overdrilled using the side of the drill guide marked with a single color ring to create a gliding hole (Step #1).



Step 1



Step 2

The other side of the drill guide can then be used (marked with two color rings) by inserting the 'top-hat' end in the already drilled gliding hole and using the standard drill bit through it to drill through the second cortex (Step #2). This standard drill is scaled in order to evaluate the appropriate screw length. Upon screw insertion, this technique will serve to lag the far cortex towards the near cortex, hence applying compression.

CAUTION

Take care not to damage the plate hole when using the Independent Lag Screw Drill Guide.

Overview

Instrumentation

Modular Handle

VariAx 2 offers a modular handle system. This is composed of two handle grip sizes (medium and large) that can be interchanged with either a bi-directional ratcheting AO-Coupling insert or a standard AO-Coupling insert.

Both handle sizes are equipped with a spin-cap to allow insertion using a two-finger technique. In order to disengage the insert from the handle, push down on the button on the distal part of the handle and pull the insert away from the handle.



CAUTION

The inserts must be removed from the handles before cleaning.

The ratcheting insert can work in three modes: clockwise ratcheting, counterclockwise ratcheting or neutral. To switch between the different modes, simply twist the distal part of the insert to the desired driving direction.

Overview

Instrumentation

Depth Measurement Options

VariAx 2 offers various options to evaluate the screw length. All drills are scaled so that the surgeon may evaluate the screw length when using the drill through the dedicated drill guides. A SpeedGuide (703891 for 2.0 drill bit and 703888 for the Speed Guide Sleeve) is also offered that allows the surgeon to drill and measure the hole depth in one step with a single instrument. For further information on the SpeedGuide, please refer to the SpeedGuide Operative Technique. Lastly, a standard depth gauge (705170) may be used either independently or through a plate hole.



Scaled Drill and Drill Guide



SpeedGuide



Depth Gauge



Taps

Taps

2.4mm and 2.7mm taps (703900 for 2.4mm screws and 703889 for 2.7mm screws) are available in the system.

⚠ CAUTION

Although all screws are self-tapping, it is recommended to use a tap if excessive resistance is felt during insertion or if the bone is dense.

Joystick for plate positioning and temporary fixation

The joystick for T8 holes can be used in any VariAx circular hole to aid in plate positioning. Additionally, they can also be used to temporarily fix the plate to the bone by inserting a K-wire with a diameter up to 1.6mm through a joystick that is already engaged in the plate hole.

After inserting the joystick tip in the circular hole, turn the knob on the upper part of the joystick clockwise to fix it in the hole. To remove the joystick, simply remove any K-wire and turn the knob counter-clockwise to disengage the tip from the hole.



Operative technique

DR Dorsal Plates


Longitudinal incision is made just ulnar to Lister's tubercle at the distal radius region.

Dissection is performed down to the extensor retinaculum. The third compartment is opened and the extensor pollicis longus is displaced radially.

The second compartment wrist extensors are subperiosteally elevated radially and the fourth compartment is subperiosteally elevated ulnarly. The terminal branches of the posterior interosseous nerve may be excised for pain reduction.

The fracture is reduced. The use of an external traction device and/or K-wire for temporary fixation may be helpful.

If necessary, adapt the plate for correct anatomical position. Removal of Lister's Tubercle might be necessary.



NOTICE

The distal lateral aspect of the plate can be bent along the groove with bending pliers to better match patient anatomy.

WARNING

- In order to reduce the likelihood of a stress riser effect and avoid reducing the fatigue properties of the implant, care should be taken to only bend the plate in between the holes.
- Only moderate bending is recommended.
- Excessive plate bending may lead to failure of the plate or the locking mechanism and should be avoided.
- Do not re-bend plates.



WARNING


- The plate bending pliers are designed to be used only in circular holes
- Always attach the bending pliers to two adjacent holes to prevent deformation of the screw holes

Operative technique

DR Dorsal Plates

The plate should be placed slightly proximal to the distal edge of the distal radius to avoid inserting screws /pegs into the joint.

The fixed angled screw trajectory provides support to the subchondral surface for a more stable construct. The distal portion of the plate is partially thickened for better screw head protection when screws are angled within the 15 degree variable angle arc.



The first pilot hole should be drilled in the oblong gliding hole.

If required, compression of the fracture site may be achieved by pulling the plate proximally with an eccentrically placed screw in the oblong hole.



Operative technique

DR Dorsal Plates

Measure the depth of the hole to determine screw length.

The screw is placed in the oblong gliding hole but not completely tightened to allow adjustment of the plate in a distal or proximal direction.

Confirm proper plate positioning by use of fluoroscopy and then tighten the first screw.

Repeat drilling, measuring, and placing of locking or non-locking screws/locking pegs into the distal holes of the plate. The position and number of screws applied depends on the type of fracture.

Place locking or non-locking screws in the proximal end of the plate.



⚠ CAUTION

Verify proper placement of screws and pegs by use of fluoroscopy to ensure that they do not penetrate the joint and are of appropriate length.

Close the incision.

⚠ CAUTION

Screw length may need to be changed after plate is fully seated on bone.

Operative technique

Lateral DR Plates

An incision is made along the radial column.

WARNING

Care must be taken to avoid injury to the superficial branch of the radial nerve.

The first dorsal compartment is released. The tendons are retracted volarly.

The brachioradialis can be elevated or the plate can be placed directly on the brachioradialis insertion.


The fracture should be reduced and stabilized with a K-wire placed from the distal radial styloid and aimed dorsal and proximal.

The plate may be slipped over the K-wire (07-40281) and placed along the radial column.

NOTICE

The plate should sit along the radial edge of the radius.

Confirm correct plate placement by use of fluoroscopy.



Operative technique

Lateral DR Plates

A non-locking screw placed in the oblong hole will compress the plate to the bone. If required, compression of the fracture site may be achieved by pulling the plate proximally with an eccentrically placed screw in the oblong hole.



A non-locking screw may also be placed in the oblong shaft-hole, instead of using the K-wire, for preliminary fixation of the shaft portion of plate.

Fill the remaining screw holes with either locking or non-locking screws, as necessary.

K-wires may also be used in conjunction with the distal screws.



The 3 in 1 K-wire bender/cutter/ inserter is used to bend K-wires distally.

It is recommended only one K-wire be placed distally at a time in order to make proper use of the K-wire Bending Pliers (64-20118).



Operative technique

Lateral DR Plates

After insertion, the tamp (64-00011) and mallet (43-09830) can be used to further insert the K-wires.



Ensure the non-locking screw in the oblong shaft hole is fully tightened.

The incision is closed.




Operative technique

Dorsal Medial DR Plates

When plate is being used in conjunction with the lateral plate, incision is made between the 4th and 5th extensor compartments.

WARNING

Care is taken to protect the superficial radial and dorsal ulnar nerve branches.



Incise the extensor retinaculum between the 4th and 5th compartments and subperiosteally expose the dorsal ulnar portion of the radius.



The fracture should be reduced and stabilized with traction and palmar flexion of the wrist or direct pressure. Ensure dorsal and ulnar fragments are reduced with K-wire if necessary.

The K-wire may be placed through the ulnar corner fragment.

The plate is slipped over the K-wire (07-40281) and placed along the dorsal ulnar corner.


Confirm correct plate placement by use of fluoroscopy.




Operative technique

Dorsal Medial DR Plates

A non-locking screw placed in the proximal oblong hole will compress the plate to the bone. If required, compression of the fracture site may be achieved by pulling the plate proximally with an eccentrically placed screw in the oblong hole.



Fill the remaining screw holes, distally to proximally, with either locking or non-locking screws, as necessary.



K-wires may also be used in conjunction with the distal screws.



Operative technique

Dorsal Medial DR Plates

The 3-in-1 K-wire bender/cutter/insertor is used to bend K-wires distally.

It is recommended only one K-wire be placed distally at a time in order to make proper use of the K-wire bending pliers (64-20118).



Ensure the non-locking screw in the oblong shaft hole is fully tightened.

The incision is closed.



Operative technique

Dorsal Lateral DR Plates

Longitudinal incision is made just ulnar to Lister's tubercle at the distal radius region.

Dissection is performed down to the extensor retinaculum. The third compartment is opened and the extensor pollicis longus is displaced radially.

The second compartment wrist extensors are subperiosteally elevated radially and the fourth compartment is subperiosteally elevated ulnarly. The terminal branches of the posterior interosseous nerve may be excised for pain reduction.

The fracture should be reduced and stabilized with traction and palmar flexion of the wrist. Ensure dorsal and radial fragments are reduced with K-wire if necessary. The K-wire may be placed through the radial corner fragment.

The plate is slipped over the K-wire (07-40281) and placed along the dorsal radial corner.




Operative technique

Dorsal Lateral DR Plates

A non-locking screw placed in the proximal oblong hole will compress the plate to the bone. If required, compression of the fracture site may be achieved by pulling the plate proximally with an eccentrically placed screw in the oblong hole.

Fill the remaining screw holes, distally to proximally, with either locking or non-locking screws.

K-wires may also be used in conjunction with the distal screws.



The 3-in-1 K-wire bender/cutter/insertor is used to bend K-wires distally.

It is recommended only one K-wire be placed distally at a time in order to make proper use of the K-wire Bending Pliers (64-20118).



Operative technique

Dorsal Lateral DR Plates

After insertion, the tamp (64-00011) and mallet (43-09830) can be used to further insert the K-wires.

Ensure the non-locking screw in the oblong shaft hole is fully tightened.


The incision is closed.



Operative technique

Volar Rim DR Plates

Incision is made through a standard volar Henry approach. Alternatively, a direct medial approach may provide superior exposure. For the direct medial approach, an extended carpal tunnel approach can be used distally and incorporated proximally to form an interval between the flexor tendons and the ulnar neurovascular bundle.



The fracture should be reduced and stabilized. The volar capsule should remain in place. Ensure the volar fragments are reduced and stabilized with a K-wire if necessary. The K-wire should be placed at the distal volar rim of the lunate facet.



NOTICE

For the Volar Rim DR Plate the K-wires are used for temporary fixation.

The plate is slipped over the K-wire (07-40281) and placed above volar capsule on the volar ulnar corner.

The vertical K-wire slots allow for distal or proximal adjustment and the horizontal K-wire slots allow for lateral or medial adjustment of the plate.



Operative technique

Volar Rim DR Plates

Ensure distal hooks penetrate the volar capsule and stabilize volar rim fragments.



With the plate in proper position, a non-locking screw is placed in the proximal oblong hole, compressing the plate to the bone.

If required, compression of the fracture site may be achieved by pulling the plate proximally with an eccentrically placed screw in the oblong hole.



CAUTION

Verify proper placement of screws and pegs by use of fluoroscopy to ensure that they do not penetrate the joint.

After stabilizing the distal fragments, fill the remaining screw holes with either locking or non-locking screws, as necessary.

The distal K-wire holes can be used for placing sutures in the volar capsule. #2-0, #3-0 XBraid S Sutures are suitable for use in the Volar Rim plate.

Temporary K-wires can now be removed.



Operative technique

Volar Rim DR Plates

Ensure the non-locking screw in the oblong shaft hole is fully tightened.

The incision is closed.




Operative technique

VariAx 2 Wrist Spanning Plates

Patient Preparation

The patient is placed in a supine position with the affected arm on a radiolucent hand fracture table.



Initial fracture reduction may be preformed by applying traction. Traction may be applied manually or with the help of a traction tower to restore the distal radius to proper length. Please see information on the Stryker Hand Traction System (07-30950).




Plate Positioning

Choose the appropriate length plate and place on the skin over the radius and over the 2nd or 3rd metacarpal until desired positioning is achieved. The wrist spanning plate is available in lengths standard (L172mm) and long (L198mm). The long plate extends more proximally and provides additional support in the radial shaft. Confirm placement with fluoroscopy. The desired placement may be marked with a marking pen.



Operative technique


VariAx 2 Wrist Spanning Plates

Exposure

An initial incision is made over the 2nd or 3rd metacarpal. If the incision is being made over the 2nd metacarpal, continue the dissection until the Extensor Carpi Radialis Longus (ECRL) tendon is identified.

! WARNING


Ensure superficial radial nerve branches are protected as necessary.



Reduce the distal radius fracture with application of K-wires or other supplemental fixation if necessary.

It is recommended to make a 2nd incision over the fractured distal radius and over the carpal bones in order to obtain reduction of articular fragments and in order to control the carpus with the center holes of the plate. This second incision should start just proximal to Lister's tubercle. The second and third dorsal compartments (ECRL and ECRB as well as EPL) should be released as necessary and minimally if possible. The VariAx 2 Wrist Spanning plate will need to be placed under the EPL.

If the plate is to be placed on the 3rd metacarpal, subperiosteal elevation of the fourth compartment is necessary to allow the plate to be placed along the dorsal aspect of the intermediate column of the distal radius. Bone graft may be used if necessary.



Operative technique

VariAx 2 Wrist Spanning Plates

The 3rd incision is made at the midline of the radius just proximal to the muscles of the abductor pollicis longus (APL) and extensor pollicis brevis (EPB) to expose the radial shaft. Incision point may need to be adjusted based on patient anatomy and fracture configuration.

WARNING


If too radial, it is possible that the superficial radial nerve may be encountered at the junction of the ECRL and the brachioradialis (BR).

Plate Bending

Bending of the plate should be avoided. However, if a slight dorsiflexion angle of the wrist is required, the wrist spanning plate may only be bent in-between the distal radius and carpal holes with the provided bending irons.

WARNING

- In order to reduce the likelihood of a stress riser effect and avoid reducing the fatigue properties of the implant, care should be taken to only bend the plate in between holes.
- Only moderate bending is recommended.
- Excessive plate bending may lead to failure of the plate or the locking mechanism and should be avoided.
- Do not re-bend plates.



Operative technique

VariAx 2 Wrist Spanning Plates

Plate Placement

Insert the plate distally to proximally under the 2nd or 4th extensor tendon compartment until the proximal end of the plate reaches the incision at the radial shaft. Blunt elevators may be utilized (705296) to create a path for the plate.

Check to ensure the wrist spanning plate is inserted with the correct orientation, with the oblong hole placed proximally.

The T8 Joystick (703927) can be used to aid in plate positioning.



Check the orientation of the wrist in neutral rotation. Check that the placement of the most distal screw holes are over the 2nd or 3rd metacarpal, and that the 4 central holes are over the carpal bones and distal radius. Use fluoroscopy as necessary.



Operative technique

VariAx 2 Wrist Spanning Plates

Check that traction is still adequately applied by either the Stryker Hand Traction System or surgical assistant.

The plate can be temporarily fixated to the bone by inserting a K-wire with a diameter up to 1.6mm through a Joystick that is already engaged in the plate hole, or by inserting an Olive Stop K-wire through any of the screw holes.

Once plate placement is confirmed, the first pilot holes should be drilled in the proximal oblong gliding hole.



Operative technique

VariAx 2 Wrist Spanning Plates

The appropriate sized non-locking screw is placed in the proximal oblong hole in the shaft of the plate and in one of the distal locking holes. Traction application to the wrist should be maintained. Do not fully tighten screw. After alignment is satisfactory, the screws are fully tightened.

A fracture reduction clamp (703822 or 702932) may be used to control the plate's position on the radial shaft and distal radius.

Fill the remaining metacarpal and carpal screw holes with either locking or non-locking screws as necessary.

Fill the remaining distal radius and radial shaft screw holes with either locking or non-locking screws as necessary.

Remove any K-wires.

Supplemental fixation can be utilized and left in based on preference.

Close all incisions.

It is recommended to remove the Wrist Spanning plate after healing is achieved.

Removal of the other VariAx 2 Distal Radius Plates is not required in general. The additional surgical trauma and the risks associated with additional anesthesia should be individually outweighed against the potential benefits for every patient. In the case of implant removal, the scar of the previous incision is (partly) re-opened and the screws and the plate are successively removed.




System Components

VariAx 2 Dedicated Wrist Tray


Top Level & 4-Level Tray Frame

consisting of:



Ref #	Description
940213	Upper Tray Wrist
1500-0006	Lower Part, 2 Level, Detachable
940347	Wrist Tray Top Layer Clip


Instruments



Ref #	Description
703921	Handle, Medium
703920	Handle, Large
703923	Handle Insert, AO, cannulated
703922	Handle Insert, AO, Ratchet, cannulated
703885	Depth Gauge for Distal Radius
703664	Screwdriver blade T8 AO
45-80040	Countersink For Screws 02.7/3.5mm AO Fitting
703663	Screwdriver Blade, AO, T8, self retaining
703675	Universal Holding sleeve
703927	Joystick for T8 screw holes
703684	Drill Guide, 2.0mm Drill, Comp/Polyaxial (T8)
703884	Drill Guide, For 2.7mm independent lag screw (T8)
703688	Drill Guide, For 2.4mm independent lag screw (T8)
703902	Drill Guide, 1.1 K-wire, Fixed Angled T8 2.0mm Drill
703896	Drill Bit, AO, Ø2.0mm x 135mm, Scaled *new
703897	Overdrill, AO, Ø2.7mm x 122mm*new
703696	Overdrill, AO, Ø2.4mm x 122mm
703899	Tap, AO, For 2.7mm Screws
703900	Tap, 2.4mm for AO screws
703888	SpeedGuide™, For 2.4/2.7 mm Screws, T8 (L = 30mm)*new
703891	SpeedGuide™ Drill, AO, Ø2.0mm (L = 30mm)

• All non-sterile plates may be ordered sterile by placing an "S" at the end of the REF Number.


VariAx 2 Dedicated Wrist Tray Drawer 1

 2.4mm and 2.7mm screws
& volar plates consisting of:




Ref #	Description
1500-0005	Drawer For Modules & Screw Racks
940348	Wrist Tray Screw Drawer Clip
940234	Screw Rack for 2.7mm Screws, T8
940235	Screw Rack for 2.4mm Screws, T8

2.4mm Locking Screw




Titanium Ref #	Length mm
656008	8
656010	10
656012	12
656014	14
656016	16
656018	18
656020	20
656022	22
656024	24
656026	26
656028	28
656030	30
656032	32
656034	34
656036	36
656038	38

2.4 Non Locking Screws



Titanium Ref #	Length mm
656108	8
656110	10
656112	12
656114	14
656116	16
656118	18
656120	20
656122	22
656124	24
656126	26
656128	28
656130	30
656132	32
656134	34
656136	36
656138	38

2.0 mm Locking Pegs



Titanium Ref #	Length mm
656616	16
656618	18
656620	20
656622	22
656624	24
656626	26

VariAx 2 Dedicated Wrist Tray

2.7mm Locking Screws



Titanium Ref #	Length mm
656308	8
656310	10
656312	12
656314	14
656316	16
656318	18
656320	20
656322	22
656324	24
656326	26
656328	28
656330	30
656332	32
656334	34
656336	36
656338	38
656340	40
656345	45
656350	50

2.7mm Non Locking Screws



Titanium Ref #	Length mm
656408	8
656410	10
656412	12
656414	14
656416	16
656418	18
656420	20
656422	22
656424	24
656426	26
656428	28
656430	30
656432	32
656434	34
656436	36
656438	38
656440	40
656445	45
656450	50

2.7mm Non Locking Partially Threaded Screws



Titanium Ref #	Length mm
656516	16
656518	18
656520	20
656522	22
656524	24
656526	26

Washer




Ref #	Description
619920	Washer

VariAx 2 Dedicated Wrist Tray

Drawer 2

Tray Content: Anatomical Dorsal, DR Fragment Specific, Wrist Spanning, Distal Ulna plates and instruments



Ref #	Description
1500-0005	Drawer For Modules & Screw Racks
940459	Wrist Tray Clip DR DU & Spanning Plates
940458	Wrist Add-on Clip DR & DU Spanning Plates
940199	Insert DU Plates & Frag Spec Instruments
940198	Insert DR & Wrist Spanning Plates

Anatomical Dorsal Plates



















Ref #	Description	Length mm	Profile Height mm
625090	Dorsal DR Plate, Narrow, Short, Right	47	2
625091	Dorsal DR Plate, Narrow, Short, Left	47	2
625094	Dorsal DR Plate, Narrow, Long, Right	58	2
625095	Dorsal DR Plate, Narrow, Long, Left	58	2
625098S*	Dorsal DR Plate, Narrow, X-Long, Right	69	2
625099S*	Dorsal DR Plate, Narrow, X-Long, Left	69	2
625092	Dorsal DR Plate, Wide, Short, Right	49	2
625093	Dorsal DR Plate, Wide, Short, Left	49	2
625096	Dorsal DR Plate, Wide, Long, Right	60	2
625097	Dorsal DR Plate, Wide, Long, Left	60	2
625100S*	Dorsal DR Plate, Wide, X-Long, Right	71	2
625101S*	Dorsal DR Plate, Wide, X-Long, Left	71	2

• All non-sterile plates may be ordered sterile by placing an "S" at the end of the REF Number.

* Available sterile only.

VariAx 2 Dedicated Wrist Tray

Anatomical Fragment Specific Plates

	Ref #	Description	Length mm	Profile Height mm
	625120	Lateral DR Plate, Short, Right	43	2
	625121	Lateral DR Plate, Short, Left	43	2
	625122	Lateral DR Plate, Long, Right	56	2
	625123	Lateral DR Plate, Long, Left	56	2
	625140	Dorsal Lateral DR Plate, Short, Right	39	2
	625141	Dorsal Lateral DR Plate, Short, Left	39	2
	625142	Dorsal Lateral DR Plate, Long, Right	49	2
	625143	Dorsal Lateral DR Plate, Long, Left	49	2
	625130	Dorsal Medial DR Plate, Short, Right	39	2
	625131	Dorsal Medial DR Plate, Short, Left	39	2
	625132	Dorsal Medial DR Plate, Long, Right	49	2
	625133	Dorsal Medial DR Plate, Long, Left	49	2
	625150	Volar Rim DR Plate, Short, Right	32	1.6
	625151	Volar Rim DR Plate, Short, Left	32	1.6
	625152	Volar Rim DR Plate, Long, Right	43	1.6
	625153	Volar Rim DR Plate, Long, Left	43	1.6

Wrist Spanning Plates

	Ref #	Description	Length mm	Profile Height mm
	625106	Wrist Spanning Plate, Standard	172	3
	625108	Wrist Spanning Plate, Long	198	3

• All non-sterile plates may be ordered sterile by placing an "S" at the end of the REF Number.

VariAx 2 Dedicated Wrist Tray

Anatomical Dorsal Plates Trials



Ref #	Description	Length mm
705890	Dorsal DR Plate Trial, Narrow Right Short	47
705891	Dorsal DR Plate Trial, Narrow, Left, Short	47
705896	Dorsal DR Plate Trial, Wide, Right, Long	60
705897	Dorsal DR Plate Trial, Wide, Left, Long	60

Anatomical Fragment Specific Plates Trials



Ref #	Description	Length mm
705920	Lateral DR Plate Trial, Right, Short	43
705921	Lateral DR Plate Trial, Left, Short	43
705930	Dorsal Med. DR Plate Trial, Right, Short	39
705931	Dorsal Med. DR Plate Trial, Left, Short	39
705940	Dorsal Lat. DR Plate Trial, Right, Short	39
705941	Dorsal Lat. DR Plate Trial, Left, Short	39
705950	Volar Rim DR Plate Trial, Right, Short	32
705951	Volar Rim DR Plate Trial, Left, Short	32

Wrist Spanning Plates Trial




Ref #	Description	Length mm
705906	Wrist Spanning Plate Trial, Standard	172

VariAx 2 Dedicated Wrist Tray

Drawer 3

Reduction Instruments, consisting of:



Ref #	Description
1500-0005	Drawer For Modules & Screw Racks
940333	Labeling Clip Core Tray Reduction
940250	Reduction Instruments 1 Insert
940251	Reduction Instruments 2 Insert
940252	Reduction Instruments 3 Insert
940350	Wrist Tray Reduction Instrument Clip


Reduction Instruments

Ref #	Description
705297	Straight Reduction Clamp, Broad
703822	Straight Reduction Clamp
702926	Repositioning forceps L130mm
702932	Repositioning forceps L143mm (Lobster Claw)
705294	Periosteal Elevator, Round Edge 6mm
705295	Periosteal Elevator, Flat Blade 13mm
705293	Periosteal Elevator, Straight Edge 6mm (optional)
705296	Periosteal Elevator, Curved Blade 13mm (optional)
45-80010	Plate Bending Pliers
700151	Hook
700664	6mm Hohmann Retractor
700665	8mm Hohmann Retractor
700667	15mm Hohmann Retractor
703938	Bending Iron for VariAx Plates
700153	Ballspike
703818	K-Wire with Stop 2mm (pack of 5)
390192	K-Wire 2.0mm, x 150mm (pack of 10)
390164	K-Wire 1.6 x 150mm (pack of 10)
390157	K-Wire 1.25mm, x 150mm (pack of 10)

VariAx 2 Wrist Add-On Tray

Top Level & 2-Level Tray Frame


consisting of:



Ref #	Description
940221	Upper Tray Wrist Add-on, 2 Level

Drawer 1













Tray Content: Anatomical Dorsal, DR Fragment Specific, Wrist Spanning, Distal Ulna plates and instruments



Ref #	Description
1500-0005	Drawer For Modules & Screw Racks
940459	Wrist Tray Clip DR DU & Spanning Plates
940458	Wrist Add-on Clip DR & DU Spanning Plates
940199	Insert DU Plates & Frag Spec Instruments
940198	Insert DR & Wrist Spanning Plates

VariAx 2 Wrist Add-On Tray

Anatomical Dorsal Plates

















	Ref #	Description	Length mm	Profile Height mm
	625090	Dorsal DR Plate, Narrow, Short, Right	47	2
	625091	Dorsal DR Plate, Narrow, Short, Left	47	2
	625094	Dorsal DR Plate, Narrow, Long, Right	58	2
	625095	Dorsal DR Plate, Narrow, Long, Left	58	2
	625098S*	Dorsal DR Plate, Narrow, X-Long, Right	69	2
	625099S*	Dorsal DR Plate, Narrow, X-Long, Left	69	2
	625092	Dorsal DR Plate, Wide, Short, Right	49	2
	625093	Dorsal DR Plate, Wide, Short, Left	49	2
	625096	Dorsal DR Plate, Wide, Long, Right	60	2
	625097	Dorsal DR Plate, Wide, Long, Left	60	2
	625100S*	Dorsal DR Plate, Wide, X-Long, Right	71	2
	625101S*	Dorsal DR Plate, Wide, X-Long, Left	71	2

• All non-sterile plates may be ordered sterile by placing an "S" at the end of the REF Number.

* Available sterile only.

VariAx 2 Wrist Add-On Tray






Anatomical Fragment Specific Plates

	Ref #	Description	Length mm	Profile Height mm
	625120	Lateral DR Plate, Short, Right	43	2
	625121	Lateral DR Plate, Short, Left	43	2
	625122	Lateral DR Plate, Long, Right	56	2
	625123	Lateral DR Plate, Long, Left	56	2
	625140	Dorsal Lateral DR Plate, Short, Right	39	2
	625141	Dorsal Lateral DR Plate, Short, Left	39	2
	625142	Dorsal Lateral DR Plate, Long, Right	49	2
	625143	Dorsal Lateral DR Plate, Long, Left	49	2
	625130	Dorsal Medial DR Plate, Short, Right	39	2
	625131	Dorsal Medial DR Plate, Short, Left	39	2
	625132	Dorsal Medial DR Plate, Long, Right	49	2
	625133	Dorsal Medial DR Plate, Long, Left	49	2
	625150	Volar Rim DR Plate, Short, Right	32	1.6
	625151	Volar Rim DR Plate, Short, Left	32	1.6
	625152	Volar Rim DR Plate, Long, Right	43	1.6
	625153	Volar Rim DR Plate, Long, Left	43	1.6

Wrist Spanning Plates

	Ref #	Description	Length mm	Profile Height mm
	625106	Wrist Spanning Plate, Standard	172	3
	625108	Wrist Spanning Plate, Long	198	3





Fragment Specific Instruments

	Ref #	Description
	43-09830	Mallet, 250g
	64-00011	Tamp
	07-40281	Zebra K-Wire 1.1mm (pack of 10)
	56-40281	K-Wire with Olive Stop, 1.4mm (pack of 5)
	64-20118	K-Wire Bending Pliers









• All non-sterile plates may be ordered sterile by placing an "S" at the end of the REF Number.

VariAx 2 Wrist Add-On Tray


Anatomical Dorsal Plates Trials

	Ref #	Description	Length mm
	705890	Dorsal DR Plate Trial, Narrow Right Short	47
	705891	Dorsal DR Plate Trial, Narrow, Left, Short	47
	705896	Dorsal DR Plate Trial, Wide, Right, Long	60
	705897	Dorsal DR Plate Trial, Wide, Left, Long	60

Anatomical Fragment Specific Plates Trials

	Ref #	Description	Length mm
	705920	Lateral DR Plate Trial, Right, Short	43
	705921	Lateral DR Plate Trial, Left, Short	43
	705930	Dorsal Med. DR Plate Trial, Right, Short	39
	705931	Dorsal Med. DR Plate Trial, Left, Short	39
	705940	Dorsal Lat. DR Plate Trial, Right, Short	39
	705941	Dorsal Lat. DR Plate Trial, Left, Short	39
	705950	Volar Rim DR Plate Trial, Right, Short	32
	705951	Volar Rim DR Plate Trial, Left, Short	32

Wrist Spanning Plates Trial

	Ref #	Description	Length mm
	705906	Wrist Spanning Plate Trial, Standard	172

Notes

Notes

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the product label and/or Instructions for Use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The Instructions for Use, Operative Techniques, Cleaning instructions, patient information leaflets and other associated labeling may be requested online at www.ifu.stryker.com or www.stryker.com. If saving the Instructions for Use, Operative Techniques, Cleaning instructions from the above mentioned websites, please make sure you always have the most up to date version prior to use.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: VariAx, Stryker. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked.

This document is not applicable to US and Canada.

Content ID: VAX-ST-72, Rev. 1, 02-2024

Copyright © 2024 Stryker

  **Manufacturer:**
Stryker GmbH
Bohnackerweg 1
2545 Selzach,
Switzerland
www.stryker.com