

BariMatt™ Plus

Pressure Redistribution Bariatric Surface



Multiple surface configurations provide the flexibility to care for a wide range of patients

Customize surface width and length with foam bolster inserts

Seamless cover minimizes infection control issues

Enhance patient comfort with high-memory visco-elastic foam

Versatile Configurations for Every Patient

The BariMatt Plus provides care for a wide range of patients. The versatile surface configurations allow for the surface to be customized to each individual patient with foam bolster inserts. The seamless cover minimizes infection control issues while the visco-elastic foam maximizes patient comfort.



Customizable width and length provides optimal comfort for each patient.

Rental Solutions Delivered by



1 800 STRYKER Option 3

Universal | Hospital | Services®

BariMatt Plus

stryker®



Medical

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Foam bolsters inserts
- Visco-elastic foam
- Seamless nylon cover
- Four size conversions:
 - 36" x 80"
 - 42" x 80"
 - 48" x 80"
 - 48" x 88"
- Seamless sleeves for insertion of bolsters to main surface

Specifications

• Model Number	2235
• Overall Length	
Retracted	80" (203.2 cm)
Expanded	88" (223.5 cm)
• Overall Width	
Retracted (both foam bolsters removed)	36" (91.4 cm)
One side bolster inserted and one bolster removed	42" (106.7 cm)
Expanded (both foam bolsters inserted)	48" (121.9 cm)
• Thickness	6" (15.2 cm)
• Weight Capacity	1000 lbs (453.6 kg)

*Stryker reserves the right to change specifications without notice.
Support surfaces are manufactured for Stryker Medical.*



Mixed Sources
Product group from well-managed
forests, controlled sources and
recycled wood or fibre
Cert no. SW-COC-1530
www.fsc.org
© 1996 Forest Stewardship Council

Copyright © 2008 Stryker
PC/SCG-7.5M-0508, 07-SMD-2184
Mkt Lit-215-122107 Rev B