

# Bedside Stands

Durable yet elegant, Stryker bedside stands are designed to meet the demanding needs of today's clinical environment. With five unique styles and a variety of laminate and hardware options, these bedside stands provide design latitude to coordinate your patient suites.



Contemporary



Traditional



Shaker



Legend™



Premier™

**Concealed casters** enhance aesthetics and help protect casters from breaking

**Removable drawers** and plastic drawer liners are extremely durable and make the stands easy to clean

**Steel drawer frames** and hinges come with a lifetime warranty

**Specialized thermoform** process created a seam-free finish on top and front panels

**Vented back** allows moisture and odors within to air out

**Raised edges** on top prevents spills and facilitates cleaning

# Bedside Stands

**stryker**



Contemporary



Traditional



Shaker



Legend™



Premier™

## Contemporary standard features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Stainless steel pull handles

## Traditional standard features

- Concealed 1.5" (4 cm) twin-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Raised drawer and door panels

## Shaker standard features

- Concealed 1.5" (4 cm) twin-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Thermoform drawer fronts and top
- Round brushed stainless or antique brass knobs
- Recessed drawer and door panels

## Legend standard features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Plastic drawers and drawer frames
- Polished chrome or brass pull handles
- Solid T-band edgemold
- Vented back

## Premier standard features

- Concealed 1.5" (4 cm) dual-wheel casters
- Finished edges
- Removable drawers and plastic drawer liners
- Vented back
- Steel drawer frames and glides
- Polished chrome or brass pull handles
- Raised spill containment top

## Optional features

- Three-drawer or drawer/door configuration
- Spill containment top (excluding Legend)
- Continuous pull handles (Premier and Contemporary models only)
- Key or slam lock on top drawer
- Choose from six standard laminates

Specifications	Contemporary	Traditional	Shaker	Legend	Premier
<b>Part #</b>	Part #4400-800/803-000	4400-880/883-000	4400-840/843-000	4400-500/503-000	4400-600/603-000
<b>Height</b>	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)	30.5" (77 cm)
<b>Width</b>	20" (51 cm)	19.25" (49 cm)	19.25 (49 cm)	18.5" (47 cm)	19.75" (50 cm)
<b>Weight</b>	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)	80 lbs (36 kg)
<b>Depth</b>	20" (51 cm)	18.5" (47 cm)	18.5" (47 cm)	18" (46 cm)	19.75" (50 cm)

- BIFMA modified static load testing up to 250 lbs.

## Custom capabilities

- Custom capabilities available - please contact your Stryker sales representative

## Warranty

- One year parts and labor
- Lifetime warranty on steel drawer frames and hinges

Mkt lit-96-02022005 Rev F  
Copyright © 2018 Stryker  
Printed in U.S.A.

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

**stryker.com**

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Legend, Premier, Stryker. All other trademarks are trademarks of their respective owners or holder.