

Emergency Relief Bed™

limited-release



In response to the urgent and growing need for patient platforms during the current COVID-19 pandemic, we have developed our limited-release Emergency Relief Bed which is designed to help emergency responders and caregivers move and position patients efficiently.

Our Emergency Relief Bed is a low-cost, readily available solution intended to quickly serve those on the front lines of the COVID-19 crisis – ranging from triage centers and pop-up areas of care, to the emergency department, and all in-patient hospital general departments.

It features a 30-degree manually adjustable backrest angle to accommodate patients in respiratory distress, including vented patients. It also features low height, a 10 cm thick mattress, and an attached IV pole.

Maximum patient weight
160 kg

30° Manually adjustable backrest
accommodates vented patients

10 cm Thick surface
included

Part number: 7600-000-900

Emergency Relief Bed specifications

- Frame dimensions: 216 cm x 87 cm x 30 cm (85" x 34.4" x 11.8")
- 30 cm low height (without mattress)
- 30° manually adjustable backrest
- 5 backrest angles
- Wheeled frame
- Caster lock brakes
- Safe Working Load: 180kg
- IV pole
- Patient weight capacity: 160 kg
- 10 cm thick foam mattress included

