

# Eye Surgery Stretcher and Surgistool Chair



## Eye Surgery Stretcher

### Access. Simple. Functionality.

**Stryker's Eye Surgery Stretcher** allows head-end clearance, giving surgeons ample leg room to manoeuvre while performing delicate eye procedures. The steady work surface and mobile design support in a efficient work-flow of surgical throughput.



### Enhanced head-end clearance

Allows the surgeon's access for either a temporal or a superior approach to surgery.



**Optional contoured or flathead piece** comfortably supports patient during procedures.

**Retractable fifth wheel** helps improve mobility.

**Metal tooth-and-groove braking system** prevents swivel and roll.

**Option of two- or three-position fold-down siderails.**

**Multiple bumper colour options** allow you to customise the stretcher for identification purposes.



Navy Grey Red

## Surgistool Chair

### The ideal surgeon's stool

**Surgistool Chair** offers customisable seat configurations, while fully adjustable backrest and armrests provide optimal support to meet the surgeon's specific positioning preferences. Easy pedal access and a stable, yet mobile, design helps to facilitate clinical performance.



**Three spindle base** provides stable control and pedal access.



**Optional padded armrests** allows versatile positioning and stability.

**Multiple seat configurations** offer personalised comfort.



Teardrop seat



Wedge seat



Saddle seat (optional)

# Eye Surgery Stretcher and Surgistool Chair



## Standard features: Eye Surgery Stretcher

- 8 cm (3") Enhanced Comfort Mattress
- 20 cm (8") Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

## Standard features: Surgistool Chair

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

## Specifications

	Eye Surgery Stretcher	Surgistool Chair
<b>Model number</b>	1089	830
<b>Overall length</b>	228.6 cm (90")	
<b>Overall width</b>	80 cm (31.5")	
<b>Weight capacity</b>	227 kg (500 lb.)	136 kg (300 lb.)
<b>Height range</b>		
High	86.4 cm (34") (to litter top)	72 cm (28.5") / 77 cm (30.5") (with riser kit)
Low	56.5 cm (22.25") (to litter top)	53 cm (21") / 58 cm (23") (with riser kit)
<b>Seat cushion style</b>		Teardrop / wedge / saddle
<b>Litter positioning</b>		
Backrest	0°– 90°	
Knee gatch	0°– 30°	
Trend./Reverse Trend.	± 18°	
<b>Patient surface</b>	66 x 221 cm (26" x 87")	
<b>Siderails</b>	33 x 139.5 cm (13" x 55")	
<b>Caster diameter</b>	20 cm (8")	

Stryker reserves the right to change specifications without notice.

## Optional features: Eye Surgery Stretcher

- |   |  |   |
|---|--|---|
| <ul style="list-style-type: none"> <li>• 10 cm (4") Ultra Comfort Mattress</li> <li>• Pre/post-op shoulder pads assembly (crank backrest only)</li> <li>• Air delivery system/ drape support</li> <li>• Crank backrest</li> <li>• Crank knee gatch</li> <li>• Defibrillator tray/ foot extender</li> <li>• Temporal or superior wristrest</li> <li>• Footboard/chartholder</li> </ul> | <ul style="list-style-type: none"> <li>• Three-sided hydraulic controls</li> <li>• Three-sided brake and steer pedals</li> <li>• Integrated patient transfer/ arm board</li> <li>• IV caddy</li> <li>• Restraint straps</li> <li>• Serving/instrument tray</li> <li>• Siderail pads</li> <li>• Surgery accessory rail</li> <li>• Three-position siderails</li> <li>• Upright oxygen bottle holder</li> </ul> | <ul style="list-style-type: none"> <li>• Oxygen bottle retaining ring</li> <li>• Serving/instrument tray holder</li> <li>• Head pads - Ultra Comfort               <ul style="list-style-type: none"> <li>• 8 cm (3"), flat</li> <li>• 8 cm (3"), concave</li> </ul> </li> <li>• IV poles               <ul style="list-style-type: none"> <li>• Permanent folding two- or three-stage at foot end</li> <li>• Removable</li> </ul> </li> <li>• Extended warranty available</li> </ul> |
|---|--|---|

## Optional features: Surgistool Chair

- 5 cm (2") riser kit
- Saddle seat
- Multi-position armrest

## Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Surgistool, Stryker. All other trademarks are trademarks of their respective owners or holders.

This material is not intended for distribution outside the EU and EFTA.



16208EN

DLS 02/2018

2017\_16208

Copyright © 2018 Stryker  
stryker.com