

GoBed® II Med/Surg Bed

Safe. Simple. Mobile.



The acceptable number of adverse events in your facility is zero. With the right technologies and processes in place, potential adverse events often can be predicted and prevented.

Extensive input and feedback from nursing professionals helped qualify the need for a med/surg bed that helps add efficiency and mobility to patient care practices. Numerous one-button or one-handed controls allow caregivers to spend their time helping patients recover quickly in a secure and comforting environment.



The GoBed II features full length siderail coverage to aid in patient safety.



A low bed height of 14.5 inches promotes proper patient positioning for safe entry and exit.



Conveniently located on the footboard, the integrated pump holder and 120-volt outlet help to maximize space at the foot end of the bed.

Available with Stryker's Chaperone® Bed Exit System and Zone Control® technology.

Single-button, in-bed scale provides repeatable weights with the patient in any position.

Complete siderail coverage and 14.5-inch low bed height contribute to a safe patient platform.

Centrally located fifth-wheel helps improve maneuverability and ease transport.

Integrated pump holder and optional 120-volt outlet conveniently located on the footboard.

GoBed II Med/Surg Bed



3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- 14.5-inch (37 cm) low-bed height
- Retractable fifth-wheel steering
- StayPut® retractable bed frame design
- 6-inch (15 cm) casters
- Trend./reverse Trend.
- Four independent electric motors
- Trend. display on footboard (only available with scale option)
- Centrally-located steer and four-wheel brake mechanism
- Four drainage bag hooks
- Eight IV pole/traction equipment sockets
- Roller bumpers
- Manual back-up for head and knee control
- Integrated pump holder
- Nurse controls on footboard and siderails
- Fixed patient controls on siderails
- Degree indicator for head elevation
- Patient restraint locations
- Electronic function lockout controls
- One-handed dampened siderail release
- Auto contour
- Photo-sensitive night light

GoBed II LX: All Standard Features Plus

- Chaperone Bed Exit System
- In-bed scale system
- CPR release

GoBed II EX: All Standard Features Plus

- Chaperone Bed Exit System
- In-bed scale system
- One-button cardiac chair
- CPR release
- Siderail communications package (includes nurse call with speakers, TV, radio, volume, room/read lights)

Specifications

Model Number	FL28C, FL28EX
Overall Length	94.25" (239.4 cm)
Overall Width	
Siderails Up	40" (101.6 cm)
Siderails Down	39" (99.1 cm)
Weight Capacity	500 lb. (227 kg)
Height Range (to litter top)	
High	29" (73.7 cm)
Low (standard with 6" casters)	14.5" (36.8 cm)
Litter Positioning	
Backrest	0°– 60° (siderails down)
Knee Gatch	0°– 28°
Trend./Reverse Trend.	±14°
Retraction	10" (25 cm)
Patient Surface	36" x 84" (91 x 213 cm)
Caster Diameter	6" (15 cm)

Electronics

Agency Approvals	UL60601-1, CSA601.1 and IEC60601-2-38
Volts	100 VAC, 120 VAC, 200 VAC, 220 VAC
Current Leakage	< 100 microamperes
Frequency	50– 60 Hz
Hospital Grade Plug	Yes

Stryker reserves the right to change specifications without notice.

Optional Features

- Chaperone Bed Exit System
- Chaperone with Zone Control
- In-bed scale system
- CPR release
- Smart TV (includes closed caption, channel up/down and mute)
- Support Surfaces: pressure ulcer prevention and management
- Bed extender with pad
- Pillow speaker interface
- Removable litter covers
- Upright oxygen bottle holder
- IV poles
 - Permanent
 - Removable
- 120-volt auxiliary outlet
- Patient helper system
- Emergency crank handle
- Footboard-mounted monitor tray