

# Gynnie OB/GYN Stretcher

stryker

Easy. Multi-functional. Convenient.

## Intuitive patient positioning

Gynnie is ideal for pelvic exams and general use, this multifunctional stretcher can outfit an entire emergency department. Gynnie can be configured with optional OB/GYN accessories for even greater functionality.

**10 cm thick Ultra Comfort mattress.**

**Retractable fifth wheel** helps improve mobility.

**Dual pedestal hydraulics** to help optimise stability.

**Permanently** attached parts to maximise caregivers' time.

**Four wheel, steel ring brakes** with dual side and head activators to help secure platform and accessibility.



## Pelvic exam conversion

Three easy steps to convert this stretcher into a pelvic exam with self-locking stirrups.



## Optional calf supports



## Glideaway siderails

Swings up and out for patient comfort and allows for zero transfer gap when lowered.



**Stryker Australia Pty Ltd**  
8 Herbert Street  
St Leonards NSW 2065  
61 2 9467 1000  
www.stryker.com.au

**Stryker New Zealand Limited**  
PO Box 17-136  
Greenlane, Auckland 1546, New Zealand  
64 9 573 1890  
www.stryker.com

JMEH\_1966\_2\_AOT  
Copyright © 2020 Stryker

# Gynnie OB/GYN Stretcher

## Standard features

- 10.2 cm (4") Ultra Comfort mattress with sealed seams
- 20.3 cm (8") Omni surface casters with wheel covers
- 66.0 cm (26") wide patient surface
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics
- Retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual-side and head activators
- Glideaway siderails
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Permanently-attached foot supports
- Pop-up push handles
- Roller bumpers
- Two IV receptacles
- Unique one-handed foot section release lever
- **Two year warranty**

## Optional features

- Hoop for drainage basin
- IV poles (permanent folding two- or three-stage or standard removable)
- Patient restraints
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder

## Specifications

|                                     |                           |
|-------------------------------------|---------------------------|
| <b>Model number</b>                 | 1061                      |
| <b>Overall length</b>               | 205.7 cm (81")            |
| <b>Overall width</b>                |                           |
| Siderails up                        | 80.0 cm (31.5")           |
| Siderails down                      | 71.1 cm (28")             |
| <b>Weight capacity</b>              | 226.8 kg (500 lb.)        |
| <b>Height range (to litter top)</b> |                           |
| High                                | 90.8 cm (35.75")          |
| Low                                 | 55.9 cm (22")             |
| <b>Litter positioning</b>           |                           |
| Backrest                            | 0° - 90°                  |
| Knee gatch                          | 0° - 90°                  |
| Trend./reverse Trend.               | ±18°                      |
| Maximum angle                       | 87.5°                     |
| <b>Patient surface</b>              | 66 x 192 cm (26" x 75.5") |
| <b>Siderails</b>                    | 33 x 99 cm (13" x 39")    |
| <b>Caster diameter</b>              | 20.3 cm (8")              |
| <b>Break away point from wall</b>   | 172.7 cm (68")            |

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local Stryker sales representative.