

stryker

# Improving outcomes in your ICU



**The ICU is a specialized environment,** staffed by highly trained caregivers, helping patients who need immediate interventions. At Stryker, we partner with caregivers to solve critical issues and ongoing problems, focusing on:

- **Early mobility**
- **Efficiency**
- **Safety and prevention**

# Altrix

## Precision temperature management system

### Prompted therapy initiation

The user-friendly graphic interface makes it easy to initiate patient therapy with just a few simple taps of the finger.

### Precise, rate controlled therapy

Constant patient and water temperature monitoring allows rapid responsiveness using thousands of power settings and steady state patient control  $\pm 0.1^{\circ}\text{C}$

Pending 510(k) clearance.  
Not available for sale.



### Individual water circuit monitoring


Simple, efficient, self-sealing hose ports eliminate water spillage and speed target temperature. All three water circuits can be individually monitored.

### Removable water reservoir

Therapy setup is simple and easy and allows for uninterrupted temperature management when adding or refilling the reservoir.

### Versatile design

Altrix is designed for versatility and speed in a safe, easily operated system that can be used in all clinical settings, for both traditional and critical applications.



**All-electric brakes**

Provide quiet operation and long-term reliability while minimizing your install costs.

**360° configurability**

Completely customizable to meet your specific service and equipment management needs.

**Designed for safety and efficiency**

Removes cables and hoses usually tethered to the wall and creates a safe, efficient work environment.

**Engineered for superior mobility**

Improves mobility within the room and provides direct access to your patient's airway during emergent events.

# Teletom

Equipment management systems



# InTouch

Critical care bed

**Zoom drive**

Motorized drive system provides mobility in the ICU and less stress on the caregiver.

**30-degree head of bed button**

Quickly adjust the bed to aid in VAP and patient safety protocols.

**One-touch electric brake**

With one touch, brake all wheels, reducing bending and stretching needed for conventional brakes.

**Integrated scale and Chaperone Bed Exit System with Zone Control**

Weigh with patient in any position, store 50 weights and track patient position with local or remote alarm and spoken patient commands.

**iBed Complete**

iBed Complete provides visual, audible and remote bed status information for protocol reminders, documentation, speaking translation and sound therapy.

# TruRize

## Clinical chair

### Caregiver specific control panel

Easily accessible, the controls are highly intuitive in a complex environment.

### Sit to stand

The Stand Assist function supports your patient to a standing position, promoting early mobility.

### Pivoting armrests

Simplify patient transfers with arms that move out of the way.

### Lighted occupant control panel

Your patient has a selection of comfort controls while full chair mobility is in your hands.

### Chair exit alarm

The integrated exit alarm system may help shorten response time to patients.

### Lateral rotation and turn assist

Isolibrium provides lateral rotation therapy up to 40 degrees to assist with both early mobility protocols and the prevention of pulmonary complications.

### Air pod technology

Patient specific pressure redistribution, envelopment and immersion can assist in pressure injury prevention.

### HeelGel technology

Pressure reduced on heels as HeelGel buckles to help prevent skin breakdown.


### Enhanced microclimate

Unique low-air-loss system channels air toward patient's skin.

# Isolibrium

## Support surface





At Stryker, we  
share your goals  
of **improving**  
**patient outcomes**  
in the ICU.



Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002 U.S.A.  
t: 269 329 2100  
toll free: 800 327 0770

[www.stryker.com](http://www.stryker.com)

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Altrix, Chaperone, HeelGel, iBed, iBed Complete, InTouch, Isolibrium, Stryker, Teletom, TruRize, Zoom, Zone Control**. All other trademarks are trademarks of their respective owners or holders.

Mkt Lit-1348 01 JUN 2016 Rev A.3  
MHD/PS 0616, SMD428

Copyright © 2016 Stryker  
Printed in U.S.A.