

LD304

Birthing Bed

Making Mum feel at home.



The LD304 birthing bed provides comfort and security for both the mother and caregiver. Simple operation, mobility and an ergonomic design are intended to provide security while positioning to Mum's comfort.

Today's labour and delivery needs have grown more challenging.

With increasing patient weights¹, routine epidurals, older mums and various birthing options, come more complex labour and delivery needs. The result? Greater physical demands on you.

Our BackSmart design offers a solution. With nearly infinite positioning and ergonomic push handles, the LD304 is easy to use and highly functional.



Convenient controls, such as single button Trendelenburg and our powered Skoocher feature, allow customised positioning with minimal effort.



An ergonomic design of our foot section allows caregivers of all heights easy visual access to the perineum.

Ergonomic design

provides comfort and support for both mother and caregiver.

Patented Lock-Rite foot section

glides off for quick and easy breakdown.

Length-adjustable footrests

provide multiple leg positioning capabilities.

Ultra Comfort mattress

with integrated, high-memory egg crate topper promotes comfort.

Patented Skoocher feature

moves mum closer to the perineal edge of the bed with the touch of a button.



Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065
61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Limited
PO Box 17-136
Greenlane, Auckland 1546, New Zealand
64 9 573 1890
www.stryker.com

1. Wang, Y., Beydoun, M.A., Liang, L., Caballero, B. and Kumanyika, S. K. (2008), Will All Americans Become Overweight or Obese? Estimating the Progression and Cost of the US Obesity Epidemic. *Obesity*, 16:2323-2330. doi:10.1038/oby.2008.351.

JMEH_1917_2_AOT
Copyright © 2020 Stryker

LD304 Birthing Bed

Standard features

- Lock-Rite slide-off foot section
- Lightweight, self-standing foot section
- Glideaway foot supports with six-way adjustability
- Single-button Trendelenburg
- Three-position siderails
- In-rail motion controls with lockout
- 44.5 cm (17.5") low bed height
- 12.7 cm (5") Ultra Comfort head and seat mattress
- 7.6 cm (3") foot mattress
- 15.2 cm (6") Omni Surface casters with wheel covers
- Locking-caster steering
- Four-wheel, steel-ring brake system with centrally located activation pedal
- 83.8 cm (33") mattress width with four motor functions (head, bed, foot, Trend.)
- Automatic pelvic tilt
- CPR release
- Manual backup
- Self-storing labour handgrips
- Plastic fluid basin
- Under-bed night light
- Wood headboard
- Permanent folding three-stage IV pole
- **Two year parts and labour warranty**

Optional features

- Patented Skoocher adjusts seat depth from 31.8 cm to 21.6 cm (12.5" to 8.5")
- Choice of removable or permanent calf supports
- Lumbar pillow and seat pillow
- Mattress overlay
- Basin liners
- 20.3 cm (8") casters
- Labour bar
- Linen for head and foot mattress

Bed Specifications

Model number	4701 (International 4701 model has 230 volt option)
Overall length	236.2 cm (93")
Overall width	
Siderails up	104.1 cm (41")
Siderails down	94.0 cm (37")
Lift capacity	
Bed	226.7 kg (500 lb.)
Foot section	136.1 kg (300 lb.)
Height range (to litter top)	
High	90.2 cm (35.5")
Low	44.5 cm (17.5")
Litter positioning	
Backrest	0° - 70°
Pelvic tilt	4°
Maximum Trend.	8°
Patient surface	84 x 206 cm (33" x 81")
Mattress thickness	
Head/seat	12.7 cm (5")
Foot	7.6 cm (3")
Seat depth	
Standard	31.8 cm (12.5")
With Skoocher	31.8 - 21.6 cm (12.5" - 8.5")
Caster diameter	
Standard	15.2 cm (6")
Optional	20.3 cm (8")
Electronics	
Agency approvals	C-US UL (120 V AC) CE (230 V AC)
Volts	120 V AC (Domestic)/230 V AC (International)
Ampere rating	10 A
Current leakage	< 300 microamperes
Frequency	60 Hz (Domestic)/ 50/60 Hz (International)
Hospital grade plug	Yes
EMC compatible	Yes
Motor function	Bead, head, foot, Trend., Skoocher (optional)

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local Stryker sales representative.