

# Michael Graves with Stryker Overbed Table



## Working together to redefine the healthcare environment

Simplified, streamlined, and incredibly lightweight, the Michael Graves with Stryker Overbed Table is designed for ease of use, mobility and easy cleaning. A large, welcoming surface provides plenty of space for personal items while obvious touch-points allow users to position the table to their comfort.



Combined with a sleek look, and less bulk, beautifully crafted hospitality handles allow caregivers and patients to more easily position the table for maximum use.



The U-shaped base fits easily around even the latest mobility solutions on beds and stretchers. A low profile allows for placement under recliners and side seating.



Intuitive activation handle allows users to move the table top easily with one hand. Simply push up or down to lift and lower.

**Dura-Edge** perimeter provides spill containment and durability over the life of the table

**Ergonomic activation handle** allows patients to easily raise and lower the table

**U-shaped base** slides around advanced mobility solutions on beds and stretchers without assistance and is compatible with most patient room equipment

**High-impact polystyrene cover** protects the base and conceals the casters

# Michael Graves with Stryker Overbed Table



## Standard features

- Concealed 1.46" (3.7 cm) dual-wheel casters
- No-grasp height adjustment handle
- Dura-Edge Spill containment top
- Hospitality handle

## Optional features

- Choose from four standard color configurations

## Specifications

**Model number** 4810-000-000

### Height

High 43.75" (111.1 cm)

Low 27.75" (70.5 cm)

**Safe working load** 20 lbs. (9 kg)

**Weight** 45 lbs (20.4 kg)

### Table top

Length 36.5" (92.7 cm)

Width 21.25" (54 cm)

### Base

Height 2" (5.1 cm)

Length 20" (50.8 cm)

Width 16" (40.6 cm)

**Caster diameter** 1.46" (3.7 cm)

## Warranty

Three years parts and labor