


Michael Graves with Stryker Patient Chair

Working Together To Redefine the Healthcare Environment



Designed to aid in nose-over-toes positioning,


the noteworthy aesthetics of the Michael Graves with Stryker Patient Chair pave the way for superior functionality. Ottoman is optional, and chair is available in two seat heights with highback, lowback, Comfort-Flex and armless configurations.



Beautifully designed extended chair arms promote independence in sitting and standing. This ergonomic design helps patients learn nose-over-toes positioning.



The practical design of the ottoman easily fits underneath the patient chair when not in use. Its height makes it functional for additional guest seating.



An optional Comfort-Flex backrest paired with standard polyurethane handrails help make positioning easy and comfortable for most body types.

Polyurethane arms provide a softer surface for grip and comfort while promoting easier ingress and egress

Steel frame offers added durability and maximum cleanability

Ottoman height allows for additional guest seating and makes it easy to store under the chair base

Four configurations provide seating options for all types of patient areas

Comfort-Flex backrest eases patients into a position that is most comfortable for their body type

Michael Graves with Stryker Patient Chair



3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Meets flammability standards of CAL 117
- Polyurethane handrails
- Metal base
- Meets ANSI/BIFMA X5.4-2005 per tests conducted in a certified test facility

Optional Features

- Choose from five standard color configurations
- CAL 133
- Ottoman
- Comfort-Flex (Highback Chair only)

Specifications

Style	Highback 18" Part # 4853-000-000	Comfort-Flex 18" Part # 4855-000-000	Lowback 18" Part # 4851-000-000
Height	43.0" (109.2 cm)	43.0" (109.2 cm)	35.0" (88.9 cm)
Width	26.5" (67.3 cm)	26.5" (67.3 cm)	26.5" (67.3 cm)
Depth			
Overall	31.5" (80.0 cm)	31.5" (80.0 cm)	30.5" (77.5 cm)
Seat	19.75" (50.2 cm)	19.75" (50.2 cm)	19.75" (50.2 cm)
Back Height from Seat	25.5" (64.8 cm)	25.0" (63.5 cm)	17.0" (43.2 cm)
Seat Height from Floor	19.0" (48.0 cm)	19.0" (48.0 cm)	19.0" (48.0 cm)
Arm Height from Floor	26.0" (66.0 cm)	26.0" (66.0 cm)	26.0" (66.0 cm)

Style	Highback 22" Part # 4854-000-000	Comfort-Flex 22" Part # 4856-000-000	Lowback 22" Part # 4852-000-000
Height	47.0" (119.4 cm)	47.0" (119.4 cm)	39.5" (100.3 cm)
Width	26.5" (67.3 cm)	26.5" (67.3 cm)	26.5" (67.3 cm)
Depth			
Overall	31.5" (80.0 cm)	31.5" (80.0 cm)	30.5" (77.5 cm)
Seat	19.75" (50.2 cm)	19.75" (50.2 cm)	19.75" (50.2 cm)
Back Height from Seat	25.5" (64.8 cm)	25.0" (63.5 cm)	17.0" (43.2 cm)
Seat Height from Floor	23.5" (60.0 cm)	23.5" (60.0 cm)	23.5" (60.0 cm)
Arm Height from Floor	30.0" (76.2 cm)	30.0" (76.2 cm)	30.0" (76.2 cm)

Style	Armless 18" Part # 4850-000-000	Ottoman Part # 4857-000-000
Height	35.0" (88.9 cm)	13" (33 cm)
Width	22.0" (55.9 cm)	18.5" (47 cm)
Depth		
Overall	28.0" (71.1 cm)	17.5" (44.5 cm)
Seat	19.75" (50.2 cm)	
Back Height from Seat	17.0" (43.2 cm)	
Seat Height from Floor	19.0" (48.0 cm)	
Arm Height from Floor	N/A	
Weight Capacity (ottoman and all models)	300 lbs (136 kg)	

Warranty

- Ten years on steel frames, arms and wood frame components
- Three years on Stryker Medical Textiles and foam
- Warranted to 300 lbs (136 kg) weight capacity