

# Prime Series Stretcher



## Mobile. Ergonomic. Configurable.

Our Prime Series Stretchers are designed to help reduce physical strain on clinicians, enhance patient comfort and help drive hospital efficiency. Transporting heavier patients is now easier with Stryker's advanced mobility; while electric functionality gives patients the independence to adjust themselves without calling for caregiver assistance. Fully equipped and configurable, our Prime Series Stretchers allows you to create the ideal solution to meet your specific needs.



### Advanced Stretcher mobility:

**Zoom with Motorised Drive**  
Virtually eliminates manual pushing.

**Big Wheel** Reduces start-up push force and steering.



**Electric-powered lift** eliminates manual pumping, allowing for quick patient access and caregiver safety. Simply touch the hydraulic pedal to raise the litter.



**Lift Assist Backrest** uses the patient's weight to raise the head of the stretcher while the drop seat design helps reduce boosting. Also available as a fully electric option.

**Recovery chair position** helps eliminate the need for unnecessary transfers by providing a comfortable seated position.



**Patient controls** The patient has control of their own comfort with the staggered control panels on the Glideaway siderails.

**IsoFlex SE Gel Support Surface<sup>1</sup>** offers pressure redistribution and comfort as a component to aid in the treatment and prevention of pressure injuries.

**Power-washable<sup>2</sup> design** with 318 kg safe working load.

<sup>1</sup>Recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by health care professional  
<sup>2</sup>1105 and 1115 models only

# Prime Series Stretcher



## Standard features

- 318 kg weight capacity
- Less than 54 cm low height
- Litter-mounted oxygen bottle holder
- Power-washable design (1105 and 1115 only)
- Center-locking brake system with dual end activators
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- Dual end brake/steer controls
- Glideaway Siderails
- 4 integrated IV receptacles
- +/- Trendelenburg
- 20 cm Omni Surface Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Thermoformed ABS base
- Integrated utility tray
- Integrated transfer board

## Warranty

Two years parts, labour and travel

## Electric functionality (30 inch only)

- Electric patient controls
- Foot end nursing controls
- Powered lift
- 120V or 240V available

## Optional features

- Big Wheel Advanced Mobility (1115)
- Zoom Stretcher with Motorised Drive System (1125)
- Support surfaces:
  - IsoFlex SE Support Surface
  - ComfortGel SE Support Surface
  - UltraComfort SE Support Surface
- Four-sided brake/steer controls (1105 and 1115)

## Specifications

	66 cm litter	76 cm litter
<b>Model number</b>	66 cm (26")	76 cm (30")
Fifth Wheel Mobility	1105	1105
Big Wheel Advanced Mobility	1115	1115
Zoom Motorised Drive System	1125	1125
<b>Overall length</b>		
1105 and 1115	215.9 cm (85")	215.9 cm (85")
1125	218.4 cm (86")	218.4 cm (86")
<b>Overall width</b>		
Siderails up	86.4 cm (34")	96.5 cm (38")
Siderails down	76.8 cm (30.25")	77.5 cm (30.5")
Siderails down (1125)	78.7 cm (31")	78.7 cm (31")
<b>Patient surface</b>	66 x 191 cm (26" x 75.25")	76 x 191 cm (30" x 75.25")
<b>Weight capacity</b>	318 kg (700 lb)	318 kg (700 lb)
<b>Height range</b>		
High	86.4 cm (34")	86.4 cm (34")
Low (1105 and 1115)	52.7 cm (20.75")	52.7 cm (20.75")
Low (1125)	58 cm (23")	58 cm (23")
<b>Litter positioning</b>		
Backrest	0°– 90°	0°– 90°
Knee gatch	0°– 40°	0°– 40°
Trend./Reverse Trend. (1105 and 1115)	±17°	±17°
Trend./Reverse Trend. (1125)	±16°	±16°
<b>Siderail length</b>	36 x 147 cm (14" x 57.75")	36 x 147 cm (14" x 57.75")
<b>Caster diameter</b>	20 cm (8")	20 cm (8")
<b>EMC</b>	IEC 60601-1-2	

- Pop-up push handles
- Electric powered knee gatch
- Lift Assist Backrest design with drop seat
- Three-sided hydraulic controls (1105 and 1115)
- Dual-ended siderail release
- IV pole
  - Permanent two- or three-stage
  - Removable

- Serving tray
- Defibrillator tray
- Defibrillator tray/foot extender/charting surface
- Chartholder
- IV caddy
- Serving tray holder/footboard
- Siderail pads
- Upright oxygen bottle holder
- Restraint straps
- Retractable cord reel (1125)

## Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker and Prime. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



PRIME-COM-3EN-20269

PRIME-Com-3EN 20269  
SDL 06/2019  
2019-20269

stryker.com

Copyright © 2019 stryker.com