

stryker

ProCeed™ Hospital Bed

An epicenter of
patient care





40% to 80%

of nurses report experiencing symptoms of psychological distress, globally.¹



78,740

nurses in the U.S. suffer on-the-job injury or illness annually.²



40-50%

of nurses report developing back pain over the course of their career.³



stryker

Designed for **durability**. Designed for **care**.
Designed for **ease**. Designed for **you**.

Your hospital beds aren't just furniture—they play a key role on a patient's journey back to health. In the face of worldwide nursing shortages and staff burn out,⁴ it's crucial that your equipment not only endure the relentless workload but also provide an intuitive, straightforward experience for the staff that rely on it.

Meet ProCeed — a hub for healthcare,
elevating experiences for nurses, staff and patients alike.
Designed with input from clinicians around the world
for durability and cleanability, ProCeed helps you deliver
efficient quality care.



Built for **durability**

You work hard for your patients and you deserve equipment that works hard for you. From cleanability to serviceability, ProCeed was built to withstand the demands of healthcare.

Cleanable

With an IPX6 rating, ProCeed's frame and electronic panel are resilient against powerful jets of water, allowing you to confidently clean your hospital bed. The support surface features welded seams to prevent water ingress.



Durable design

With a 260 kg safe working load, the bed has been rigorously strength tested. The entire frame, including wheels, siderails, headboard, footboard and accessories, have been comprehensively tested for impact resistance.⁵

Simple to service

ProCeed is designed for ease of serviceability and can be part of your ProCare service plan, so you can keep your fleet performing as it should.

ProCare Services

Built for **patient and care team safety**

We know safety is a top priority for you. With a low bed height, bed status indicators, and full coverage siderails, ProCeed was created with both staff and patient safety in mind.

Ergonomically designed, full coverage folding siderails

ProCeed's ergonomically designed, full coverage siderails help keep patients in bed when they need to be, and fold down and tuck away when they don't.



Status indicators on siderail

With indicators on the siderail, nurses can see from a distance the status of the battery, brakes, AC power, bed height, bed angle and more.



Low height

ProCeed's low bed height of 34 cm not only offers an ergonomic position for nurses and staff to assist patients, but also supports patient mobility and reduces their risk of injury from falls.



Built for **efficiency and ease of use**

A nurse's job is anything but easy. That's why ProCeed was thoughtfully designed to ease the burden of everyday tasks like x-rays, patient transportation and repositioning.

Radiolucent backrest

With ProCeed's radiolucent backrest, staff can take x-rays without transferring the patient to another surface, saving valuable time and energy.

Headboard options

Choose between a headboard that stays stationary when the frame tilts, or a mobile headboard that tilts, raises and lowers as the bed moves into different positions.



Fifth wheel

With fifth wheel mobility, an extra wheel at the core of the bed offers heightened control, reduces start-up push force and steering helping to reduce caregiver back injuries.



Customer Support Services

Technical Support

Our Technical Support comprises a team of professionals available to help with your ProCeed needs.

Stryker's ProCare Services

With ProCare Services, our trained technicians help ensure your equipment is ready to perform when it's needed. ProCare is a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

Flex Financial Program

Through Stryker's Flex Financial business, we provide a range of smart alternatives designed to fit your organisation's needs. We collaborate with our customers for assistance in finding the financial solution that suits your needs.*

*Program availability will vary depending on country and product portfolio. Please check with your local sales representative for full details and availability

References: **1:** James Buchan; Howard Catton. Recover to rebuild: Investing in the nursing workforce for health system effectiveness. March 2023. ICN. [https://www.icn.ch/sites/default/files/2023-07/ICN_Recover-to-Rebuild_report_EN.pdf] **2:** Illness, Injuries, and Fatalities. 11 February 2022. U.S. Bureau of Labor Statistics. 15 May 2024. [[https://www.bls.gov/iif/snapshots/osn-registered-nurses-2016-20.htm#:~:text=Nonfatal%20occupational%20injuries%20and%20illnesses,\(FTE\)%20workers%20was%20390.6.](https://www.bls.gov/iif/snapshots/osn-registered-nurses-2016-20.htm#:~:text=Nonfatal%20occupational%20injuries%20and%20illnesses,(FTE)%20workers%20was%20390.6.)] **3:** Rayhan A. Tariq; Jemima S. George; George Ampat; Tammy J. Toney-Butler. Back Safety. 22 July 2023. National Library of Medicine. 15 May 2024. [<https://www.ncbi.nlm.nih.gov/books/NBK519066/#:~:text=Of%20all%20healthcare%20workers%2C%20nurses,%5B23%5D%20%5B24%5D.>] **4:** Jester R. Editorial - Global shortage of nurses - Rebecca Jester for May 2023 issue. Int J Orthop Trauma Nurs. 2023 May;49:101018. doi: 10.1016/j.ijotn.2023.101018. Epub 2023 Mar 24. PMID: 37041090; PMCID: PMC10038672. **5:** test data on file

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Flex Financial, ProCare, ProCeed. All other trademarks are trademarks of their respective owners or holders. ProCeed : Class I. The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

