

# ProForm<sup>®</sup> Support Surface



ProForm support surface is used to assist in the prevention and treatment of pressure injuries. ProForm's foam design helps to promote immersion and envelopment, helping to enhance pressure redistribution.

## Specifications

<b>Model number</b>	2815
<b>Overall length</b>	
Standard	80"/84" (203.2 cm/213.4 cm)
Behavioral health	80"/84" (203.2 cm/213.4 cm)
<b>Overall width</b>	35" (88.9 cm)
<b>Thickness</b>	6" (15.2 cm)
<b>Safe working load</b>	500 lb. (226.8 kg)

## Warranty

Three-year, non-prorated warranty on mattress and cover

## Features

- Non-powered design
- Fire barrier\*
- Latex free
- Dartex top cover (behavioral health only)

## Optional features

- Fire barrier\* (standard only)
- Nylon top cover (standard only)
- Dartex top cover (standard only)

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

**stryker.com**

**ProForm is a non-powered, pressure redistributing support surface** through its multi-layered foam.

**360 watershed** prevents fluid intrusion

**High density** siderail foam provides a supportive crib to assist with ingress and egress

**Covered foam crib** increases patients surface area.

\*Fire retardant inner cover meets flammability standard CAL TB129, Boston IX-11, 16 CFR 1632, 16 CFR 1633, and CGSB/CAN 2-4.2 Method 27.7-M-77.

Stryker reserves the right to change specifications without notice.

Dartex is a registered trademark of Dartex Coating Limited Corporation.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: ProForm, Stryker. All other trademarks are trademarks of their respective owners or holder.