

SideKick® Disinfecting Wipes



Compatible. Versatile. Convenient.

SideKick is a cleaning and disinfecting wipe that is convenient and effective at disinfecting healthcare surfaces with a two-minute wet contact time. Its broad-spectrum efficacy is tuberculocidal, bactericidal, fungicidal and virucidal. SideKick wipes are the only cleaning and disinfecting solution that is 100% compatible¹ for use on our beds, stretchers, patient support surfaces², furniture, chest compression systems, defibrillators, emergency care cots, cot fasteners and chairs – guaranteed.



Compatible

Tested 100% compatible on over 200 materials. Use SideKick wipes with confidence knowing you will not degrade or discolor our beds, stretchers, patient support surfaces², furniture, chest compression systems, defibrillators, emergency care cots, cot fasteners and chairs.¹



Versatile

SideKick wipes' broad-spectrum efficacy is tuberculocidal, bactericidal, fungicidal and virucidal. It was developed to clean and disinfect a variety of healthcare surfaces.²



Convenient

SideKick is a premoistened wipe formulated with a two-minute wet contact time. The wipes are ready to use, meaning there is no mess or potential for chemical-mixing error.

SideKick disinfecting wipes are proven effective for the following:²

Bactericidal

- Mycobacterium bovis, BCG (tuberculosis) (TB)
- Escherichia coli-with ESBL (E. coli)
- Carbapenem-Resistant Klebsiella pneumoniae (CRKP)
- Methicillin-Resistant Staphylococcus aureus (MRSA)
- Salmonella enterica
- Vancomycin-Resistant Enterococcus faecium (VRE)
- Acinetobacter baumannii
- Enterobacter aerogenes
- Listeria monocytogenes
- Pseudomonas aeruginosa
- Serratia marcescens
- Staphylococcus aureus

Fungicidal

- Trichophyton interdigitale (ringworm)

Virucidal

- Avian Influenza A Reassortant Virus (H3N2)
- Hepatitis B Virus* (HBV)
- Herpes Simplex II Virus (HSV 2)
- Rhinovirus Type 14 (RV-14)
- Human Immunodeficiency Virus* (HIV-1)
- Rotavirus Strain WA
- Hepatitis C Virus* (HCV)
- Swine Influenza A Virus (H1N1)
- Influenza A Virus Strain A2/Hong Kong

SideKick Disinfecting Wipes



Standard features

- Compatible on Stryker's products¹
- Broad-spectrum efficacy
 - Tuberculocidal
 - Bactericidal
 - Fungicidal
 - Virucidal
- 2-minute wet contact time
- Neutral pH of 7.6
- Alcohol quat formulation
- For use on many healthcare surfaces²
- EPA reg number: 70144-2-92728

Active ingredients

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride: **0.154%**
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride: **0.154%**
- Isopropanol: **21.000%**
- Other ingredients: **78.692%**

Specifications

	Standard wipes	Large wipes
Model number	2060-000-001	2060-000-002
Wipe dimensions	6" x 10"	9" x 12"
Count per canister	100	65
Weight of canister	1 lb. 7.8 oz. (676.6g)	1 lb. 11.9 oz. (790.6g)

Store in a cool, well-ventilated area.

Stryker reserves the right to change specifications without notice.

Always refer to the SideKick wipes label for additional information and instructions for use. References and/or data for claims are on file at Stryker and available upon request.

***HIV-1, HBV and HCV Precautions:** This product kills HIV-1, HBV, and HCV on pre-cleaned hard, non-porous surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

1. Compatibility testing refers to the validation of Product surface materials performed by Stryker to not degrade, discolor or crack with exclusive use of SideKick for daily device cleaning and disinfection for the Product's stated expected life. Subject to footnote 2 below. Products refer to beds, stretchers, patient support surfaces, furniture, chest compression systems, defibrillators, emergency care cots, cot fasteners and chairs manufactured by and/or for Stryker.

2. Cleaning and disinfection claims exclude the Equilibrium support surface cover. Equilibrium by Dartex® and Dartex® are registered trademarks of Dartex Coatings, Ltd.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: SideKick, Stryker. All other trademarks are trademarks of their respective owners or holder.

Mkt Lit-1537 12 DEC 2017 Rev D
Copyright © 2018 Stryker
Printed in the U.S.A.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

stryker.com/sidekick