

SV2 Hospital Bed

Customisable to fit your needs

The next standard of care.

Creating a hospital bed like our SV2 takes collaboration. By working together with medical professionals like you all over the world, we are delivering a high-quality bed you can customise to your specific needs.

The SV2 is a feature-rich, value-driven bed that includes electronic controls with convenient, one-touch button positioning for both you and your patients, a low bed indicator, mobile headboard, optional fifth wheel, optional under-bed lighting and much more.



Low bed indicator

Clearly indicates to caregivers when the SV2 is in its lowest position, helping to reduce patient injury from falls.



Fifth wheel mobility

Enhances maneuverability, helping to reduce caregiver back injuries.



Smart egress

One-touch button brings bed surface to its exiting height, helping to reduce patient falls.

Fold-down siderails

Simplify patient exiting, while enhancing safety by fully covering and minimising entrapment gaps. Allowing ease of use for caregivers by tucking away with only one hand.

Double regression

Allows backrest and knee-gatch regression to help reduce both sacral pressure and displacement of the patient.

Mobile headboard

Easily adjustable bed height for a more ergonomic positioning for the caregivers.

"Always-on" underbed light

Offers patients better sight during the night when exiting the bed, helping to reduce falls.

Making a difference. Reducing our impact.

When we make products like our SV2 bed, we seek ways to use less and waste less – because we consider the wellbeing of our planet to be as significant as the health of the individuals who live on it.

SV2 Hospital Bed



Stryker's Regional Headquarters
 Herikerbergweg 110
 1101 CM, Amsterdam
 The Netherlands

patientcare.eu.stryker.com

Standard features

- Battery backup
- Linak DC motors (for quiet operation)
- Low bed indicator
- Manual CPR on both sides of backrest
- Electric lift
- Electric legrest
- Electric backrest
- Lockable and tuck-away siderails
- Single-wheel casters
- Removable litter covers for cleaning
- Removable headboard and footboard
- Dual Linak siderail controls
- Ten (10) years of life expectancy
- Four corner brake, neutral, steer pedals
- One anti-static caster to ground frame
- Crash bumpers
- Easy adjustable foot prop for knee-bend positioning and foot elevation
- Electro-static painted metal frame
- Angle indicators for backrest
- Angle indicators for Trendelenburg and Anti-Trendelenburg (bed angle)
- IPX4 liquid splash protection
- IEC 60601-1:2012 (edition 3.1)
- IEC 60601-2-52:2009 + A1:2015 (edition 1.1)
- IEC 60601-1-2:2014 (edition 4.0)

Optional features

- 5th wheel
- Nurse control unit
- Patient pendant
- Under bed light
- X-ray translucent backrest and cassette holder
- Bed extender with linen tray
- dual-wheel casters

Packages

SV2 package plus

Includes all standard features, plus:

- Nurse pendant
- Dual wheel casters
- 5th wheel

SV2 package complete

Includes all standard features, plus:

- Under bed light
- Bed extender with linen tray
- Dual wheel casters
- Patient pendant
- X-ray translucent backrest and cassette holder
- Nurse pendant
- 5th wheel

Specifications

	SV2 Hospital Bed
Model number	7500
Safe working load	250 kg
Note: Safe working load indicates the sum of the patient, mattress and accessory weight	
Maximum patient weight	215 kg
Product weight	150 kg
Overall product size	
Length	220 cm (±1 cm)
Length with bed extender	251 cm (±1 cm)
Width	99 cm (±1 cm)
Product height (without mattress)	
High	75 cm (±1 cm)
Low	37 cm (±1 cm)
Under product clearance	15 cm
Caster size (single and dual)	15 cm
Product angle indicator	0° - 15°
Backrest angle indicator	0° - 90°
Backrest angle	0° - 65°
Knee gatch angle	0° - 30°
Trendelenburg/Reverse Trendelenburg	0° - ±12°
Electronics	100 - 240 VAC, 50 - 60 Hz

Specifications

Mattresses

	7002-2-012 basic foam	7002-2-014 basic foam	7002-5-012 visco foam
Model number	7002-2-012 basic foam	7002-2-014 basic foam	7002-5-012 visco foam
Model with flame barrier (meets EN 597-1, EN 597-2, BS6807, UNI9175)	7002-2-512 basic foam	7002-2-514 basic foam	7002-5-512 visco foam
Model with flame barrier (meets BS7177)	N/A	N/A	7002-5-712 visco foam
Length	200 cm	200 cm	200 cm
Width	86 cm	86 cm	86 cm
Height	12 cm	14 cm	12 cm
Weight	7.1 kg	8.2 kg	8.6 kg

Specifications

Bolsters for bed extender

	7002-4-018 with 12 cm mattress	7002-4-020 with 14 cm mattress
Model number	7002-4-018 with 12 cm mattress	7002-4-020 with 14 cm mattress
Model with flame barrier	7002-4-518 with 12 cm mattress	7002-4-520 with 14 cm mattress

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Active Sensor Technology, IsoAir, SV2. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



2017-14264EN
 DFS 05/2017

Copyright © 2017 Stryker