

# Symmetry<sup>®</sup> II Bariatric Chair

Comprehensive Patient Room  
Seating Solutions



Sleek and sturdy for your patients' needs

Durability and design unite to offer a high-performing Bariatric Chair. Transitional styling and clean lines provide a classic look that is easily integrated in all lounge and patient room environments.

## Standard Features

- 750 lbs dynamic weight capacity
- Removable seat and slip-covered cushions facilitate maintenance, cleaning and enhance durability
- High-density Ultracell foam promotes maximum comfort and resiliency
- Perma-Mesh™ suspension with 2" steel frame provides additional support
- Accentuated lumbar design supports lower back
- Arched, open-arm design with extended arm caps permits easy entry and exit
- Limited lifetime warranty

## Optional Features

- Choose from a variety of fabrics and stains
- CAL 133

## Specifications

Part Number	3500-000-420
Overall Width	35" (89 cm)
Overall Height	37" (94 cm)
Depth	
Overall	25" (64 cm)
Seat	18" (46 cm)
Weight Capacity	750 lbs (340 kg)
Weight	70 lbs (32 kg)
Seat Width	31.5" (80 cm)
Seat Height From Floor	18.5" (47 cm)
Arm Height From Floor	26.5" (67 cm)
Back Height from Seat	18.5" (47 cm)

## Warranty

Limited lifetime warranty

Stryker reserves the right to change specifications without notice.  
Furniture is manufactured exclusively for Stryker Medical.