

# Symmetry II

## Recliner Bed



### Versatility & resilience for quality care.

With Stryker's Recliner Bed, expect durability without compromising comfort. Three highly-versatile positions (chair, lounge and full flat-bed) will help guests feel right at home in any patient room.

#### Standard features

- 8 cm (3") back swivel casters
- 8cm (3") fixed front locking casters
- CAL 117

#### Optional features

- Choose from a variety of fabrics and stains
- CAL 133

#### Warranty

- Limited lifetime warranty

#### Specifications

Model	3500-000-650
Overall width	88.9 cm (35")
Overall height	90.2 cm (35.5")
Overall depth	94.0 cm (37")
Seat depth	53.3 cm (21")
Lounge position depth	167.6 cm (66")
Flat-bed position depth	215.9 cm (85")
Weight	59.0 kg (130 lbs)
Seat width	68.6 cm (27")
Seat height from floor	44.5 cm (17.5")
Arm height from floor	66.0 cm (26")
Weight capacity	136.1 kg (300 lb.)

**Easily converts** to lounge or sleep positions.

**High-resiliency polyurethane foam** and fiber-wrap cushioning provide maximum comfort and lasting shape.

Equipped with **front-locking casters**.

**Limited lifetime warranty.**



**Lounge Position**



**Full Flat-bed Position**

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local sales representative.

**Stryker Australia Pty Ltd**  
 8 Herbert Street  
 St Leonards NSW 2065  
 61 2 9467 1000  
 www.stryker.com.au

**Stryker New Zealand Limited**  
 PO Box 17-136  
 Greenlane, Auckland, New Zealand  
 64 9 573 1890  
 www.stryker.com