

TruRize®

Clinical Chair



The next step in mobility.

Early patient mobility is a crucial component during the recovery process. Its importance to a patient's health and recovery cannot be emphasized enough.*

TruRize® helps to promote early patient mobility. It gives caregivers the assurance they need to carefully mobilise their patients.

Key Features

- Simple, user-friendly technology that equips caregivers with the tools to help promote efficient patient handling
- Powered-pivoting design guides the movement of TruRize® into positions that help promote patient comfort
- An industry first integrated chair exit alarm system may help to shorten response time to patient falls*



A



B



C



D



E



F



- A Pivoting armrests
- B Caregiver control panel
- C Occupant control panel

- D One-touch brakes
- E Integrated Foley bag hooks
- F Vertical transportation handles

Stand Assist Function supports the patient to a standing position, helping to promote early patient mobility.

* Early Mobilization in the Intensive Care Unit: A Systematic Review
Joseph Adler, PT, DPT, CCS and Daniel Malone, PhD, MPT, CCS

* Powered, clinical chair intended for hospital use only.



Standard Features

- Powered chair
- Dual-sided Caregiver specific control panels
- Dual-sided lighted occupant control panel with lockout controls
- Six positions
(Stand assist, Recline, Upright 2, Upright 1, Transfer, Trend)
- Height Adjustability
- Integrated chair exit alarm system complaint with IEC 60601-1-8 medium priority alarm
- Pivoting armrests
- One-touch brake
- Vertical handles for maneuvering
- Two integrated Foley bag hooks
- Two line management hooks
- 3.5" (8,9 cm) Casters
- Four fabric color options
(Meets CAL 117 flammability)

Optional Features

- Battery Backup
- Nurse Call Cable (6.35 mm / 1/4" jack)
- CAL 133 Fire Barrier

Specifications

Model Number	3333-000-000
Overall Width (Upright 1)	86 cm (34")
Overall Height (Upright 1)	135 to 155 cm (53" to 61")
Overall Length (Full-Flat)	185 cm (73")
Seat	
Width	66 cm (26")
Depth	52.1 cm (20.5")
Height from floor (Upright 1)	48 to 66 cm (19" to 26")
Arm Height from Floor	71 to 89 cm (28" to 35")
Weight Capacity	158 kg (350 lbs.)
Electronics	
Agency Approvals	IEC60601-1
Volts	100-240 VAC
Ampere Rating	4 A
Frequency	50/60 Hz
Battery (optional)	(2) 12V DC lead acid batteries

Warranty

- Two years for parts

Stryker reserves the right to change specifications without notice.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: TruRize Clinical Chair and Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

OTHMED-BRO-18

RRD 2016/06

2016-9734

Copyright © 2016 Stryker

Website: patientcare.eu.stryker.com