

stryker



TruRize

Clinical Chair

Truth in mobility.

The next step in mobility.

Early patient mobility is crucial to the recovery process hence its importance to a patient's health must be emphasised.¹

With TruRize, Stryker has created a powered clinical chair with a new way of thinking about mobility.

TruRize helps to promote early patient mobility, giving caregivers the assurance and help they need to get their patients out of bed early.

It's only a chair by de inition: TruRize is the next step forward in patient mobility.

Key features

- Simple, user-friendly technology helps to equip caregivers with the tools to help improve patient handling.
- Powered-pivoting design guides the movement of TruRize into positions that help promote patient comfort.
- Integrated chair exit alarm system may help shorten response time to patient falls.



Vertical transportation handles allow caregivers of a wide range of heights to easily push, manoeuvre and configurate TruRize.

Lighted occupant control panel is designed to allow patients to control their own comfort while limited full chair mobility to the hands of the caregiver.



Caregiver specific control panel is easily accessible and easy to use, the control panel offers simplicity in a complex environment.

Vertical handles are designed for optimal ergonomic manoeuvring.

Pivoting armrests move to aid in patient transfer.

Height adjustable to accommodate a wide range of patients.

Integrated Foley bag hooks, located on both sides of TruRize, position drainage bags below the waist.

One-touch brakes, located on both sides of the chair, simultaneously lock all four dual-wheels.

stryker

Functionality & configurability.



Stand assist function is designed to support patients to a standing position, helping to promote early patient mobility.

Early patient mobility

Long-term patient immobility puts a patient's recovery and overall health at risk. TruRize is designed to promote early patient mobility by supporting the patient to a standing position.

Patient comfort

From its external design to the functionality of its features, TruRize has been carefully engineered with patient comfort in mind. As a clinical chair with an integrated exit alarm system, TruRize may help shorten response time to patients when they fall.¹

¹Integrated chair exit alarm system that is IEC 60601-1-8 Medium Priority Alarm compliant.

Caregiver confidence


Caregiver injuries can occur as a result of manual lifting, moving or shifting of patients.² TruRize has been designed to help promote better patient handling.

Full-flat transfer position is intended to aid in patient transfers with pivoting armrests to aid in the process. Simply align TruRize with the patient's bed for easy transfers.



Six set positions

The convenient caregiver-specific control panel allows caregivers to choose from six set positions with a simple push of a button.



- A** Integrated chair exit alarm
- B** Stand assist function
- C** Height adjustability
- D** Full-flat transfer position
- E** Integrated Trendelenburg position
- F** Occupant control lock
- G** Brake indicator



1. Trendelenburg position



2. Full-flat transfer



3. Recline position



4. Upright one



5. Upright two



6. Stand Assist

Customer support services

Technical support

Stryker technical support comprises a team of professionals available to help with your TruRize needs. Contact via phone at 1 800 667 558 or email at ssptechservices@stryker.com.

ProCare service from Stryker

Product reliability is just the beginning of the Stryker commitment to lifetime customer satisfaction. Stryker utilises advanced metrics to ensure equipment uptime and assist customers in achieving their patient and caregiver goals.

ProCare service solutions are customisable and scalable to any budget. All ProCare offerings are supported by our talented factory-trained Stryker technicians, ensuring all equipment and documentation is maintained to the highest standards.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex financial program

Our financial programs provide a range of smart alternatives designed to fit your organisation's needs. We offer flexibility beyond a cash purchase with payment structures that can be customised to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia Pty Ltd

8 Herbert Street
St Leonards NSW 2065
61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Limited

PO Box 17-136
Greenlane
Auckland
New Zealand
64 9 573 1890
www.stryker.com

JMEH_1938_3_AOT
Copyright © 2020 Stryker

1. Berry A, Beattie K, Bennett J, Chaseling W, Cross Y, Cushway S, Hassan A, Jones S, Longhurst E, Moore R, Phillips D, Plowman E, Scott J, Smith K, Thomas L and Elliott D (2017). Physical Activity and Movement: a Guideline for Critically Ill Adults. Agency for Clinical Innovation NSW Government ISBN 978-1-74187-976-6
2. Safe Work Australia, Priority industry snapshot: Health care and social assistance, June 2018.