



# Tray layout card – Gamma4 Indication tray



Insert layer



Ref #	Description	Qty	Ref #	Description	Qty
<b>A</b> 1420-1010	Indication Tray Insert	1	<b>Optional</b>		
<b>1</b> 1420-0300	Extraction Shaft	1	<b>4</b> 2351-0240	One Shot Device	1
<b>2</b> 1420-1020	Free Space Mat	1	<b>Spare part</b>		
<b>3</b> 1320-0234	Set Screwdriver, Flexible Shaft	1	1500-0040	Full Size Tray Lid	1
			<b>Kit</b>		
			1420-0000	Indication Kit	1

# Tray layout card – Gamma4 Indication tray

Base layer



Ref #	Description	Qty	Ref #	Description	Qty
<b>B</b> 1420-1005	Indication Tray Base	1	<b>8</b> 1420-0220	Precision Sleeve	1
<b>1</b> 2351-6000	Opening Reamer Handle	1	<b>9</b> 1420-0160	Lag Screw Sleeve	1
<b>2</b> 1420-0050	Opening Sleeve	1	<b>10</b> 1420-0260	Lag Screw Driver (two parts)	1
<b>3</b> 1420-0055	Opening Guide	1	<b>11</b> 1420-0240	Lag Screw Reamer	1
<b>4</b> 1420-0080	Opening Reamer	1	<b>Spare part</b>		
<b>5</b> 1420-0110	Targeting Sleeve (two parts)	1	1420-0112	Targeting Sleeve Knob	1
<b>6</b> 1420-0100	Proximal Targeting Arm	1	<b>Kit</b>		
<b>7</b> 1420-0105	Nail Holding Screw	2	1420-0000	Indication Kit	1

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Gamma4, Precision Sleeve, Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia  
 8 Herbert Street,  
 St Leonards NSW 2065  
 T: 61 2 9467 1000 F: 61 2 9467 1010

Stryker New Zealand  
 PO Box 17-136  
 Greenlane Auckland 1546  
 T: 64 9 573 1890 F: 64 9 573 1891