

Count Sheet

IMN Basic Instruments Tray: 23560580-T

Tray Insert

Ref. #	Description	Location	Qty	Check box
2351-0101	Self-retaining Screwdriver Sleeve, Short	Insert Layer	1	
2351-0100	Screwdriver Bit, Short	Insert Layer	1	
2351-0105	Screwdriver Bit, Medium	Insert Layer	1	
2351-0040	Ball Tip Screwdriver	Insert Layer	1	
2351-0106	Self-retaining Screwdriver Sleeve, Medium	Insert Layer	1	
2351-0180	Extraction Shaft	Insert Layer	1	
1806-0150	T2 Strike Plate	Insert Layer	1	
2351-0110	Screwdriver Bit, Long	Insert Layer	1	
2351-0070	Tissue Protection Sleeve, Long	Insert Layer	1	
2351-0160	Freedhand Depth Gauge, Short	Insert Layer	1	
2351-4280	Locking Drill Sleeve, Long	Insert Layer	1	
2351-0170	Freedhand Depth Gauge, Long	Insert Layer	1	
2351-4290	Locking Trocar, Long	Insert Layer	1	
2351-0340	Screw Scale	Insert Layer	1	
2351-0150	Guided Depth Gauge	Insert Layer	1	
Total Instrument Count			15	
2355-5010	Basic IMN Insts - Tray Insert	Insert Layer	1	

Tray Base

Ref. #	Description	Location	Qty	Check box
2351-0420	X-Ray Ruler	Base Layer	1	
2351-0020	Reduction Rod	Base Layer	1	
2351-0140	Quick-lock Delta Handle, Modified Trinkle	Base Layer	1	
2351-0030	Guide Wire Handle	Base Layer	1	
2351-0010	Curved Awl	Base Layer	1	
2351-0011	Awl Plug	Base Layer	1	
2351-0380	Guide Wire Pusher	Base Layer	1	
1806-0022	Guide Wire Ruler	Base Layer	1	
2351-0400	Compression Screwdriver	Base Layer	1	
2351-0060	Slotted Hammer	Base Layer	1	
1806-0110	Universal Rod	Insert Layer	1	
Total Instrument Count			10	
2355-5005	Basic IMN Insts - Tray Base	Base Layer	1	

Hospital Signature _____

Count Sheet

IMN Basic Instruments Tray: 23560580-T

Additional Items

Ref. #	Description	Location	Qty	Check box
Total Instrument Count				

Hospital Signature _____

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/ or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate related entities own, use or have applied for the following trademarks or service marks: SOMA, Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia
 8 Herbert Street,
 St Leonards NSW 2065 T: 61 2
 9467 1000
 F: 61 2 9467 1010

Stryker New Zealand
 PO Box 17-136 Greenlane
 Auckland 1546
 T : 64 9 573 1890
 F: 64 9 573 1891

stryker.com