


**stryker**

**Quality you can trust,**

customisation you need




**Surgeon iD**

Mandible reconstruction

Surgeon iD allows you to create your own pre-bent mandible reconstruction plate. Plates are designed to meet your individual needs and help reduce intraoperative adaptation time. Surgeon iD provides you with an on-site, customised and pre-bent solution.

- Select your plate design features on personalised design session
- 2.0mm or 2.8mm plate profile height
- Planned on normative mandible models
- Available in small, medium and large sizes
- Utilise sizer tool to determine appropriate plate size for each patient

**A customised solution,**  
ready when you are



### **Personalised design session**

- Select profile height, length and run of your plate
- Determine screw hole placement
- Customised plate strength by increasing bar widths

## Craniomaxillofacial

**This document is intended solely for the use of healthcare professionals.**

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Surgeon iD and Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.

CMF-FL-124EN\_Rev. None\_23991  
SDL 08/2020  
2020-26576  
Copyright © 2020 Stryker  
stryker.com



Stryker Leibinger  
GmbH & Co. KG  
Bötzingen Straße 41  
79111 Freiburg  
Germany



0197