

Power-LOAD

Stetcher fastener system.



Key features

- Lifting arms
- Head-end lock LED indicators
- Control panel
- Stretcher release handles
- Foot-end lock LED indicator
- Linear transfer system
- Inductive charging
- Operation guide

Optional features

- Wheel guide

Warranty

- One year parts, labour, and travel
7-year service life.

Specifications*

Model number	6390
Length	
Overall length	95 in (241 cm)
Minimum length	85 in (216 cm)
Width	24.5 in (62 cm)
Weight	
Total weight	211.5 lb (96.5 kg)
Floor plate assembly	16.5 lb (7.5 kg)
Anchor assembly	23 lb (10.5 kg)
Transfer assembly	67 lb (30.5 kg)
Trolley assembly	105 lb (48 kg)
Maximum weight capacity*	700 lb (318 kg)
Minimum operators required	
Occupied stretcher	2
Unoccupied stretcher	1
Recommended loading height	22 in to 36 in (56 cm to 91 cm)
Battery	12V, 5 Ah lead acid battery (6390-001-468)

*Maximum weight capacity represents patient weight. Safe working load of 870 lb (395 kg) represents the sum of the stretcher total weight and patient.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The Power-LOAD stretcher fastener system is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822.

The yellow and black colour scheme is a proprietary trademark of the Stryker Corporation.

Power-LOAD is currently not available, specifications subject to change without notice. Pending 510(K) clearance, not available for sale. The device has not been authorised as required by the rules of the Federal Communications Commission. The device is not, and may not be, offered for sale or lease, or sold or leased, until authorisation is obtained.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Patents pending.

For more information contact your local sales representative

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065

Stryker New Zealand Ltd
515 Mt Wellington Highway
Mt Wellington, Auckland 1060

61 2 9467 1000
www.stryker.com.au

64 9 573 1890
www.stryker.com



Mkt Lit-1243 JULY 2016 REV A.2
Copyright © 2016 Stryker