

Power PRO XT

Powered ambulance stretcher.



Key features

- Automatic fastener shut-off when loaded
- Manual back-up release handle
- Automatic high-speed retract
- Battery-powered hydraulic lift system
- Preventative maintenance hour meter
- Battery pack power level indicator
- Adjustable load height with jog function
- High visibility powder-coated frame
- Lightweight, rugged aluminum construction
- Lower lifting bar
- Dual height lift bar
- Lift-capable safety bar
- Integrated bumper system
- Retractable head section
- Highly visible controls
- Floor-mounted safety hook
- One-hand release, fold down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Oversized wheels with sealed caster and wheel bearings
- Sealed bolster mattress
- Shock or flat leg positioning
- Built-in pull handle
- X-frame guards
- Power washable
- SMRT power system (12V DC, 120V AC and 240V AC options available)
- X-restraints

Optional features

- Heavy duty two- or three-stage IV poles (patient right or left)
- Permanent or removable O2 bottle holders (head end, foot end or fowler)
- Sealed flat mattress
- Grey mattress
- Dual wheel locks
- Head extension
- Pillow
- Equipment hook
- Backrest storage pouch
- Head end storage flat
- Defibrillator platform
- Base storage net
- Knee gatch
- SMRT charger mounting bracket
- Cobalt blue restraints
- PowerLOAD compatibility
- PerformanceLOAD compatibility

Warranty

- Two-year parts, labour and travel
- One-year soft goods
- Three-year X-frame components
- Three-year limited powertrain
- Lifetime on all welds*

Extended warranties available.

* 7-year service life.

For more information contact your local sales representative

Stryker Australia Pty Ltd 8 Herbert Street
St Leonards NSW 2065

Stryker New Zealand Ltd 515 Mt Wellington Highway
Mt Wellington, Auckland 1060

61 2 9467 1000
www.stryker.com.au

64 9 573 1890
www.stryker.com

Specifications

Model number	6506
Height¹ (infinite height positioning between lowest and highest position)	
Highest position	41.5 in (105 cm)
Lowest position	14 in (36 cm)
Length	
Standard	81 in (206 cm)
Minimum	63 in (160 cm)
Width	23in (58 cm)
Weight²	125 lb (57 kg)
Wheels	
Diameter	6 in (15 cm)
Width	2 in (5 cm)
Articulation	
Backrest	0-73°
Shock position	+ 15
Optional knee gatch	30°
Maximum weight capacity³	700 lb (318 kg)
Minimum operators required	
Occupied stretcher	2
Unoccupied stretcher	1
Recommended fastener system	
Power-LOAD	Model 6390
Floor mount	Model 6370
Wall mount	Model 6371
Recommended loading height⁴	Up to 36 in (91 cm)

¹ Height measured from bottom of mattress, at seat section, to ground level.

² Stretcher is weighed with one battery pack, without mattress and restraints.

³ 700 lb weight capacity with an unassisted lift capacity of 500 lb.

⁴ Can accommodate load decks up to 36 in. Load height can be set between 26 in and 36 in.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The Power-PRO XT is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822. The Power-PRO XT is designed to be compatible with competitive stretcher fastener systems.

The yellow and black colour scheme is a proprietary trademark of the Stryker Corporation. Power-LOAD Compatibility option currently not available, specifications subject to change.

Pending 510(K) clearance, not available for sale.

The device has not been authorised as required by the rules of the Federal Communications Commission. The device is not, and may not be, offered for sale or lease, or sold or leased, until authorisation is obtained.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Patents pending.

