

stryker

SMRT
power system
Powered ambulance stretcher.



SMRT Power system.

For Stryker powered ambulance stretchers

Stryker managed recharging technology

The SMRT Power System will give your Power-PRO stretcher superior performance. This is a tough, professional-grade system designed for the high-pressure world of EMS. The SMRT Power System eliminates the time-consuming and elaborate charging protocols and "tuning" usually required for maintaining reliable, high performance battery-powered systems.


With the SMRT Power System you get:

- Faster charging — approximately 2 hours
 - Zero memory—repeated partial charges do not diminish SMRT Paks
 - Automatic power management extends battery life
 - Charge SMRT Paks when and where needed—convenient in-ambulance and wall outlet charging options
 - Sealed durable construction for all-weather toughness
 - Two SMRT Paks, one SMRT Charger and a power cord.
- Keep one SMRT Pak in the SMRT Charger and the other in your Power-PRO stretcher



Change a SMRT Pak with one hand, in seconds

To remove the SMRT Pak from your stretcher, press the red release button and slide the SMRT Pak out. This can be done with one hand. Slide in a fresh SMRT Pak until you hear a click.



Result

SMRT Power technology automatically protects against damage from over-charging and high temperatures that dramatically shorten battery life in unpredictable EMS applications. Under harsh test conditions, SMRT Power outperforms and outlasts comparable batteries.

Actual battery performance and life experienced by individual EMS services may differ greatly due to variables including environmental factors, patient size, call frequency and ambulance stretcher condition.



SMRT charger automatic charging

SMRT Charger automatically evaluates the condition of each SMRT Pak and cannot overcharge or overheat SMRT Paks.

SMRT Pak (2)

Powerful and dependable. Long run times give you confidence on every call.



Accessories



12v DC cable (automotive)



240V AC - Australia (AS-3112 without ground pin)



SMRT charger mounting bracket



Customer support services

Technical support

Stryker technical support comprises a team of professionals available to help with your InTouch needs. Contact via phone at 1 800 667 558 or email at ssptechservices@stryker.com.

ProCare service from Stryker

Product reliability is just the beginning of the Stryker commitment to lifetime customer satisfaction. Stryker utilises advanced metrics to ensure equipment uptime and assist customers in achieving their patient and caregiver goals.

ProCare service solutions are customisable and scalable to any budget. All ProCare offerings are supported by our talented factory-trained Stryker technicians, ensuring all equipment and documentation is maintained to the highest standards.

All ProCare agreements provide:

- Stryker-authorised service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex financial program

Our financial programs provide a range of smart alternatives designed to fit your organisation's needs. We offer flexibility beyond a cash purchase with payment structures that can be customised to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Reference

1. Testing conducted by Intertek, an independent global provider of product and certification test services. www.intertek.com
2. A "Simulated EMS Call" consists of four ambulance stretcher extension/retraction cycles performed during a representative EMS service call.
3. Each test cycle consists of a full charge followed by a sequence of simulated EMS Calls until pack voltage dropped below 18.8 V. To simulate the expected power conditions and battery management practices, an over-charge cycle was included with every EMS Call.

Stryker reserves the right to change specifications without notice.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia Pty Ltd

8 Herbert Street
St Leonards NSW 2065

61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Ltd

515 Mt Wellington Highway
Mt Wellington, Auckland 1060

64 9 573 1890
www.stryker.com