

**XPS**

expandable  
patient surface



# XPS – expandable to fit a wide variety of patients and environments.

## XPS is there when you need it


XPS provides an expanded patient surface area that can be easily retrofitted to compatible cots. XPS is adjustable with seven locking positions and includes a wider mattress that reduces transfer gap and is designed with your patient's comfort in mind. XPS helps address the growing obesity trend and supports a variety of patients and environments.

Cots with XPS meet current tip stability standards for both patient and provider safety<sup>1</sup>. Engineered for durability and cleanability, with a 6061-T6 aluminum core and over-mold design.

1. Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL. BS-EN 1789 for Power-PRO XT, Power-PRO TL and Performance-PRO XT.


### Standard features

- Increased patient surface area
- Adjustability for patients and environments
- Seven locking positions
- Integrated into cot — always there
- Easily retrofitted to existing compatible cots
- Enhanced safety built into release handles
- Durable aluminum over-mold design
- Mattress design reduces transfer gap
- Enhanced patient comfort
- Compliant with tip stability and dynamic crash certifications<sup>1</sup>



### Seven locking positions

Between 10 and 52 degrees



### Increased surface area

From 23 to 33 inches



### Option or retrofit XPS compatible cots

Power-PRO XT model 6500/6506

Power-PRO TL model 6550

Performance-PRO XT model 6085/6086



## Safety

Compliant with tip stability and dynamic crash certifications



## Durability

Engineered for durability with an aluminum core


## Always there

Integrated into cot for accessibility on every call



## Release handles

Enhanced safety built into release handles



## Reduced transfer gap

Mattress design reduces transfer gap



## XPS specifications

<b>Includes</b>	
2 XPS siderails XPS mattress	
<b>Weight<sup>1</sup></b>	8 lb (4 kg)
<b>Seven locking positions</b>	Between 10° – 52°
<b>Surface area expansion</b>	
<b>Locking angle</b>	<b>Width – total surface area</b>
10°	23 in (58 cm)
17°	25 in (64 cm)
24°	27 in (69 cm)
31°	29 in (74 cm)
38°	30 in (76 cm)
45°	32 in (81 cm)
52°	33 in (84 cm)
<b>Height</b>	10 in (25 cm)
<b>Length</b>	30 in (76 cm)
<b>XPS mattress width<sup>2</sup></b>	23 in (58 cm)
<b>XPS compatible cots<sup>3</sup></b>	
Power-PRO XT	Model 6500/6506
Power-PRO TL	Model 6550
Performance-PRO XT	Model 6085/6086

1. Additional weight compared to standard siderails not including mattress
2. Width is measured at widest point (standard bolster mattress 19 in / 48 cm)
3. Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL. BS-EN 1789 for Power-PRO XT, Power-PRO TL and Performance-PRO XT.

Stryker reserves the right to change specifications without notice.  
All numbers rounded to nearest whole value.

## Technical Support

Our Technical Support comprises a team of professionals available to help with your XPS needs. Contact via phone at 1 800 STRYKER or email at [medicaltechnicalsupport@stryker.com](mailto:medicaltechnicalsupport@stryker.com)

## Stryker's ProCare Services

Every day, you count on your medical equipment to perform at its best. With ProCare Services, our people help to ensure your equipment is ready to perform when it's needed and make it easier to get the most from your investment. When an issue arises, we promise that we'll solve it — performing repairs quickly and correctly.

ProCare isn't just a service program. It's a partnership you can count on to give you one less thing to worry about, and one more reason to feel confident you're doing all you can for your clinicians, staff and patients.

All ProCare agreements provide:

- Stryker-authorized service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time\*

\*Based on the provisions of the Service Agreement and the location of the product.

## Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures that can be customized to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks:

**Flex Financial, Performance-PRO, Power-PRO, ProCare, Stryker, XPS.**

All other trademarks are trademarks of their respective owners or holders.

3800 E. Centre Avenue  
Portage, MI 49002 U.S.A.  
t: 269 329 2100  
f: 269 329 2315  
toll free: 800 327 0770

[ems.stryker.com](http://ems.stryker.com)