

stryker

XPS

Expandable patient surface.



XPS

Expandable patient surface.

The all-in-one solution for a variety of patients and environments.

XPS is there when you need it.

XPS provides an expanded patient surface area that can easily be retrofitted with compatible stretchers. XPS is adjustable with 7 locking positions and includes a wider mattress that reduces transfer gap and is designed with patient comfort in mind. This solution helps address growing obesity trends and supports a variety of patients and environments. Stretchers with XPS meet current tip stability standards for both patient and provider safety¹. Engineered for durability and cleanability, with a 6061-T6 aluminium core and over-mold design.

- Increased patient surface area
- Adjustability for patients and environments
- 7 locking positions
- Integrated into stretcher – Always there
- Easily retrofitted to existing compatible stretchers
- Enhanced safety built into release handles
- Durable aluminium over-mold design
- Mattress design reduces transfer gap
- Enhanced patient comfort
- Compliant with tip stability and dynamic crash certifications³

Safety

Compliant with tip stability and dynamic crash certifications.



Release handles

Enhanced safety built into release handles.



7 locking positions

Between 10° and 52°



Increased surface area

From 58 cm to 84 cm



Option or retrofit

Power-PRO XT model 6500/6506



stryker

Durability

Engineered for durability with an aluminium core.

Always there

Integrated into stretcher for accessibility on every call.

Reduced transfer gap

Mattress design reduces transfer gap.





Customer support services

Technical support

Stryker technical support comprises a team of professionals available to help with your InTouch needs. Contact via phone at 1 800 667 558 or email at ssptechservices@stryker.com.

ProCare service from Stryker

Product reliability is just the beginning of the Stryker commitment to lifetime customer satisfaction. Stryker utilises advanced metrics to ensure equipment uptime and assist customers in achieving their patient and caregiver goals.

ProCare service solutions are customisable and scalable to any budget. All ProCare offerings are supported by our talented factory-trained Stryker technicians, ensuring all equipment and documentation is maintained to the highest standards.

All ProCare agreements provide:

- Stryker-authorised service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex financial program

Our financial programs provide a range of smart alternatives designed to fit your organisation's needs. We offer flexibility beyond a cash purchase with payment structures that can be customised to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Reference

1. Additional weight compared to standard siderails not including mattress
2. Width is measured at widest point (Standard bolster mattress 19 in / 48 cm)
3. Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL.
BS-EN 1789 for Power-PRO XT.

Stryker reserves the right to change specifications without notice.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065

61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Ltd
515 Mt Wellington Highway
Mt Wellington, Auckland 1060

64 9 573 1890
www.stryker.com