

## XPS

Expandable patient surface.



### Key features

- Increased patient surface area
- Adjustability for patients and environments
- 7 locking positions
- Integrated into stretcher – always there
- Easily retrofitted to existing compatible stretchers
- Enhanced safety built into release handle
- Durable aluminum over-mold design
- Mattress design reduces transfer gap
- Enhanced patient comfort
- Compliant with tip stability and dynamic

### Specifications

<b>Weight<sup>2</sup></b>	8 lb (4kg)
<b>7 Locking positions</b>	Between 10°-52°
<b>Surface area expansion</b>	
<b>Locking angle</b>	Width-total surface area
<b>10°</b>	23 in (58 cm)
<b>17°</b>	25 in (64 cm)
<b>24°</b>	27 in (69 cm)
<b>31°</b>	29 in (74 cm)
<b>38°</b>	30 in (76 cm)
<b>45°</b>	32 in (81 cm)
<b>52°</b>	33 in (84 cm)
<b>Height</b>	10 in (25 cm)
<b>Length</b>	30 in (76 cm)
<b>XPS mattress width<sup>3</sup></b>	23 in (58 cm)
<b>XPS compatible stretchers</b>	
<b>Power-PRO XT</b>	Model 6506

<sup>1</sup>Certified to IEC 60601-1 for Power-PRO XT and Power-PRO TL.  
<sup>2</sup>BS-EN 1789 for Power-PRO XT, Power-PRO TL, and Performance-PRO XT.  
<sup>3</sup>Additional weight compared to standard siderails not including mattress  
<sup>3</sup>Width is measured at widest point (Standard bolster mattress 19 in / 48 cm)

Stryker reserves the right to change specifications without notice.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information contact your local sales representative

**Stryker Australia Pty Ltd** 8 Herbert Street  
St Leonards NSW 2065

**Stryker New Zealand Ltd** 515 Mt Wellington Highway  
Mt Wellington, Auckland 1060

61 2 9467 1000  
www.stryker.com.au

64 9 573 1890  
www.stryker.com