



Genuine Accessories

McGRATH™ MAC EMS

Video Laryngoscope

McGRATH™ MAC EMS

Video Laryngoscope






**Thank you for choosing Physio-Control
as your partner in helping save lives and
improve patient care.**

The McGRATH™ MAC EMS video laryngoscope is shipped with one battery and the Instructions for Use. Disposable blades are shipped separately. In this catalog you will find an overview of blades and additional accessories available for the McGRATH™ MAC EMS. If you have any questions about accessories or disposables, please contact your Physio-Control representative.

Disposable Laryngoscope Blades


McGRATH™ MAC Disposable Laryngoscope Blades are high performance, low-cost, single-use sterile blades. They are used for routine as well as many difficult airways, from pediatric to adult patients. They have a Macintosh-like curvature and are therefore easy to use with limited training. McGRATH™ MAC Blades are made of robust optical polymer and blade lenses are coated with hydrophilic optical surface treatment.



McGRATH™ MAC 2 Disposable Laryngoscope Blade

Designed for pediatric patients older than eight weeks or at least 4.5 kg. Packed individually and shipped in boxes of 10.


Box of 10 11996-000414



McGRATH™ MAC 3 Disposable Laryngoscope Blade

Designed for adult patients. Packed individually and shipped in boxes of 10.

Box of 10 11996-000415




McGRATH™ MAC 4 Disposable Laryngoscope Blade

Designed for large adult patients. Packed individually and shipped in boxes of 10.


Box of 10 11996-000416

McGRATH™ X blade™

The McGRATH™ X blade™ is a sterile disposable single-use blade that extends the range of proprietary disposable laryngoscope blades for the McGRATH™ MAC EMS system. The X blade™ contains a sweeping acute curvature for extreme airways, complementing the existing range of McGRATH™ MAC Disposable Laryngoscope Blades which are utilized in routine to difficult airways.



Replacement Batteries



The McGRATH™ MAC EMS utilizes 3.6V Lithium power cells, which provide up to 250 minutes of usage. A representation of minutes remaining is displayed on screen. The on/off power switch is integrated into the battery pack, so when a new battery is fitted, the on/off function is renewed, keeping the product young.

McGRATH™ 3.6V EMS Battery

The McGRATH™ MAC EMS is shipped with one 3.6V EMS Battery. Replacement batteries are shipped individually.
11996-000394

McGRATH™ MAC EMS

Video Laryngoscope





Physio-Control is now part of Stryker.

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 800.895.5896 (Canada) or visit our website at www.physio-control.com

Physio-Control Headquarters

11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support

P. O. Box 97006
Redmond, WA 98073
Toll free 800 442 1142
Fax 800 426 8049

Physio-Control Canada

Physio-Control Canada Sales, Ltd.
45 Innovation Drive
Hamilton, ON
L9H 7L8
Canada
Toll free 800 895 5896
Fax 866 430 6115