


stryker



Full flexibility
for better care

LIFEPAK® 15 monitor/defibrillator

Versatile. Intuitive. Powerful.

Take control of patient care with one of the most adaptable emergency response devices found in any area of the hospital. Flexing up for critical events, this industry-leading, easy-to-use monitor/defibrillator delivers immediate insights for quick action and is ready when you need it most. Discover a fully loaded monitoring solution with escalating energy to 360J, advanced CPR capabilities and enhanced visibility that helps make your job easier and improve patient outcomes.

✕ ✕ Flexes up



- Ideal for anywhere in the hospital—from crash carts to advanced monitoring
- Standardize equipment and streamline staff training
- Adapt to protocols and care needs across departments with highly configurable options



Easy to use

- One button, one job—quickly activate therapies for all levels of providers
- Respond confidently with easy access to a simplified keypad
- Receive CPR data for effective post-event review through standard defibrillation pads



Powerful

- Treat difficult-to-defibrillate patients with escalating energy up to 360J
- Enhance CPR quality through physiological feedback using EtCO₂ trending
- Speed up care team activations with simple acquisition and transmission of 12-lead reports to other facilities or departments



Dependable

- Highest user satisfaction rating in all individual categories (MD Buyline, Q3 2019)
- Track device readiness with virtual device monitoring and automatic alerts for status changes
- Service your device with customized plans and onsite biomed training



Maximum flexibility at a glance

Easily manage STEMI patients with ST-segment trend monitoring and 12-lead ECG transmissions via the LIFENET® System

Dual, long-lasting lithium-ion batteries deliver six hours of monitoring—additional option for AC power

Durable, shock-absorbing handle and double-layered screen withstands falls from bed height or drops in transit



Keep a consistent compression rate without the need for extra external hardware or expensive pads with a clear CPR metronome

Large screen with high resolution maximizes visibility and includes a third channel for monitoring waveforms

Easily collect patient and CPR performance data for accurate documentation and quick analysis in CODE-STAT™ data review software and service

Help identify return of spontaneous circulation (ROSC) and gauge CPR quality with waveform ranges as low as 0-20 mmHg (consistent with AHA Guidelines) using Integrated Medtronic Microstream™ EtCO₂ monitoring

Reinforced cable connector gives you the confidence for secure therapy delivery

Reduce the need for extra handhelds with options for advanced monitoring including NIBP, EtCO₂, SpO₂, CO, Met, 12-lead ECG, temperature and invasive pressure

100mm printer offers clear visibility of waveforms and patient summary data



A legacy
of trust

strykeremergencycare.com

For a brief summary of important disclosure and safety information regarding LIFEPAK devices, visit our website at <https://www.strykeremergencycare.com/disclosure-and-safety-information>

Please consult Operating Instructions at www.strykeremergencycare.com or call 800.442.1142 for complete list of indications, contraindications, warnings, cautions, potential adverse events, safety and effectiveness data, instructions for use and other important information.

Emergency Care

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, LIFENET, LIFEPAK, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

For further information, please contact Stryker at 800 442 1142 (U.S.), 800 668 8323 (Canada) or visit our website at strykeremergencycare.com.

Copyright © 2021 Stryker. 3301599_J



Manufactured by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052
U.S.A.
Toll free 800 442 1142
strykeremergencycare.com

Distributed in Canada by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2 Canada
Toll free 800 668 8323