



stryker

2023 ENT
product
catalog

Balloon dilation

Nasal implants

Cryotherapy

Bioresorbables

Navigation

Visualization

Surgical instruments

We are Stryker

**Together with our customers,
we are driven
to make healthcare better**

**Our mission is at the heart of what we do
and believe.**

We collaborate with our customers to develop innovative products and services that ultimately improve the lives of patients. We understand that in today's challenging healthcare environment, our customers need solutions that improve quality and efficiency in the operating room and, ultimately, the patient experience. These challenges are creating opportunities for us to make healthcare better, for our customers and the patients they serve.

Stryker is one of the world's leading medical technology companies, and, together with our customers, we are driven to make healthcare better. Stryker offers a diverse array of innovative products and services in ENT, Orthopaedics, Medical and Surgical, Neurotechnology and Spine that help improve patient and hospital outcomes. Stryker is active in over 100 countries around the world.

Follow @StrykerENT



Balloon dilation

XprESS ENT dilation system	03
Audion ET dilation system	04

Nasal implants

LATERA absorbable nasal implant system	05
--	----

Cryotherapy

ClariFix cryotherapy device	06
-----------------------------	----

Bioresorbable dressings

NasoPore bioresorbable nasal dressing	07
HemoPore nasal dressing	08
XeroGel nasal/epistaxis pack	08
OtoPore bioresorbable ear dressing	09

Navigation

Stryker ENT navigation system and accessories	10
Stryker ENT navigation system instruments	12

Visualization

FocESS image solutions	14
FocESS sinusscopes	15


Surgical instruments

ESSx microdebrider	16
MiniFESS turbinate forceps	17
MiniFESS Blakesley thru-cut forceps	18
MiniFESS Takahashi forceps	19
MiniFESS sickle knife	20
MiniFESS Blakesley forceps	21
MiniFESS seekers	22
MiniFESS Light Seeker sinus confirmation tool	23
MiniFESS surgical instrumentation	24
Reinforced anesthesia needle	25
Suction and irrigation tools	26

Balloon dilation

XprESS[®]

ENT dilation system



Product	Ordering no.	Description
XprESS LoProfile system	LPLF-105	5 x 20 mm, 1 system
XprESS LoProfile system	LPLF-106	6 x 20 mm, 1 system
XprESS LoProfile system	LPLF-107	7 x 20 mm, 1 system
XprESS LoProfile system	LPLF-205	5 x 8 mm, 1 system
XprESS LoProfile system	LPLF-206	6 x 8 mm, 1 system
XprESS Ultra system	ULF-105	5 x 20 mm, 1 system
XprESS Ultra system	ULF-106	6 x 20 mm, 1 system
XprESS Ultra system	ULF-205	5 x 8 mm, 1 system
XprESS Ultra system	ULF-206	6 x 8 mm, 1 system

1 system includes: 1 balloon dilation device, 1 inflation syringe, 1 extension line, 1 bending tool,
1 PathAssist LED Light Fiber (5110-001)

Balloon dilation

Audion[™]

ET dilation system



Product	Ordering no.	Description
Audion ET dilation system	AET-100	1 system

1 system includes: 1 balloon dilation device, 1 inflation syringe, 1 extension line, 1 inflation lock

Nasal implants

LATERA[®]

absorbable nasal implant system



Product	Ordering no.	Description
LATERA implant system	LATSYS20	20 mm, 1 system
LATERA implant system	LATSYS24	24 mm, 1 system

1 system includes: 1 implant delivery device, 2 absorbable nasal implants

Cryotherapy

ClariFix[®]

cryotherapy device



Product	Ordering no.	Description
ClariFix system	CFX-2000	1 system

1 system includes: 1 ClariFix cryotherapy device, 2 cryogen canisters and 2 sterile canister caps

NasoPore®

bioresorbable nasal dressing



Product	Ordering no.	Description
NasoPore dressing	5400-010-004	Standard, 4 cm, 8-pack
NasoPore dressing	5400-010-008	Standard, 8 cm, 8-pack
NasoPore dressing	5400-020-004	Firm, 4 cm, 8-pack
NasoPore dressing	5400-020-008	Firm, 8 cm, 8-pack
NasoPore FD dressing	5400-020-108	Fast degrading, 8 cm, 8-pack
NasoPore dressing	5400-030-004	Extra firm, 4 cm, 8-pack
NasoPore dressing	5400-030-008	Extra firm, 8 cm, 8-pack
NasoPore dressing	5400-212-004S	Standard, 4 cm, 2-pack
NasoPore dressing	5400-222-004S	Firm, 4 cm, 2-pack

HemoPore® and XeroGel®

hemostatic dressing with chitosan



Product	Ordering no.	Description
Hemopore nasal dressing		
HemoPore dressing	5400-020-208	8 cm, 8-pack
XeroGel nasal/epistaxis pack		
XeroGel pack	XG-102	4.0 L x 2.4 W x 0.3 cm thick, box of 2 individual sterile units
XeroGel pack	XG-108	4.0 L x 2.4 W x 0.3 cm thick, box of 8 individual sterile units

Bioresorbable dressings

OtoPore® bioresorbable ear dressing




Product	Ordering no.	Description
OtoPore dressing	5400-010-000	Cylinder, standard, 2 x 1 cm, 8-pack
OtoPore dressing	5400-020-000	Cylinder, firm, 2 x 1 cm, 8-pack
OtoPore dressing	5400-020-100	Square, firm, 2.0 L x 4.0 W x 0.15 cm thick, 8-pack


Navigation

Stryker ENT navigation system and accessories

Scopis ENT software



Scopis ENT software with TGS (Target Guided Surgery)



Electromagnetic navigation unit

Product	Ordering no.	Description
Stryker ENT navigation system		
Scopis ENT software	8000-020-001	Basic software for navigation in ENT applications
Scopis ENT software with TGS (Target Guided Surgery) technology	8000-020-002	Premium software for planning and navigation in ENT applications with advanced planning features and augmented reality technology
Electromagnetic navigation unit	8000-010-003	High-performance navigation PC with electromagnetic measuring technology


Stryker ENT navigation system and accessories (cont.)



Field generator




Field generator mounting arm



Headrest

Product	Ordering no.	Description
Navigation platform accessories		
Field generator	8000-010-004	Accessory used for the electromagnetic position measurement of patient trackers and navigated instruments
Field generator mounting arm	8000-010-005	Flexible mounting arm to position the field generator
Universal mounting arm	8000-010-008	Smaller arm for flexible positioning of the field generator
Pro equipment cart	8000-030-002	Includes: four shelves and a storage drawer
Medical keyboard	8000-030-010	U.S./international keyboard, wired, number pad, USB
Mouse, wired	8000-030-020	Wired, USB, three buttons, scroll function
Mouse, wireless	8000-030-021	Wireless, USB, scroll function
Headrest - adapter A package	8000-011-005	TGS universal headrest, mounting arm, adapter A, 0.25" x 2.0"
Headrest - adapter B package	8000-012-005	TGS universal headrest, mounting arm, adapter B, 0.25" x 2.5"
Headrest - adapter C package	8000-013-005	TGS universal headrest, mounting arm, adapter C, 0.25" x 3.25"
Field generator stand	8000-014-005	Stand with mounting arm to position the field generator

Stryker ENT navigation system instruments



Instrument clamp, forceps



Instrument clamp, 2-6 mm



Instrument clamp, 6-10 mm



Instrument clamp, 10-16 mm

Product	Ordering no.	Description
Navigation instruments		
Patient tracker electromagnetic	8000-040-001	Disinfectable, used to locate the patient within the EM field
Registration pointer electromagnetic	8000-050-003	Disinfectable, used for surface registration
Navigation tool extension cable	8000-050-011	Extension cable for navigation tools
Instrument clamp, forceps	8000-060-010	For attaching the universal tracker to rigid conventional surgical instruments
Instrument clamp, 2-6 mm	8000-060-011	For attaching the universal tracker to rigid conventional surgical instruments between 2-6 mm
Instrument clamp, 6-10 mm	8000-060-012	For attaching the universal tracker to rigid conventional surgical instruments between 6-10 mm
Instrument clamp, 10-16 mm	8000-060-013	For attaching the universal tracker to rigid conventional surgical instruments between 10-16 mm
Patient tracker tabs	8000-100-001	Adhesive pads for fixation of patient trackers electromagnetic, 100 per unit

Stryker ENT navigation system instruments (cont.)



TGS guidewire Pointer and suction tubes Instrument clamp sphere Instrument clamp universal

Product	Ordering no.	Description
Navigation instruments, single use		
TGS guidewire	8000-060-009	Allows navigation of XprESS LoProfile system
Navigation instruments, maximum 10 uses		
Patient tracker electromagnetic	8000-040-002	Used to locate the patient within the EM field
Precision pointer electromagnetic	8000-050-001	1.5 mm, malleable
Suction tube, Frazier electromagnetic	8000-050-005	3 mm, malleable, straight suction
Suction tube, Eicken electromagnetic	8000-050-006	3 mm, malleable, curved suction
Endoscope tracker electromagnetic	8000-060-001	Navigation suitable for 0°, 30° and 45° endoscopes with 4 mm diameter and 180 mm length
Calibration body electromagnetic	8000-060-002	Calibration of suitable conventional surgical instruments
Universal tracker electromagnetic	8000-060-006	Allows navigation of rigid conventional surgical instruments in combination with an instrument clamp
Instrument clamp sphere	8000-060-020	Adapter to track rigid conventional surgical instruments between 2-4.6 mm
Instrument clamp universal	8000-060-021	3-6 mm, universal tracker


FocESS[®] image solutions



FocESS HD wireless camera

Product	Ordering no.	Description
FocESS HD wireless camera	WCS-100	Includes: 1 receiver, 1 battery charger and 2 batteries
Monitor, 23"	OPM-100-23	23" HD monitor
Monitor, 32"	OPM-100-32	32" HD monitor

FocESS[®] sinusscopes



Product	Ordering no.	Description
FocESS sinuscope, 0°	ES-0	3 mm, 0°, autoclavable
FocESS sinuscope, 30°	ES-30	3 mm, 30°, autoclavable
FocESS sinuscope, 45°	ES-45	3 mm, 45°, autoclavable
FocESS sinuscope, 70°	ES-70	3 mm, 70°, autoclavable
FocESS sterilization tray	TRA-102	Sterilization tray for FocESS sinuscopes; autoclavable; holds 3 sinuscopes and 6 light post adapters

ESSx[®] microdebrider



Product	Ordering no.	Description
Microdebrider individual components		
Irrigation cassettes	5290-075-000	10-pack
ESSx microdebrider	5290-601-200	Handpiece
Microdebrider blades		
Aggressive blade	5290-011-000	4.0 mm, 11 cm, 5-pack
Sharp edge blade	5290-340-000	4.0 mm, 8 cm, 5-pack
Aggressive serrated blade, 40°	5290-435-100	3.5 mm, 11 cm, 5-pack
Aggressive serrated blade, 40°	5290-440-100	4.0 mm, 11 cm, 5-pack
Jaguar blade	5290-530-000	3.5 mm, 8 cm 5-pack
Jaguar blade	5290-540-000	4.0 mm, 8 cm, 5-pack
Aggressive serrated blade	5290-628-000	2.5 mm, 8 cm, 5-pack
Aggressive serrated blade	5290-635-000	3.5 mm, 8 cm, 5-pack
Aggressive serrated blade, 60°	5290-635-100	3.5 mm, 8 cm, 5-pack
Aggressive blade	5290-638-000	3.5 mm, 8 cm, 5-pack
Aggressive serrated blade, 60°	5290-640-100	4.0 mm, 11 cm, 5-pack
Aggressive serrated blade	5290-645-000	4.0 mm, 8 cm, 5-pack
Hooded rounded bur	5290-647-000	3.0 mm, 8 cm, 5-pack
Aggressive blade	5290-648-000	4.0 mm, 8 cm, 5-pack
Angled aggressive blade	5290-740-100	4.0 mm, 11 cm, 5-pack
Convex angled aggressive blade	5290-740-200	4.0 mm, 11 cm, 5-pack
Angled aggressive blade, 12°	5290-744-100	4.0 mm, 11 cm, 5-pack
Angled aggressive blade, 60°	5290-760-100	4.0 mm, 11 cm, 5-pack

MiniFESS[®]

turbinate forceps



MiniFESS turbinate forceps

Product	Ordering no.	Description
MiniFESS turbinate forceps	CS-101	110 mm working length, 5 x 19 mm jaws

MiniFESS[®]

Blakesley thru-cut forceps



MiniFESS Blakesley thru-cut forceps

Product	Ordering no.	Description
MiniFESS Blakesley thru-cut forceps	CS-102	110 mm working length, straight, 2 x 5 mm jaws
MiniFESS Blakesley thru-cut forceps	CS-103	110 mm working length, 45°, 2 x 5 mm jaws

MiniFESS® Takahashi forceps



Product	Ordering no.	Description
MiniFESS Takahashi forceps	CS-104	Straight, 2.5 x 8 mm jaws

MiniFESS® sickle knife



Product	Ordering no.	Description
MiniFESS sickle knife	CS-105	Curved, 10 mm

MiniFESS®

Blakesley forceps



Product	Ordering no.	Description
MiniFESS Blakesley forceps	CS-106	Straight, 110 mm working length, 2.5 x 7 mm jaws
MiniFESS Blakesley forceps	CS-107	45°, 110 mm working length, 2.5 x 7 mm jaws
MiniFESS Blakesley forceps	CS-108	90°, 110 mm working length, 2.5 x 7 mm jaws

MiniFESS®

seekers



Product	Ordering no.	Description
MiniFESS sphenoid seeker/freer	FS-100	Sphenoid seeker/freer
MiniFESS maxillary seeker	MS-120	120° angle
MiniFESS maxillary seeker	MS-135	135° angle
MiniFESS maxillary seeker set	MS-255	Includes: both 120° angle and 135° angle

MiniFESS® Light Seeker®

sinus confirmation tool



Product	Ordering no.	Description
MiniFESS Light Seeker	IS-100	78°, 1.5 mm ball tip, includes standard ACMI connection plus 1 Wolf and 1 Storz connector adapter
Portable red LED light source	PRLS-100	Portable, non-sterile, reusable light source to be used with the MiniFESS Light Seeker, includes 2 batteries and 1 battery charger

MiniFESS®


surgical instrumentation



Product	Ordering no.	Description
MiniFESS surgical instrument set	CST-1000	Includes: 3 Blakesley forceps (straight, 45°, 90°), 2 Blakesley thru-cut forceps (straight, 45°), 1 Takahashi forceps, 1 turbinate forceps, 1 sickle knife, and 1 sterilization tray
MiniFESS seeker set	PASS-100	Includes: 2 maxillary seekers (120° and 135° angle), 1 sphenoid seeker/freer and 1 MiniFESS Light Seeker
MiniFESS sterilization tray	TRA-101	Autoclavable; holds 8 finger-grip and 6 seeker-style instruments

RAN


(reinforced anesthesia needle)



Product	Ordering no.	Description
Reinforced anesthesia needle	RAN-027-2	27 G, 3.5", 2-pack (only sold in a bundle configuration)
Reinforced anesthesia needle	RAN-027-5	27 G, 3.5", 5-pack

Suction


and irrigation tools



Product	Ordering no.	Description
Cliq suction pump	LCP-AS-1001D	Includes: 800 mL collection canister, microbial filter, intermediate tubing and suction tubing
Cyclone system	CC-100	Includes: 1 Cyclone device, 1 syringe, 1 bending tool

Unmatched support

Stryker provides the technology, training and service you need to help you improve your patient outcomes. With Stryker, you gain access to a comprehensive suite of solutions to help make your vision of patient care a reality. From technology to training, from reimbursement tools to patient education, Stryker is there to help.



Reimbursement

Today's healthcare environment is ever-evolving, and keeping up is critical for your practice. That's why we offer a variety of services and support tools to help you navigate the intricacies of reimbursement.

Contact our reimbursement hotline

Talk live with an expert to address your questions.

P 844 842 4232 | ENT-reimbursement@stryker.com



Flex Financial

Our payment options deliver smart alternatives designed to fit your business needs. We offer flexibility beyond a cash purchase – enabling you to secure the solutions you need for your practice, patients and staff. No other company provides strategic financial consulting combined with our breadth of products, meeting both your financial and equipment needs.

Contact your Sales Representative for more information about our flexible payment options.



Virtual training

Whether you are in the office, the OR or at home, our MedEd team is by your side. We've expanded our MedEd program in an effort to help you help your patients. We are committed to providing superior training on your schedule.



Live training

We understand there is still a demand for hands on training and peer to peer interaction. That's why we've created a combination of SITE visits, Area Labs, and virtual trainings to meet your needs.

ENT ordering information

To order, contact Customer Service | P 800 253 3210 | F 877 790 7772 | inst.stryker.cs@stryker.com

ENT

See Instructions for Use - visit: ent.stryker.com/physician-resources/instructions-for-use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Audion, Building Blocks, ClariFix, Cyclone, ESSx, FocESS, HemoPore, LATERA, Light Seeker, MiniFESS, NasoPore, OtoPore, Scopis, Stryker, XeroGel and XprESS. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

4278-018 rG 06/2022 | Copyright © 2022 Stryker

Stryker ENT
3600 Holly Lane North,
Suite 40
Plymouth, MN 55447
O 866 620 7615
F 763 463 1599

ent.stryker.com