

EasyClip[®]



Low profile. Dynamic compression.

EasyClip


Superelastic staple fixation system

Constant dynamic compression

EasyClip allows for active compression across the osteotomy site, avoiding pin tract infections and stress risers.¹


Compressive force²

| | |
|------------------|---------|
| EZ08, EZ10, EZ12 | 12 lbs. |
| EZ15, EZ18 | 17 lbs. |
| EZ20, EZ25 | 34 lbs. |



Patented forceps³

Reuseable instrumentation allows for simple removal, repositioning or adjusting staple sizes intra-operatively.



Forefoot

- 1 Proximal phalange
- 2 Distal metatarsal
- 3 Metatarsophalangeal joint

EasyClip size

- EZ 08, 10
- EZ 10, 12
- EZ 15, 18

Midfoot

- 4 Navicular cuneiform joint
- 5 Lisfranc joint
- 6 1st TMT joint

- EZ 15, 18
- EZ 15, 18
- EZ 18, 20

Rearfoot

- 7 Talonavicular joint
- 8 Calcaneocuboid joint
- 9 Subtalar joint

- EZ 20, 25
- EZ 20, 25
- EZ 20, 25

Visit footankle.stryker.com or email us at footandankle@stryker.com.

References

- Julie A. Neumann, MD, Kathleen D. Reay, MD, Kendall E. Bradley, BS, MS, and Selene G. Parekh, MD, MBA. "Staple Fixation for Akin Proximal Phalangeal Osteotomy in the Treatment of Hallux Valgus Interphalangeus." *Foot & Ankle International* 2015, Vol. 36(4) 457-464
- Data on File
- U.S. Pat. 8,137,351

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: EasyClip, Stryker. All other trademarks are trademarks of their respective owners or holders.

Manufactured by:
Stryker GmbH
Bohnackerweg 1
CH-2545 Selzach
Switzerland

www.stryker.com