

SALVATION™ 2

Midfoot Nail



SALVATION™ 2
Comprehensive Limb Salvage Solutions

MORE SIMPLE
MORE STABLE
MORE COMPREHENSIVE



 **WRIGHT™**
FOCUSED EXCELLENCE

SALVATION™ 2

Midfoot Nail

Superior Strength*

The Midfoot Nail is an 8mm intramedullary device that provides the strength of a nail with the ease of a beam.



Fixed Device with **8mm** of Controlled Compression

Cross-locking technology provides anti-migration capabilities while also being capable of providing up to 8mm of controlled compression across the fusion sites.

Simple Application

The cannulated technique and streamlined instrumentation allows for simple application in the most complex pathologies.



*Data on file.



1023 Cherry Road
Memphis, TN 38117
800 238 7117
901 867 9971
www.wright.com

56 Kingston Road
Staines-upon-Thames
Surrey TW18 4NL
United Kingdom
+44 (0)845 833 4435

161 Rue Lavoisier
38330 Montbonnot Saint Martin
France
+33 (0)4 76 61 35 00

™ and ® denote Trademarks and Registered Trademarks of Wright Medical Group N.V. or its affiliates.
©2019 Wright Medical Group N.V. or its affiliates. All Rights Reserved.

AP-011370A_14-Mar-2019