


# ToeTac<sup>®</sup> Xpress

## Hammertoe fixation system



# ToeTac Xpress

## Hammertoe fixation system

### Cannulated

- ▶ 0° cannulated option can be used with a k-wire to temporarily fixate the MTP joint\*

### PEEK


- ▶ Biocompatible, non-resorbable polymer
- ▶ Radiolucent implant allows for visualization of the fusion site

### Sterile

- ▶ Compact kit includes all instruments needs for procedure
- ▶ Sterile, single-use kits may reduce the risk of cross-contamination and post-operative infection, improving patient safety.

### Joint preparation

- ▶ Includes unique, sterile, cup and cone reamers designed for PIP joint preparation



Visit [footankle.stryker.com](http://footankle.stryker.com) or email us at [footandankle@stryker.com](mailto:footandankle@stryker.com).

\* 10° non-cannulated option also available.

#### References

1. Siegel, Geoffrey, MD, et. al. Cost analysis and surgical site infection rates in total knee arthroplasty comparing traditional vs. single-use instrumentation. The Journal of Arthroplasty 30 (2015) 2271-2274.

### Foot & Ankle

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, ToeTac. All other trademarks are trademarks of their respective owners or holders.

**Manufacturer:**  
Stryker GmbH  
Bohnackerweg 1  
2545 Selzach,  
Switzerland

[stryker.com](http://stryker.com)

This document is only approved for use in the USA.