

# Sonopet iQ

Ultrasonic Aspirator System

Bone cutting tips

## 12cm iQ Micro Claw

5500-25B-315  
4/box

**L** 4.8in (121mm) **CL** 1.6mm **CW** 2mm



## 12cm iQ Claw

5500-25B-318  
4/box

**L** 4.8in (121mm) **CL** 1.7mm **CW** 2.8mm



## 12cm iQ 120 Claw

5500-25B-317  
4/box

**L** 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm



## 12cm iQ Apex 360

5500-25B-312  
4/box

**L** 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm



## 20cm iQ Apex 360

5500-25B-302  
4/box

**L** 7.9in (200mm) **CL** 2.5mm **CW** 3.15mm



## 18cm iQ Micro Claw

5500-25B-303  
4/box

**L** 7.2in (182.8mm) **CL** 1.6mm **CW** 2mm



## 18cm iQ Apex Micro Claw

5500-25B-319  
4/box

**L** 7.2in (182.8mm) **CL** 1.6mm **CW** 2mm



## 11cm iQ Knife

5500-25B-107  
4/box


**L** 4.2in (106mm) **CL** 12.5mm **CW** 0.83mm



## 11cm iQ Apex Knife

5500-25B-114  
4/box

**L** 4.2in (106mm) **CL** 12.5mm **CW** 0.83mm



**L** = Length **CL** = Cutting length **CW** = Cutting width

**12cm iQ Micro**

5500-25S-309

4/box

**L** 4.7in (120mm) **OD** 1.6mm **ID** 1.20mm**12cm iQ Standard**

5500-25S-307

4/box

**L** 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Large**

5500-25S-308

4/box

**L** 4.7in (112.1mm) **OD** 2.49mm **ID** 2.0mm**20cm iQ Standard**

5500-25S-301

4/box

**L** 8.1in (205mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Barracuda**

5500-25S-310

4/box

**L** 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm**12cm iQ Serrated Large**

5500-25S-316

4/box

**L** 4.6in (115.7mm) **OD** 2.51mm **ID** 2.0mm**L** = Length **OD** = Outer diameter **ID** = Inside diameter**Neurosurgical**

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilisation (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Sonopet iQ and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

D0000011657 Rev. AA  
G55/PS  
Copyright © 2019 Stryker  
Printed in USA

Distributed by:  
Stryker Australia  
8 Herbert Street,  
St Leonards NSW 2065  
T: +61 2 9467 1000  
F: +61 2 9467 1010

Stryker New Zealand Limited  
PO Box 17-136  
Greenlane, Auckland 1546  
T: 64 9 573 1894  
F: 64 9 5731895

[www.StrykerMedEd.com](http://www.StrykerMedEd.com)