

stryker



Driving  
performance,

**delivering  
versatility**

**Spetzler-Malis and SILVERGlide  
Non-Stick Bipolar Forceps**

# Your **superior experience** starts with our people

## Supporting you throughout your procedures

### People

Our people stand out. They're dedicated, highly skilled experts that care deeply and listen carefully. When you need a partner in improving care, we're by your side.

### Culture

We create a dynamic, supportive environment that enables amazing people to do amazing things. Count on us to do what's right and deliver—every time.

### Product innovation

Innovation is in our DNA—from Dr. Homer Stryker to the researchers, engineers and surgeons creating solutions that are reshaping neurotechnology today.

### Medical education


Innovation and skill come together in our educational outreach to providers. Each year, we help surgical teams around the world master the latest medical technologies and care practices.

### Research and development


We amplify our investments in R&D by listening to our customers and collaborating with them to create solutions that are simply superior.

### Outcomes

We share your values and your commitment to innovative care that heals effectively and efficiently. With you, we're driven to make healthcare better.



All of our non-stick bipolar forceps **resist tissue buildup** and provide **excellent visualisation**. Regardless of your preference for single-use Spetzler-Malis or reusable SILVERGlide, you can expect **precision, control** and a responsive extension of your hands and eyes from us.



# Disposable Spetzler-Malis

## Standard Bipolar Forceps


**Flexible lightweight cord and remote connection plug** helps achieve balance

**Spring force** elevates control and tactile feedback

**Single use** ensures excellent tip alignment and coagulation performance

**Sleek bayonet design** improves visualisation

**Non-stick technology** enhances performance



### 7" (18cm)

- 6760-180-005 0.5mm tip
- 6760-180-010 1.0mm tip
- 6760-180-015 1.5mm tip

### 8" (20cm)

- 6760-200-005 0.5mm tip
- 6760-200-010 1.0mm tip
- 6760-200-015 1.5mm tip


### 9" (23cm)

- 6760-230-005 0.5mm tip
- 6760-230-010 1.0mm tip
- 6760-230-015 1.5mm tip

## Dual Irrigating Bipolar Forceps

**Transparent fine tubing** allows for better observation

**Dual irrigation design** maintains moisture of the tips and vessels



### 7" (18cm)

- 6770-180-005 0.5mm tip
- 6770-180-010 1.0mm tip
- 6770-180-015 1.5mm tip

### 8" (20cm)


- 6770-200-005 0.5mm tip
- 6770-200-010 1.0mm tip
- 6770-200-015 1.5mm tip

### 9" (23cm)

- 6770-230-005 0.5mm tip
- 6770-230-010 1.0mm tip
- 6770-230-015 1.5mm tip

## Slim Bipolar Forceps

**Slim profile** provides access in narrow fields



### 8" (20cm)


- 6780-200-005 0.5mm tip
- 6780-200-010 1.0mm tip
- 6780-200-015 1.5mm tip

### 9" (23cm)

- 6780-230-005 0.5mm tip
- 6780-230-010 1.0mm tip
- 6780-230-015 1.5mm tip

# SILVERGlide Non-Stick

## ProSeries Bipolar Forceps




**Standard two-pin round connector**  
works with most standard foot-switching bipolar electro-surgical generators

**Silver plating**  
reduces collateral heating, "hot spots" and sticking for performance

**Ultra-thin tines**  
allow better access in surgery

**Lightweight handle**  
for easy manipulation and balance

**Guide-stop**  
minimises scissoring and increases control



### 5.7" (14.5cm) bayonet

- 6730-145-004 0.4mm
- 6730-145-007 0.7mm
- 6730-145-010 1.0mm
- 6730-145-012 1.2mm
- 6730-145-015 1.5mm

### 6.3" (16cm) bayonet

- 6730-160-004 0.4mm
- 6730-160-007 0.7mm
- 6730-160-010 1.0mm
- 6730-160-012 1.2mm
- 6730-160-504, 45° up 0.4mm

### 7.5" (19cm) bayonet

- 6730-190-004 0.4mm
- 6730-190-007 0.7mm
- 6730-190-010 1.0mm
- 6730-190-012 1.2mm
- 6730-190-015 1.5mm
- 6730-190-020 2.0mm
- 6730-190-504, 45° up 0.4mm
- 6730-190-512, 45° up 1.2mm
- 6730-190-612, 45° down 1.2mm

### 8.6" (22cm) bayonet

- 6730-220-004 0.4mm
- 6730-220-007 0.7mm
- 6730-220-012 1.2mm
- 6730-220-015 1.5mm
- 6730-220-020 2.0mm
- 6730-220-107, 15° up 0.7mm
- 6730-220-307, 30° up 0.7mm
- 6730-220-507, 45° up 0.7mm
- 6730-220-607, 45° down 0.7mm
- 6730-220-512, 45° up 1.2mm
- 6730-220-612, 45° down 1.2mm

### 7.5" (19cm) straight bayonet

- 6735-190-007 0.7mm
- 6735-190-010 1.0mm
- 6735-190-310, 30° up 1.0mm
- 6735-190-320, 30° up 2.0mm



### 9.4" (24cm) bayonet

- 6730-240-004 0.4mm
- 6730-240-007 0.7mm
- 6730-240-012 1.2mm
- 6730-240-015 1.5mm
- 6730-240-020 2.0mm
- 6730-240-507, 45° up 0.7mm
- 6730-240-607, 45° down 0.7mm
- 6730-240-807, 80° up 0.7mm
- 6730-240-907, 80° down 0.7mm



### 11" (28cm) bayonet

- 6730-280-012 1.2mm
- 6730-280-015 1.5mm

## Keyhole Maximum Access Bipolar Forceps



### 5.5" (14cm) keyhole bayonet

- 6750-140-004 0.4mm
- 6750-140-007 0.7mm

### 6.7" (17cm) keyhole bayonet

- 6750-170-004 0.4mm
- 6750-170-007 0.7mm
- 6750-170-012 1.2mm

### 7.8" (20cm) keyhole bayonet

- 6750-200-004 0.4mm
- 6750-200-007 0.7mm
- 6750-200-012 1.2mm
- 6750-200-015 1.5mm

### 9.2" (23cm) keyhole bayonet

- 6750-230-007 0.7mm
- 6750-230-012 1.2mm
- 6750-230-015 1.5mm
- 6750-230-107, 15° up 0.7mm
- 6750-230-307, 30° up 0.7mm

### 10.2" (25cm) keyhole bayonet

- 6750-250-007 0.7mm
- 6750-250-012 1.2mm
- 6750-250-015 1.5mm
- 6750-250-507, 45° up 0.7mm
- 6750-250-607, 45° down 0.7mm
- 6750-250-807, 80° up 0.7mm
- 6750-250-907, 80° down 0.7mm

### 11.7" (28cm) keyhole bayonet

- 6750-280-012 1.2mm
- 6750-280-015 1.5mm

# SILVERGlide Non-Stick

## Precision Bipolar Forceps



### 6.1" (15.5cm) SP bayonet

6720-155-004 0.4mm  
6720-155-007 0.7mm  
6720-155-012 1.2mm  
6720-155-015 1.5mm

### 8.7" (22cm) SP bayonet

6720-220-004 0.4mm  
6720-220-007 0.7mm  
6720-220-012 1.2mm  
6720-220-015 1.5mm  
6720-220-020 2.0mm

### 7.0" (18cm) SP straight

6725-180-007 0.7mm  
6725-180-010 1.0mm  
6725-180-410, 30° up 1.0mm  
6725-180-420, 30° up 2.0mm

### 9.7" (24cm) SP bayonet

6720-240-004 0.4mm  
6720-240-007 0.7mm  
6720-240-015 1.5mm  
6720-240-020 2.0mm

### 7.5" (19cm) SP bayonet

6720-190-004 0.4mm  
6720-190-007 0.7mm  
6720-190-010 1.0mm  
6720-190-015 1.5mm  
6720-190-020 2.0mm  
6720-190-512, 45° up 1.2mm  
6720-190-612, 45° down 1.2mm

## Standard Bipolar Forceps



### 4.0" (10cm) jewelers

6715-100-004 0.4mm  
6715-100-007 0.7mm

### 7.5" (19cm) cushing

6710-190-004 0.4mm  
6710-190-007 0.7mm  
6710-190-010 1.0mm  
6710-190-012 1.2mm  
6710-190-015 1.5mm  
6710-190-020 2.0mm

### 5.0" (12.5cm) adson

6715-125-004 0.4mm  
6715-125-007 0.7mm  
6715-125-012 1.2mm

### 7.8" (20cm) scoville-greenwood

6710-200-015 1.5mm

### 6.0" (15cm) semkin

6715-150-012 1.2mm

### 8.3" (21cm) cushing

6710-210-007 0.7mm  
6710-210-012 1.2mm  
6710-210-015 1.5mm  
6710-210-020 2.0mm

## Irrigating Bipolar Forceps



Hollow tine design  
brings irrigation to the  
tip of the forceps

Separate fluid and  
electrical connectors  
accept a variety of cord-tubing sets



### 7.5" (19cm) bayonet

6740-190-007 0.7mm  
6740-190-012 1.2mm

### 8.7" (22cm) bayonet







6740-220-007 0.7mm  
6740-220-012 1.2mm

### 10.2" (25cm) bayonet

6740-250-007 0.7mm  
6740-250-012 1.2mm

## What sets SILVERGlide apart?

The elimination of tissue buildup allows you to work without interruption. These non-stick characteristics continue to perform use after use. The SILVERGlide alloy's level of thermal diffusivity is 37 times higher than that of stainless steel commonly used in forceps.

	Standard forceps	SILVERGlide Non-Stick	
"Hot spots" from sparks occur during coagulation			"Hot spots" start to occur during coagulation
The "hot spots" cause localised sticking			Heat is diffused rapidly to prevent sticking
Tissue seared to the forceps tips is ripped away			The forceps can be removed without tissue sticking

## Accessories

### SILVERGlide Bipolar Forceps

#### SILVERGlide Bipolar Forceps sterilisation trays

6700-002-000 Two forceps sterilisation tray  
6700-005-000 Five forceps sterilisation tray  
6700-010-000 10 forceps sterilisation tray (forceps not included)

#### Other SILVERGlide accessories

6700-000-002 Polishing cloths  
6700-000-001 Tip guards  
6702-000-000 Disposable cable  
6701-000-000 Reusable cable



### Disposable Spetzler-Malis Bipolar Forceps cord and irrigation sets

#### Rotary cord and tubing sets

6790-100-003 Bipolar irrigator integrated tubing and cord set, disposable  
6790-100-004 Bipolar irrigator integrated tubing and cord set, disposable, unitised plug

#### Gravity style cord and tubing sets

6790-100-001 Integrated tubing and bipolar cord set, disposable  
6790-100-002 Integrated tubing and bipolar cord set, disposable, unitised plug

#### Cord only sets

6790-100-005 Bipolar cord set, disposable  
6790-100-006 Bipolar cord set, disposable, unitised plug



## Malis Bipolar Generator and Irrigator system

### Malis Bipolar Generator and Irrigator system

6700-313-500 Malis Bipolar Generator (includes power cord, does not include foot pedal)  
6700-313-600 Malis Bipolar Irrigator (does not include power cord)  
0996-851-011 Universal power cord  
6700-313-004 Bipolar generator/irrigator interconnect cable  
6700-313-700 Bipolar generator cart



**Irrigator**  
Optional integrated irrigator

**Oversized color screen**  
for easy comprehension

**Programmable user profiles**  
store information for up to eight surgeons

**CUT and COAG knobs**  
facilitate quick setting adjustments

**Independent bipolar outputs**  
enhance functionality

**DUAL WAVE technology**  
provides confidence during cutting and coagulating

### Foot pedals

6700-313-003 Single foot pedal (COAG only, does not include required adapter cable)  
6700-313-304 Single foot pedal adapter cable  
6700-313-002 Dual foot pedal



**Whether you choose reusable or disposable,  
our non-stick bipolar forceps act as a specialised  
extension of your hands and eyes.**

For more information, contact your Neurosurgical Sales Representative.

## Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilisation (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Malis, SILVERGlide, Spetzler and Stryker. All other trademarks are trademarks of their respective owners or holders.

D0000053188 Rev. AA\_SSP\_NZ  
G55/PS  
Copyright © 2020 Stryker  
Printed in USA

Stryker Australia  
8 Herbert Street  
St.Leonards  
NSW 2065 Australia  
Ph: +61 2 9467 1000  
Fax: +61 2 9467 1010

Stryker New Zealand  
515 Mt. Wellington Highway  
Auckland 1060 New Zealand  
Ph: +64 9 573 1890  
Fax: +64 9 573 1891