

stryker



SpineMap Go™ augmented fluoroscopy and SpineMap® 3D navigation



**What's right
for you?**

Solutions to suit your style

From technology in our personal lives and at work, we all have a preferred way of doing things. Some are satisfied with a more utilitarian approach while others want all the bells and whistles. That's why we offer multiple solutions to meet your needs for:

- Individualized surgical workflows
- Preferred imaging devices
- Standard hospital rooms to fully integrated suites as well as outpatient surgical centers
- Various budget parameters

Choosing between augmented 2D or 3D technology is not about choosing what's right or wrong. Instead, it's about finding the visualization solution that's best for you. Here's a quick look at some core features.



	SpineMap Go	SpineMap 3D
Views / Imaging		
Fluoroscopic C-arm	✓	
3D imaging device or pre-op CT		✓
Biplanar fluoro views	✓	
AP/Lateral views	✓	✓
Pre-op CT/MRI merge		✓
Image manipulation	✓	✓
Image browser	✓	
Planning and placement		
Virtual screw visualization	✓	✓
Virtual screw preplanning/dropping		✓
Virtual screw families		✓
Measuring tool		✓
Annotation		✓
Positioning assistance		
C-arm guidance	✓	
Registration methods		
Automatic registration	✓	✓
Instant re-registration	✓	
Anatomical registration		✓
A.I.M. registration (with SpineMask Tracker)		✓
Basic operation		
Live camera for setup	✓	✓
Remote software control	✓	✓
Compatibility with:		
3D image set		✓
C-arm image set	✓	
NavLock tool tracker adapters	✓	✓
Noninvasive SpineMask Tracker	✓	✓
Stryker's navigated spine instruments	✓	✓
NAV3/3i Platform and NavSuite3 kit	✓	✓

Learn more

For a demo or more information, call your sales representative or **800 253 3210**. You may also learn more at **agt.stryker.com**.

Advanced Guidance Technologies

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: NAV3, NAV3i, NavLock, NavSuite, SpineMap, SpineMap Go, SpineMask and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

D0000010309 Rev. AA
G55/PS
Copyright © 2019 Stryker
Printed in USA

Stryker
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 323 2742
toll free: 800 253 3210

agt.stryker.com