

Sage

**AirTAP® LC**

**Lift Compatible Patient Repositioning System**

**Advancing mobility**  
with convenience  
and confidence



# Sage AirTAP<sup>®</sup> LC

## Lift Compatible Patient Repositioning System

### Promoting safety and mobility

We understand the importance of early mobility. To achieve your mobility protocols and comply with guidelines, you need a readily available solution that stays under your patient and fits seamlessly into your work flow. That's why we combined our innovative air assist technology with lift compatibility to promote safety and mobility.



**M<sup>2</sup> Microclimate Body Pad**  
Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

**Hi-tech fabric**  
Promotes a healthy microclimate for the skin and is compatible with low air loss technology  
Easy to clean surface



**Hose Protection Sleeve**

Helps protect hose from environmental contamination



**Point of care power switch**

Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside

**HEPA**  
Equipped filter



**Prevalon Air Pump**  
Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA filter is tested and certified to perform to HEPA standards and is **99.97% efficient in processing 0.3 μm particles.**

**Quick connect valve**  
Provides an easy, secure connection and a quick release



**30° Body and Anchor Wedge System**

Helps sustain a 30° side-lying position and reduces pressure by offloading sacrum  
Minimize manual activity as wedges slide into a turn, reducing the need for a log roll  
Minimizes patient migration down the bed  
Reduces the need for boosting, minimizing shear and friction



**Integrated head support**

Ensures patient's head maintains neutral position while inflated

**Easy roll handles**  
Encourage proper posture for healthcare worker

**Black positioning handles**  
Encourage proper ergonomic position for the healthcare worker when preparing to reposition patient



**Intuitive strap system**



Enables multiple lifting positions

## 3-in-1 mobility

- 1 In-bed positioning**  
Addresses risk factors of pressure injury
- 2 Lateral transfer**  
Utilizes air-assisted technology to safely transfer patients
- 3 Vertical transfer**  
Enables patients to safely transfer to the bedside chair



# Simple interventions. Extraordinary outcomes.

We are your partner for proven prevention. Our market-leading products solve real healthcare problems and are backed by clinical evidence. Our products allow you to deliver essential patient care with confidence by addressing risk factors that can lead to infections, skin injury, and caregiver injuries.

**We are driven to solve real problems and make healthcare better for you and your patients.**




Healthcare workers are critical to delivering essential care and improving patient outcomes. Keeping them safe should be a priority for every hospital. **The goal of the Journey to Zero program is to achieve zero harm for healthcare workers** – which helps improve care and makes hospitals safer for everyone. The program includes a range of repositioning and transfer devices along with educational services, assessments, and utilization tracking to help promote a culture of safety.

**We're committed to actively partnering with your facility to reduce the risk of patient handling injuries and improving outcomes.**

**Let's start the journey to zero – together.**

## Product ordering



37 in x 65 in  
Weight limit: 600 lb / 272 kg

### Sage AirTAP LC Lift Compatible Patient Repositioning System

- (1) Lift Compatible Glide Sheet
- (1) 30° Body and Anchor Wedge System
- (1) M<sup>2</sup> Microclimate Body Pad

5 systems/case  
Reorder #7287

### Sage AirTAP LC Lift Compatible Patient Repositioning System Glide Sheet

- (1) Lift Compatible Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

10 systems/case  
Reorder #7289

**30° Body and Anchor  
Wedge System**  
5 pair/case  
Reorder #7295

**M<sup>2</sup> Microclimate Body Pad**  
36" x 51"  
30 pads/case (6 bags of 5)  
Reorder #7250

Reference 1. Data on file, Sage Products LLC, IESTRP-CC007.3 testing conducted by APC Filtration Inc. Brentford, ON Canada. June 23-24, 2021.

Zero healthcare worker injuries and zero doubt messages are not guarantees and are aspirational in nature.