

stryker

Sage
HalfMATS™
Mobile Air Transfer System

Easily transfer patients
for OR specialty procedures



Sage

HalfMATS

Mobile Air Transfer System



Promotes easy transfer and positioning in the OR

We understand specialty procedures can have unique challenges. You shouldn't have to worry when you transfer surgical patients from surface to surface in the OR. That's why we designed HalfMATS to fit seamlessly into your workflow and address your transfer needs during specialty procedures so you can achieve the best possible outcomes.

Imaging capability


- Acceptable for use in MRI*, X-Ray, CT Scan, and ultrasounds¹



Operating room specialty procedures


Lithotomy positioning

- Gynecology
- Urology
- Cystoscopy



Orthopedic positioning

- Hip fracture
- Total hip replacement



Cradle inflation

- Cradles patient with 2-phase inflation and raises patient in one smooth motion

Quick Connect Valve

- Provides an easy, secure connection and a quick release

*Prevalon HalfMATS is MRI safe by rationale. The device is made from all non-metal materials; therefore MR safety testing was not performed. Compatibility tests did not show artifacts. Based on rationale, the HalfMATS is electronically non-conductive and non-magnetic.

Microclimate Management Pad

- Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

Easy Grab Handles

- Positioned along outer edges

Patient placement positioning (Breakline)


- Ensures proper placement of the patient to promote workflow efficiency

Prevalon Air Pump

- Easily fits into the Prevalon Air Pump Cart for transport or in cabinets and shelves in patient room

HEPA Equipped Filter

- Tested and certified to perform to HEPA standards and is 99.97% efficient in processing 0.3 μm particles²



Head placement indicator

- Provides easy guide to ensure proper head placement



Integrated head support

- Inflation provides comfortable support for the head and neck



Hose Protection Sleeve

- Helps protect hose from environmental contamination



Point of Care Power Switch

- Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside






Simple interventions. **Extraordinary outcomes.**

We are your partner for proven prevention. Our market-leading products solve real healthcare problems and are backed by clinical evidence. Our products allow you to deliver essential patient care with confidence by addressing risk factors that can lead to infections, skin injury, and caregiver injuries.

We are driven to solve real problems and make healthcare better for you and your patients.

Product **ordering**



HalfMATS **Mobile Air Transfer System**

(1) Mobile Air Transfer HalfMATS - 39" x 49"
(1) Microclimate Management Pad
Weight limit: 500 lbs / 227 Kg

10 systems/case
Reorder #3230



Microclimate Management Pad

23" x 36"
30 pads/case (6 bags of 5)
Reorder #7550



Prevalon Air Pump 120V
1/case
Reorder #7455



Prevalon Air Pump Cart
1 cart/case
Reorder #7475



Air Pump Bed Adapter
1/case
Reorder #7485



HEPA Equipped Replacement Filter
4/case
Reorder #7465



Hose Protection Sleeve (HPS)
50/case
Reorder #7460

References: 1. Data on file, Sage Products LLC. ECO16-0426 12.1.16. 2. Data on file. Sage Products LLC. MSC081 1660 PCS.