

Sage

PrimaFit™

External Urine Management for the Female Anatomy

1-step connection

Suction tubing adapter designed for easy connection to hospital tubing

Accommodates a range of body types


Long, flexible tubing designed to allow multiple suction tubing angles

Gentle adhesive pad

Keeps device in place and can be easily and gently removed

Helps keep skin dry

Ultra-soft fabric diverts urine away from skin




Flexible fit

Contours to patient anatomy and maintains shape for the duration of use

Stays in place

Soft, tapered end gently fits in perineal area



Helps promote early catheter removal

Address the #1 CAUTI risk factor, the indwelling catheter¹

PrimaFit is an external urine management system designed to help manage urinary incontinence in patients with female anatomy and promote early catheter removal, addressing the #1 risk factor of catheter-associated urinary tract infections (CAUTI). PrimaFit serves as an alternative to incontinence pads, diapers, or, at times, catheters while transitioning patients to independent continence.

When properly positioned, PrimaFit's ultra-soft fabric diverts urine away from the patient's skin. Urine is then absorbed into the system's core and suctioned into a collection canister, helping keep the skin dry and addressing a risk factor or incontinence-associated dermatitis (IAD).

Approximately

75%

of hospital-acquired UTIs are associated with an indwelling urinary catheter.¹



(1) PrimaFit, adult size/package

40 devices/case
Reorder #5400-X

References: 1. Centers for Disease Control and Prevention, Catheter-associated Urinary Tract Infections (CAUTIs). Available at: https://www.cdc.gov/hai/ca_uti/uti.html. Accessed August 24, 2017.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: PrimaFit, Sage, and Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

PrimaFit is a Class I medical device CE.

This document is not intended for distribution outside EU, UK and EFTA.

SAGE-UM-PSTR-862197_REV-0_en_us © 2024 Stryker



Sage Products LLC
3909 Three Oaks Road
Cary, IL 60013
800 323 2220

stryker.com