


Aleutian[®]

Anatomically-Narrow Lordotic Interbody System




The Aleutian Anatomically-Narrow (AN) Lordotic Interbody System has a full range of anatomically-designed PEEK interbodies for multiple spinal applications. The system offers a comprehensive solution for degenerative disc disease (DDD) up to a Grade I spondylothesis. This offers a solution to restore and maintain disc height and achieve sagittal balance in the lumbar spine. For all indications, the Aleutian implants are intended to be used with supplemental internal fixation appropriate for the implanted level.


Aleutian Anatomically-Narrow Lordotic Interbody System

Interbody design

8.5 x 24mm




8.5 x 28mm



Bilateral Aleutian AN Lordotic Interbodies pictured with the Mesa Degenerative Spinal System


Tantalum marker placement



- 8.5 x 24 and 8.5 x 28mm footprints available in posterior heights ranging from 4–12mm and lordosis angles of 6, 12, and 18°
- Posterior bone graft funnel provision
- Inserter built to support an insert and rotate technique
- Is designed to accommodate anatomic structure of endplates
- Manufactured of biocompatible PEEK polymer

Inserter design features

Inserter designed to facilitate in-situ rotation



Spine division

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. The products depicted are CE marked according to the Medical Device Directive 93/42/EEC. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Aleutian, K2M, Mesa, Stryker. All other trademarks are trademarks of their respective owners or holders.

ALEANL-SS-1_25379
Copyright © 2020 Stryker



K2M, Inc.
600 Hope Parkway SE
Leesburg, VA 20175
t: 571 919 2000
www.stryker.com