

stryker

Knotilus⁺

PEEK knotless suture anchor



Secure. Streamlined. Versatile.

Knotilus⁺

A streamlined, cross-specialty knotless anchor portfolio.

The Knotilus⁺ platform includes full procedural offerings with implants and corresponding instruments to facilitate procedures in the shoulder, hip and foot and ankle.¹




Secure

The anchor's cortical barb-lock technology is designed to optimize compression fixation in bone, supporting secure, knotless repairs.⁴

Streamlined

Knotilus⁺ is designed to improve procedural efficiency following insertion of the implant.²


The knob on the inserter handle is designed to provide tactile clicks to signal disengagement of the implant from the inserter with just four turns.^{2,3}



Versatile

Procedure-specific implants and instrumentation meet indication-specific needs.

Knotilus⁺ facilitates multiple stitch configurations such as mattress stitch, simple stitch and luggage stitch.



5 implant offerings



Knotilus+ 2.4x8.9mm PEEK suture anchor, hip

Knotilus+ 2.9x12.5mm PEEK suture anchor, hip

Knotilus+ 2.4x11.3mm PEEK suture anchor

Knotilus+ 2.9x12.5mm PEEK suture anchor

Knotilus+ 2.9x15.5mm PEEK suture anchor

19 instruments offerings

Available in two lengths: crown tip, fork tip, fish mouth and offset tip guides



Knotilus+ fork tip guide, hip



Knotilus+ crown tip guide, hip



Knotilus+ crown tip guide



Knotilus+ offset tip guide



Knotilus+ fishmouth tip guide



Drills to facilitate use in medium and hard bone

Part number	Description
Implants	
3911924511	Knotilus+ 2.4x11.3mm PEEK suture anchor
3911929512	Knotilus+ 2.9x12.5mm PEEK suture anchor
3911929515	Knotilus+ 2.9x15.5mm PEEK suture anchor
Instrumentation	
3911124511	Knotilus+ 2.4x11.3mm drill
3911224511	Knotilus+ 2.4x11.3mm drill, hard bone
3911129512	Knotilus+ 2.9x12.5mm drill
3911229512	Knotilus+ 2.9x12.5mm drill, hard bone
3911129515	Knotilus+ 2.9x15.5mm drill
3911229515	Knotilus+ 2.9x15.5mm drill, hard bone
3911600501	Knotilus+ crown tip guide
3911600502	Knotilus+ fishmouth tip guide
3911600503	Knotilus+ offset tip guide
3911000501	Knotilus+ bullet tip obturator, cannulated
3911000502	Knotilus+ pencil tip obturator

Additional part numbers coming soon.

References:

1. Knotilus+ IFU (P53789)
2. DHD14627
3. TR22906
4. TR22978

Sports Medicine

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Knotilus and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.