




stryker®

Sports Medicine

CHAMPION® *SlingShot* SUTURE PASSERS




E F F I C I E N T



The **Champion Slingshot** has been optimized for use in the shoulder and allows the surgeon to both pass and retrieve suture using a single portal. Its ergonomic thumb-slide design helps to eliminate additional procedural steps and potentially reduce operating time.

Optimized for the Shoulder

Ergonomically designed thumb-slide mechanism helps to simplify the suture passing and retrieval steps. The small Nitinol Arms allow for the ability to quickly retrieve the suture, while the razor sharp needle tip helps to minimize potential damage to the labral tissue.



Versatility

With a 45° Left, 45° Right, and 70° Up tip option, the Champion SlingShot is specifically designed to address various labral pathologies of the Shoulder.



PART NUMBER	DESCRIPTION
CAT02854	Champion SlingShot, 45° Left
CAT02855	Champion SlingShot, 45° Right
CAT02862	Champion SlingShot, 70° Up

5670 Greenwood Plaza Blvd. Ste. 200
Greenwood Village, CO 80111
t: 201 831 5000
www.stryker.com
www.sportsmedicine.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Champion, Stryker. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2015 Stryker
Printed in USA
1000-901-822 8/15

www.sportsmedicine.stryker.com

