



Cannulas

Product offerings



Stryker's comprehensive cannula offering includes **flexible, telescoping and rigid body** cannulas.

Dri-Lok



Adjustable inflow/outflow valve

- Transparent **rigid** body cannula available in threaded and non-threaded options
- Color coded sizes (5.0, 6.5 and 8.0mm) to accommodate a wide range of instrumentation and procedures

GateWay



- **Flexible** silicone cannula available in 6 lengths and 2 diameters, designed for easy insertion and improved function¹
- Increased stiffness of the distal flange is designed to reduce the risk of unintentional removal²
- Double dam designed to manage fluid
- Transparent color designed for improved visualization

Transport



- **Telescoping** cannula featuring in-situ length adjustment to accommodate most patient sizes
- Helps maintain optimal positioning and accounts for swelling in the joint space

Expanded size guide

Length 1, 2, 3:
60mm–80mm

Length 4, 5, 6:
80mm–110mm

Length 7, 8, 9:
110mm–140mm

Available in a variety of **lengths, thread patterns and diameters.**

GateWay silicone cannulas

Part number	Description
3910080020	8mm x 20mm GateWay silicone cannula
3910080030	8mm x 30mm GateWay silicone cannula
3910080040	8mm x 40mm GateWay silicone cannula
3910080050	8mm x 50mm GateWay silicone cannula
3910012030	12mm x 30mm GateWay silicone cannula
3910012040	12mm x 40mm GateWay silicone cannula

Transport cannulas

Part number	Description
CAT00222	Transport cannula, 08mm, length 1, 2, 3
CAT00223	Transport cannula, 08mm, length 4, 5, 6
CAT00224	Transport cannula, 08mm, length 7, 8, 9
CAT00792	Transport obturator, reusable
CAT01515	Stick, cannulated
CAT01514	Non-cannulated

Dri-Lok cannulas

Part number	Description
3910075500	5.0mm x 75mm threaded cannula
3910075501	5.0mm x 75mm non-threaded cannula
3910075502	5.0mm x 75mm Dri-Lok cannula distal thread
3910090500	5.0mm x 90mm Dri-Lok threaded cannula
3910090502	5.0mm x 90mm Dri-Lok cannula distal thread
3910075651	6.5mm x 75mm non-threaded cannula
3910075652	6.5mm x 75mm Dri-Lok cannula distal thread
3910090650	6.5mm x 90mm Dri-Lok threaded cannula
3910090652	6.5mm x 90mm Dri-Lok cannula distal thread
3910075800	8.0mm x 75mm threaded cannula
3910075801	8.0mm x 75mm non-threaded cannula
3910075802	8.0mm x 75mm Dri-Lok cannula distal thread
3910090800	8.0mm x 90mm Dri-Lok threaded cannula
3910090802	8.0mm x 90mm Dri-Lok cannula distal thread

Dri-Lok obturators

Part number	Description
3910075006	6.5mm x 75mm cannulated obturator
3910090006	6.5mm x 90mm cannulated obturator
3910075008	8mm x 75mm cannulated obturator
3910090008	8mm x 90mm cannulated obturator

References:

1. Improved function over rigid body cannula
2. Stryker DHD13873 Rev A 2019

Sports Medicine

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: GateWay, Dri-Lok, Transport and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

1000903189 Rev A

Copyright © 2019 Stryker

Stryker
Sports Medicine
5670 Greenwood Plaza Blvd.
Ste. 200
Greenwood Village, CO 80111
t: 866 596 2022
www.sportsmedicine.stryker.com