


# **XBraid TT, Force Fiber and XBraid S**



# XBraid TT



XBraid TT is a smooth, low profile, 100% ultra-high molecular weight polyethylene suture tape designed to increase tendon-to-bone contact in soft tissue repairs. The custom design of the XBraidTT allows the suture to lie flat against the tissue to increase the suture footprint.

## Features and benefits


- XBraid TT is a tie-able tape with comparable knot stacks to #2 suture.
- Designed with lubricity in mind, XBraid TT easily slides after anchor insertion and during knot tying.
- Easy to open packaging designed to help with efficiency in the operating room.
- Custom braid technology provides a suture tape that is specifically designed to not fold under tension.
- Ultimate tissue pull through load: 206N<sup>2</sup>

Product number	Description
3910900010	1.0mm XBraid TT, blue/black, suture tape, CT-3 (taper, 1/2c, 22mm), 28" suture
3910900044	1.0mm XBraid TT, blue/black, Passing Loop suture tape, non-needled, 24" suture with 7" loop
3910900130	1.0mm XBraid TT, blue, Passing Loop suture tape, non-needled, 24" suture with 1.5" loop
3910900017	1.2mm XBraid TT, violet co-braid, 100% UHMWPE
3910900045	1.2mm XBraid TT, white/blue, non-needled, 46" suture
3910900046	1.2mm XBraid TT, blue/black, non-needled, 46" suture
3910900047	1.2mm XBraid TT, white/black, non-needled, 46" suture
3910900048	1.2mm XBraid TT, white/blue, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900129	1.2mm XBraid TT, white/blue, Passing Loop suture tape, non-needled, 24" suture with 1.5" loop
3910900132	1.2mm DualBraid TT, white/blue+white/black, dual braid suture tape, non-needled, 40" suture
3910900002	1.4mm, white/blue, suture tape, non-needled, 40" overall length, 9" tape working length
3910900003	1.4mm, white/black, suture tape, non-needled, 40" overall length, 9" tape working length
3910900004	1.4mm XBraid TT, white/blue, suture tape, non-needled, 46" suture
3910900005	1.4mm XBraid TT, white/blue, suture tape, CT (taper, 1/2c, 40mm), 40" suture
3910900006	1.4mm XBraid TT, blue/black, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900009	1.4mm XBraid TT, white/blue, suture tape, ST-1 (straight taper needle, 60mm), on both ends, 40" suture
3910900011	1.4mm XBraid TT, blue/black, suture tape, CT (taper, 1/2c, 40mm), 40" suture
3910900012	1.4mm XBraid TT, white/blue, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900013	1.4mm XBraid TT, black/white, suture tape, non-needled, 40" suture
3910900043	1.4mm XBraid TT, white/blue, Passing Loop suture tape, non-needled, 24" suture with 5" loop
3910900064	1.4mm XBraid TT, white/blue, WhipStitch suture tape, ST-1 (straight taper needle, 60mm), 20" suture with 18" loop
3910900124	1.4mm XBraid TT, blue/black, WhipStitch suture, ST-1 (straight taper needle, 60mm), 18" suture loop
3910900125	1.4mm XBraid TT, white/black, suture tape, DA MO-6 (taper, 1/2c, 26mm), 40" suture
3910900128	1.4mm XBraid TT, white/blue, Passing Loop suture tape, non-needled, 24" suture with 1.5" loop
3910900134	1.4mm DualBraid TT, white/blue+white/black, suture tape, non-needled, 40" suture
3910900014	1.5mm XBraid TT, white/blue, suture tape, non-needled, 36" suture
3910900055	1.5mm XBraid TT, white/black, suture tape, non-needled, 40" suture

# XBraid TT



Product number	Description
3910900135	1.8mm DualBraid TT, white/blue+white/black, suture tape, non-needed, 40" suture
3910900142s	Double pack, 1.8mm XBraid TT, white/blue, suture tape, 9" flat, non-needed, 40" suture; 1.8mm XBraid TT, blue, suture tape, 9" flat, non-needed, 40" suture
3910900056	1.8mm XBraid TT, white/blue, suture tape, 9" flat, non-needed, 40" suture
3910900056s	1.8mm XBraid TT, white/blue, suture tape, non-needed, 40" suture with 9" flat length
3910900057	1.8mm XBraid TT, white/black, suture tape, 9" flat, non-needed, 40" suture
3910900058	1.8mm XBraid TT, white/blue, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900059	1.8mm XBraid TT, white/black, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900069	1.8mm XBraid TT, black/white, suture tape, da ST-1 (straight taper needle, 60mm), 40" suture
3910900018	2.0mm XBraid TT, violet co-braid, 100% UHMWPE
3910900019	2.0mm XBraid TT, violet, 100% UHMWPE
3910900007	2.2mm XBraid TT, white/blue, suture tape, CTX (taper, 1/2c, 48mm), 40" suture
3910900008	2.2mm XBraid TT, white/black, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900015	2.2mm XBraid TT, white/blue, suture tape, MO-6 (taper, 1/2c, 26mm), 40" suture
3910900016	2.2mm XBraid TT, white/black, suture tape, CTX (taper, 1/2c, 48mm), 40" suture
3910900144	Double pack, 2.2mm XBraid TT, white/black, suture tape, non-needed, 40" suture; 2.2mm XBraid TT, black, suture tape, non-needed, 40" suture



# XBraid S

The XBraid S Suture line offers high-strength, round sutures in various sizes and colors. Sutures are composed of 100% ultra-high molecular weight polyethylene and designed to increase tendon-to-bone contact in soft tissue repairs.

## Features and benefits

- Designed with lubricity in mind, XBraid S easily slides after anchor insertion and during knot tying.
- Easy to open packaging designed to help with efficiency in the operating room.

Product number	Description
3910900039	#4-0 XBraid S, white/blue, suture, BB (taper, 3/8c, 17mm), 18" suture
3910900063	#4-0 XBraid S, white/blue, WhipStitch suture, BB (taper, 3/8c, 17mm), 10" loop length
3910900122	#4-0 XBraid S, white/blue, WhipStitch suture, ST-1 (straight taper needle, 60mm), 10" suture loop
3910900038	#3-0 XBraid S, white/blue, suture, BB (taper, 3/8c, 17mm), 18" suture
3910900036	#2-0 XBraid S, white/blue, suture, CT-3 (taper, 1/2c, 22mm), 18" suture
3910900037	#2-0 XBraid S, white/black, suture, non-needed, 18" suture
3910900062	#2-0 XBraid S, white/blue, WhipStitch suture, CT-3 (taper, 1/2c, 22mm), 20" loop length
3910900072	#2-0 XBraid S, white/blue, XPass suture, non-needed, 50" suture (one end stiffened 12")
3910900121	#2-0 XBraid S, white/blue, WhipStitch suture, ST-1 (straight taper needle, 60mm), 20" suture loop
3910900127	#2-0 XBraid S, white/blue, suture, non-needed, 36" suture
3910900042	#2-0 XBraid S, white/blue, Passing Loop suture, non-needed, 24" suture with 1.5" loop
3910900034	#0 XBraid S, white/blue, suture, CT-3 (taper, 1/2c, 22mm), 39" suture
3910900035	#0 XBraid S, white/black, suture, non-needed, 39" suture
3910900041	#0 XBraid S, white/blue, Passing Loop suture, non-needed, 24" suture with 1.5" loop
3910900061	#0 XBraid S, white/blue, WhipStitch suture, CT-3 (taper, 1/2c, 22mm), 20" loop length
3910900110	#0 XBraid S, white, shuttle loop suture, ST-1 blunt (straight blunt point needle, 60mm), 22" suture with 1" Brummel center loop
3910900111	#0 XBraid S, white, shuttle loop suture, non-needed, 22" suture with 1" Brummel center loop
3910900120	#0 XBraid S, white/blue, WhipStitch suture, ST-1 (straight taper needle, 60mm), 20" suture loop
3910900033	#1 XBraid S, white/black, suture, CT-3 (taper, 1/2c, 22mm), 46" suture
3910900126	#1 XBraid S, white/blue, suture, non-needed, 46" suture
3910900131	#1 DualBraid S, blue+blue/black, suture, non-needed, 46" suture
3910900023	#2 XBraid S, blue/black, suture, MO-6 (taper, 1/2c, 26mm), 39" suture
3910900025	#2 XBraid S, white/blue, suture, MO-6 (taper, 1/2c, 26mm), 39" suture
3910900026	#2 XBraid S, white/black, suture, CT taper, 1/2c, 40mm), 39" suture
3910900027	#2 XBraid S, white/blue, suture, non-needed, 9" suture
3910900028	#2 XBraid S, white/black, suture, non-needed, 39" suture
3910900029	#2 XBraid S, white/blue, suture, ST-1 (straight taper needle, 60mm) on both ends, 39" suture
3910900030	#2 XBraid S, white/blue, suture, CT (taper, 1/2c, 40mm), 39" suture
3910900040	#2 XBraid S, white/blue, Passing Loop suture, non-needed, 24" suture with 1.5" loop
3910900060	#2 XBraid S, white/blue, WhipStitch suture, ST-1 (straight taper needle, 60mm), 20" suture with 18" loop


# XBraid S



Product number	Description
3910900070	#2 XBraid S, white/blue, XPass suture, non-needled, 50" suture (one end stiffened 12")
3910900071	#2 XBraid S, white/black, XPass suture, non-needled, 50" suture (one end stiffened 12")
3910900076	#2 XBraid S, white/blue, XPass suture with loop, non-needled, 26" suture (one end stiffened 12", other end .75")
3910900112	#2 XBraid S, white/blue, XPass shuttle loop suture, non-needled, 32" suture with 1" brummel center loop
3910900113	#2 XBraid S, blue/black, XPass shuttle loop suture, non-needled, 32" suture with 1" center loop
3910900116	#2 XBraid S, blue, suture, non-needled, 39" suture
3910900123	#2 XBraid S, blue/black, WhipStitch suture, ST-1 (straight taper needle, 60mm), 18" suture loop
3910900133	#2 DualBraid S, blue+blue/black, dual braid suture, non-needled, 39" suture
3910900143	Double pack, #2 XBraid S, white/blue, suture, non-needled, 39" suture; #2 XBraid S, white/black, suture, non-needled, 39" suture
3910900052	#5 XBraid S, white/black, suture, CTX (taper, 1/2c, 48mm), 39" suture
3910900117	Needle XBraid S, MO-6 (taper, 1/2c, 26mm), 10mm loop
3910900118	Needle XBraid S, CTX (taper, 1/2c, 48mm), 10mm loop
3910900119	Needle XBraid S, ST-1 (straight taper needle, 60mm), 10mm loop



# Force Fiber



Force Fiber is a coreless, surgical suture designed with a low knot profile and better knot security than comparative sutures with a core. Unlike polyblend and some polyester sutures, Force Fiber is composed of a proprietary combination of materials that allow for a unique braid construction, designed to contribute to its strength. It is a non-absorbable suture composed of ultra-high molecular weight, polyethylene fiber allowing for a durable and strong solution.

## Features and benefits

- 53% better tensile strength than the leading competitor.<sup>1\*</sup>
- Designed to offer greater tensile strength than polyester and some polyblend sutures, Force Fiber is suitable/optimal for ligation of soft tissue for orthopedic applications.
- Lubricious for smoother, more efficient performance.
- Blue co-braid design for increased visibility in the joint.
- Available with taper, cutting, or K-point needles.

Product number	Description
3910900022	#2 Force Fiber, white/blue, suture, KC-7 (k-point, 1/2c, 38.9mm), 38" suture
3910900024	#2 Force Fiber, white/blue, suture, C-7 (taper, 1/2c, 38.9mm), 38" suture
3910900021	#2 Force Fiber, white/blue, suture, non-needed, 38" suture
3910900020	#2 Force Fiber, white/blue, suture, HC-5 (taper, 1/2c, 25.9mm), 38" suture
Cat01585	#2 Force Fiber, white/green, suture, non-needed, 36" suture
Cat01586	#2 Force Fiber, white/black, suture, non-needed, 36" suture
3910900051	#5 Force Fiber, white/blue, suture, non-needed, 38" suture
3910900050	#5 Force Fiber, white/blue, suture, K-60 (3/8c, sternotomy, 49.2mm), 38" suture



**References:**

1. Barber FA et al: Sutures and Suture Anchors – Update 2006. Arthroscopy 22:1063, 2006. \*Arthrex
2. Test Report: DHFD12758

## Sports Medicine

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Force Fiber, Stryker and XBraid. All other trademarks are trademarks of their respective owners or holders.