

ICONIX


All-Suture Anchor Platform



Features & Benefits


IntelliBraid Technology

The IntelliBraid Technology creates a bunching effect using targeted compression zones within the implant sheath for optimal fixation.



Minimal Bone Removal

The ICONIX disposable drills feature a unique Self Centering (SC) Technology, to help ensure accurate pilot hole placement.



Instrumentation

Drill Guide Options (2.3mm & 1.4mm):

- Straight Short
- Straight
- 12°
- 12° Slant
- 25°
- 25° Slant
- Disposable Kit Options




Anchor Options XBraid TT and Force Fiber


2.3mm Self-Punching Anchors



2.3mm Anchors



1.4mm Anchors



Part Number	Description
ICONIX	
3910-500-512	ICONIX 1 WITH INTELLIBRAID TECHNOLOGY, 1.4MM ANCHOR WITH 1 STRAND #2 FORCE FIBER (BOX OF 5)
3910-500-522	ICONIX 2 WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 2 STRANDS #2 FORCE FIBER (BOX OF 5)
3910-500-525	ICONIX 25 WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 2 STRANDS #5 FORCE FIBER (BOX OF 5)
3910-500-532	ICONIX 3 WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 3 STRANDS #2 FORCE FIBER (BOX OF 5)
ICONIX TT	
3910-500-312	ICONIX 1 TT WITH INTELLIBRAID TECHNOLOGY, 1.4MM ANCHOR WITH 1.2MM XBRAID TT (BOX OF 5)
3910-500-322	ICONIX 2 TT WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 2.0MM XBRAID TT (BOX OF 5)
ICONIX Needles	
3910-500-212	ICONIX 1 NEEDLES WITH INTELLIBRAID TECHNOLOGY, 1.4MM ANCHOR WITH 1.4MM XBRAID TT (BOX OF 5)
3910-500-412	ICONIX 1 NEEDLES WITH INTELLIBRAID TECHNOLOGY, 1.4MM ANCHOR WITH 1 STRAND #2 FORCE FIBER (BOX OF 5)
3910-500-222	ICONIX 2 NEEDLES WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 2.3MM XBRAID TT (BOX OF 5)
3910-500-422	ICONIX 2 NEEDLES WITH INTELLIBRAID TECHNOLOGY, 2.3MM ANCHOR WITH 2 STRANDS #2 FORCE FIBER (BOX OF 5)
ICONIX SPEED	
3910-500-921	ICONIX SPEED WITH 1 STRAND 1.2MM & 1 STRAND 2.0MM XBRAID TT
3910-500-931	ICONIX SPEED TRIPLE LOADED 1.2MM XBRAID TT
3910-500-920	ICONIX SPEED DOUBLE LOADED 2.0MM XBRAID TT
3910-500-922	ICONIX SPEED DOUBLE LOADED #2 FORCE FIBER

5670 Greenwood Plaza Blvd. Ste. 200
Greenwood Village, CO 80111
t: 201 831 5000
www.stryker.com
www.sportsmedicine.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks; ICONIX, XBraid, Stryker. All other trademarks are trademarks of their respective owners or holders.

1000902283 Rev C
Copyright © 2018 Stryker
Printed in USA