



Tray layout card – Gamma4 Indication Tray

Insert layer





Ref #	Description	Qty.
A 1420-1010	Indication Tray Insert	1
1 1420-0300	Extraction Shaft	1
2 1420-1020	Free Space Mat	1
3 1420-0192	Anti-Rotation Sleeve	1
4 1420-0190	Anti-Rotation Clip	1
5 1420-2100	Intermediate Targeting Sleeve	1
6 1420-3020	RC Inserter	1
7 1420-3010	RC Guide	1

Ref #	Description	Qty.
Optional		
8 1320-0175	Final Impactor	1
9 2351-0240	One Shot Device	1
Spare Part		
1500-0040	Full Size Tray Lid	1

Tray layout card – Gamma4 Indication Tray



Standard base layer



Ref #	Description	Qty	Ref #	Description	Qty
B 1420-1005	Indication Tray Base	1	9 1420-0160	Lag Screw Sleeve	1
1 2351-6000	Opening Reamer Handle	1	10 1420-0260	Lag Screw Driver (two parts)	1
2 1420-0050	Opening Sleeve	1	11 1420-0270	Set Screw Driver	1
3 1420-0055	Opening Guide	1	12 1420-0240	Lag Screw Reamer	1
4 1420-0080	Opening Reamer	1			
5 1420-0110	Targeting Sleeve (two parts)	1	Spare part		
6 1420-0100	Proximal Targeting Arm	1	1420-0112	Targeting Sleeve Knob	1
7 1420-0105	Nail Holding Screw	2	Indication kit		
8 1420-0220	Precision Sleeve	1	1420-0000	Indication Kit	1

Tray layout card – Gamma4 Indication Tray

Plus base layer



Ref #	Description	Qty	Ref #	Description	Qty
B 1420-1005	Indication Tray Base	1	9 1420-0160	Lag Screw Sleeve	1
1 2351-6000	Opening Reamer Handle	1	10 1420-0260	Lag Screw Driver (two parts)	1
2 1420-0050	Opening Sleeve	1	11 1420-0270	Set Screw Driver	1
3 1420-0055	Opening Guide	1	12 1420-0240	Lag Screw Reamer	1
4 1420-0080	Opening Reamer	1			
5 1420-0110	Targeting Sleeve (two parts)	1	Spare part		
6 1420-0120	Plus Targeting Arm	1	1420-0112	Targeting Sleeve Knob	1
7 1420-0125	Plus Nail Holding Screw	2	Indication kit		
8 1420-0220	Precision Sleeve	1	1420-0005	Plus Indication Kit	1

This document is intended solely for healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area. Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Gamma4, Precision Sleeve, Stryker. All other trademarks are trademarks of their respective owners or holders.



Manufactured by:
 Stryker GmbH
 Bohnackerweg 1
 2545 Selzach, Switzerland
 Phone: +41 (0) 32 641 66 66
 Fax: +41 (0) 32 641 66 60