

# Count sheet

## T2 Alpha Femur GT Antegrade tray

The instrument and implant banks must be ordered separately. This sheet contains all products you will need for the case.

### Base

Ref #	Description	Location	Qty.	Check box
2353-3103	Targeting arm femur GT	Tray	1	
2351-0240	One shot device	Tray	1	
2353-3104	Nail holding screw tibia / femur GT	Tray	2	
2351-0221	Recon lag screwdriver rod	Tray	1	
2355-5030	Drills and pins insert	Tray	1	

### Opening 13mm insert

Ref #	Description	Location	Qty.	Check box
2351-6113	Opening reamer	One tray level	1	
2351-6000	Opening reamer handle	One tray level	1	
2351-6213	Opening reamer sleeve	One tray level	1	
2351-6313	Opening reamer trocar ø13	One tray level	1	
2351-0145	Quick-lock T-handle, modified trinkle	One tray level	1	

### Femur antegrade GT / PF indication insert

Ref #	Description	Location	Qty.	Check box
2351-0210	Recon lag screw drill	One tray level	1	
2351-0200	Recon K-wire sleeve	One tray level	2	
2351-0190	Recon tissue protection sleeve	One tray level	2	
2353-3301	Recon knob	One tray level	1	
2351-0225	Recon lag screwdriver bit	One tray level	1	
2351-0250	External compression device femur	One tray level	1	

\* Please note that the instrument bank and implant set are separate part numbers and must be ordered separately.

Hospital signature \_\_\_\_\_

